



LAVISION

WE COUNT ON PHOTONS

NanoStar

family of intensified
12 bit camera systems
for imaging and
spectroscopy.

NanoStar incorporates
a high resolution
fast-scan CCD,
very fast readout,
and low noise.

SPECIAL

Fast electronic
CCD-shutter on-chip
(superior to MCP
bracket-pulsing
or MCP gating)

- ▶ **Image Intensifier:** gate width to 5 ns, repetition rate to 2 MHz, UV – NIR with single photon sensitivity
- ▶ **CCD:** 12 bit, 12.5 MHz readout, 1280 x 1024 pixels or 1340 x 1024 pixels, horizontal and vertical binning
- ▶ **Interface:** PCI board connected with fiber optic cable (10-1500 m)
- ▶ **DaVis software:** image acquisition and analysis, camera control, synchronization and control of external devices, programming of all features in 'C' syntax
- ▶ **NEW Double frame mode:** 2 full frames with interframing time down to 500 ns



Applications

- ▶ Combustion Research
- ▶ Spray Diagnostics
- ▶ Fluorescence Imaging and Spectroscopy
- ▶ Engine Diagnostics
- ▶ Flow Analysis
- ▶ Raman Imaging
- ▶ Plasma Diagnostics
- ▶ Laser Ablation

- ▶ **Ultrafast Readout:** 8 frames /s (1280 x 1024 pixels)
High System Dynamic: > 2000 : 1 @ single shot
- ▶ **Compact Design:** all controllers built in camera head
- ▶ **Test Certificate:** detection efficiency (counts / photoelectron), gain curve, single or double shot linearity, dark noise
- ▶ **Options:** spectrograph, stepping motor, chiller, lenses, stereoscope, filters, additional analogue input, software packages for device control (spectrographs, lasers, stepping motors)

LAVISION GMBH

ANNA-VANDENHOECK-RING 19 / D-37081 GOETTINGEN / GERMANY
E-MAIL: INFO@LAVISION.DE / WWW.LAVISION.DE
TEL. +49-(0)5 51-90 04-0 / FAX +49-(0)5 51-90 04- 100

LAVISION INC.

301 W. MICHIGAN AVE. / SUITE 403 / YPSILANTI, MI 48197 / USA
E-MAIL: SALES@LAVISIONINC.COM / WWW.LAVISIONINC.COM
PHONE: (734) 485 - 0913 / FAX: (240) 465 - 4306

NanoStar



NanoStar is a turn-key system.

It consists of camera head including control units (for CCD and image intensifier), PC, PCI interface and TTL-I/O board. The

DaVis software integrates these devices into a single system. This allows immediate implementation and synchronization to the application.

The **image intensifier** is proximity focused with single stage MCP. Different photo cathodes are available with 18 or 25 mm diameter.

The **DaVis software** controls all camera functions (gate width, delay, gain, binning, windowing) and performs accurate timing. By use, the user can quickly acquire a data series. Measurement and evaluation series can also be user programmed with 'C' syntax Command Language.

Note:

Above data is for NanoStar camera family. For the specific models see corresponding technical data sheets.

Data provided by LaVision is believed to be true. However, no responsibility is assumed for possible inaccuracies or omissions. All data are subject to change without notice.

► A. General System Specifications

Image intensifier gate	_____	down to 5 ns
System dynamic	_____	> 2000:1 @ single shot
Sensitivity	_____	typical: > 80 counts / photoelectron
Spectral range	_____	190 – 900 nm
Double frame	_____	2 full frames with interframing time down to 500 ns
On-chip CCD shutter	_____	down to 100 ns electronic CCD shutter

► B. CCD Chip and Control Unit

Sensor	_____	SuperVGA
Number of pixels	_____	1280 x 1024 pixels
Pixel size	_____	6.7 x 6.7 μ m
Full-well capacity	_____	25,000 electrons
Frame rate	_____	8 frames/s
Readout rate	_____	12.5 MHz
Readout noise	_____	< 2 counts RMS @ 12.5 MHz
A/D converter	_____	12 bit @ 12.5 MHz
Binning	_____	horizontal 1-8, vertical 1-32
Cooling type	_____	2-stage Peltier, forced air (optional liquid)

► C. Image Intensifier and Control Unit

Design	_____	proximity focused, single stage MCP
Photocathode	_____	18 mm or 25 mm, S20, S25, GaAs or GaAsP
Gate width	_____	5 ns – 1000 s
Delay	_____	0 ns - 1000 s, with minimum 45 ns intrinsic delay
Jitter	_____	< 0.5 ns (< 5 ns for gates \geq 100 ns)
Repetition rate	_____	to 2 MHz (3 kHz for gates < 20 ns)
Phosphor	_____	P43 or P46 (decay to 10% in 1 ms or 0.3 μ s)
Coupling to CCD	_____	lens optics

► D. Personal Computer, PCI Interface and TTL-I/O Board

state-of-the-art PC, monitor, TTL-I/O interface, PCI interface with fiber optic cable.

► E. DaVis software

Operating System	_____	Windows XP
Data acquisition	_____	Camera, timing, user programmable
Data processing	_____	buffer, column, row or pixel addressing
Command Language	_____	'C' syntax for system and user functions
Communication	_____	RS 232 (GPIB optional), TTL-I/O board

► F. Options

Spectrograph	_____	various types with adapter to camera head
Device control	_____	spectrograph, stepping motor, laser
Optics	_____	lenses (e.g. UV-Nikkor), telescopes, filters, long distance microscopes (Questar), macro lenses
Upgrades	_____	LIF, Rayleigh, Raman, Spray, ...
Analogue input	_____	simultaneous recording (e.g. single shot control)

LA VISION GMBH

ANNA-VANDENHOECK-RING 19 / D-37081 GOETTINGEN / GERMANY
E-MAIL: INFO@LAVISION.DE / WWW.LAVISION.DE
TEL. +49-(0)5 51- 90 04- 0 / FAX +49-(0)5 51- 90 04- 100

LA VISION INC.

301 W. MICHIGAN AVE. / SUITE 403 / YPSILANTI, MI 48197 / USA
E-MAIL: SALES@LAVISIONINC.COM / WWW.LAVISIONINC.COM
PHONE: (734) 485 - 0913 / FAX: (240) 465 - 4306