


# FlowMaster PIV system in Application


## Large scale digital PIV applied to landslide generated impulse waves

**Reference:** Hermann M. Fritz, Laboratory of Hydraulics, Hydrology and Glaciology (VAW), Swiss Federal Institute of Technology (ETH),  
**Source:** 10th Int. Symp. on Applications of Laser Techniques to Fluid Mechanics 38p1, Lisbon, Portugal, July 2000


Digital large scale PIV was applied to an extremely unsteady three phase flow consisting of granular matter, air and water. Areas of interest up to 0.8 m by 0.8 m are investigated in the impulse wave generation zone.



Impact experiments were conducted in a rectangular prismatic water wave channel (LxWxH: 11m, 0.5m, 1m). Landslides were modelled with an artificial granular material (PP-BaSO<sub>4</sub>).




Experimental setup: LaVision's FlowMaster 2 CCD camera (Kodak -ES1.0), and a twin cavity Nd:YAG laser were used for PIV.




Velocity vector field with flow reattachment and characteristic saddle point above slide shoulder. hill slope angle= 45°, stillwater depth= 0.45 m, slide mass= 27 kg, impact velocity= 2.5 m/s

## PIV-investigation of the in-cylinder tumble flow in an IC-engine

**Reference:** W. Hentschel, B. Block, Research and Development, Volkswagen AG, Wolfsburg, Germany; **Sources:** Application of laser-optical diagnostics for the support of direct-injection gasoline combustion process development 4. Int. Symp. für Verbrennungsdagnostik, Baden-Baden, May 2000; PIV-investigation of the in-cylinder tumble flow in an IC-engine. 8th. Int. Conf. on Anemometry, Rom, Sept.1999




**PIV setup on the transparent engine.**  
Arrangement of the LaVision FlowMaster 3S PIV camera (left) and the light sheet optic (bottom)



### PIV setup on the glass ring engine

**Nd:Yag double-pulse laser:** 532 nm, ca. 2 x 20 - 50 mJ, pulse separation: 4 - 8  $\mu$ s, **light sheet width:** ca. 1 mm  
**seeding:** oil droplets, **double-image CCD camera:** resolution 1024 x 1280 pixels, f = 85 mm / aperture 8 / 11, with 532 nm interference filter, **evaluation:** 2D-cross correlation  
**interrogation cells:** 32 x 32 pixels = 0,8 x 0,8 mm<sup>2</sup> → 32 x 40 vectors, **field of view:** 32 x 26 mm<sup>2</sup>, in the cylinder centre




**Intake and compression flow**  
PIV-measurement at a 4-valve SI-engine, central cross section

# Flow Field Measurement of Unsteady Flows in Radial Diffusors using PIV


**Reference:** U. Lohmann, Energy and Power Station Technology, University of Essen, Germany

The operational characteristic of a blower is not exclusively depending on the shape of the rotor, also the stator has great influence on the maximum efficiency. Instantaneous PIV-measurements can be used for optimization of the diffuser flow and thus optimization of the blower performance.



Experimental setup using LaVision's FlowMaster 3S PIV system

Setup of blower and diffuser. Aperture angle of the diffuser can be varied.




Diffuser flow field at design operation (left image, averaged). Velocity at the end of the diffuser is reduced about a factor of 2 compared to the entry region of the diffuser. Most of the kinetic energy is converted into pressure. Diffuser flow field with high restriction and high aperture angle of diffuser (right image). Increased turbulence and backward flow field areas are visible.


# Highly-Resolved Three-Dimensional Velocity Measurements via Dual-Plane Stereo Particle Image Velocimetry (DSPIV) in Turbulent Flows

**Reference:** John A. Mullin and Werner J.A. Dahm, Laboratory for Turbulence & Combustion (LTC), Department of Aerospace Engineering, University of Michigan **Source:** AIAA conference 2002-0290, Reno, NV, USA, Jan. 2002


A frequency-based dual-plane stereoscopic particle image velocimetry (DSPIV) technique is used to obtain fully-resolved simultaneous measurement of all nine components of the instantaneous velocity gradient tensor field  $\nabla u(x,t)$  at the small scales of a turbulent flow. The technique is based on two essentially independent stereo PIV systems that provide three-component velocity measurements in two differentially-spaced light sheets of different colors.



Asymmetric forward-forward scatter arrangement of the four FlowMaster 3S cameras, showing one camera in each color-pair oriented at 20-degrees and the other at 30-degrees relative to the sheet-normal.



Schematic indicating layout of coflowing turbulent jet facility used in this study.



Three-component velocity vectors in the green (left) and red (right) light sheet planes obtained in the centerline of the coflow and turbulent jet at Reynolds Number 6000. Laser sheet thickness: 800  $\mu\text{m}$ ; sheet separation: 400  $\mu\text{m}$

LA VISION GMBH

ANNA-VANDENHOECK-RING 19 / D-37081 GOETTINGEN / GERMANY  
E-MAIL: INFO@LAVISION.DE / WWW.LAVISION.DE  
TEL. +49-(0)5 51-90 04-0 / FAX +49-(0)5 51-90 04-100

LA VISION INC.

301 W. MICHIGAN AVE. / SUITE 403 / YPSILANTI, MI 48197 / USA  
E-MAIL: SALES@LAVISIONINC.COM / WWW.LAVISIONINC.COM  
PHONE: (734) 485-0913 / FAX: (248) 465-4306