

METALS

and Related Products

Reference Materials
Catalogue 2008



Excellence through measurement

LGC Standards office contact details:

Bulgaria

LGC Standards Sp. z o.o.
Офис 102
Бул. "Г.М. Димитров" 52А
1125 София, България

Tel: +359(0) 2 971 4955
Fax: +359(0) 894 780 154
Email: bg@lgcstandards.com

Czech Republic

LGC Standards Sp. z o.o.
Organizační složka
Hněvkovského 65
Brno
617 00
Česká republika

Tel: +420 543 529 205
Fax: +420 543 529 205
Email: cz@lgcstandards.com

Finland

LGC Standards AB
Albanoliden 5, 4 trappor
Box 1737
SE-501 17 Borås
Sweden

Tel: +46 (0)33 20 90 60
Fax: +46 (0)33 20 90 79
Email: se@lgcstandards.com

France

LGC Standards SARL
6, Rue Alfred Kastler
B.P. 83076
67123 Molsheim
France

Tel: +33 (0)3 88 04 82 82
Fax: +33 (0)3 88 04 82 90
Email: fr@lgcstandards.com

Germany

LGC Standards GmbH
Mercatorstr. 51
46485 Wesel
Germany

Tel: +49 (0)281 9887 0
Fax: +49 (0)281 9887 199
Email: de@lgcstandards.com

Hungary

LGC Standards GmbH
Magyarországi Közvetlen
Kereskedelmi Képviselet
H-2000 Szentendre
Körte u. 4.
Hungary

Tel/Fax: +36 (06) 26 314 891
Email: hu@lgcstandards.com

Ireland

4 Manor Lodge
Magherafelt
County Derry
BT45 6QL
Ireland

Tel: +44 (0)28 7930 0078
Mobile: 07921 872595
Email: patrick.henry@lgcstandards.com

India

LGC Promochem
PO Box 8061
No. 142, 3rd Floor, 5th Cross
Rajmahal Vilas Extension
Bangalore 560 080
India

Tel: +91-080-2361 4774
Fax: +91-080-2361 3859
Email: in@lgcpromochem.com

Italy

LGC Standards S.r.l.
Via Venezia, 23
20099 Sesto San Giovanni
Italy

Tel: +39 02 2412 6830
Fax: +39 02 2412 6831
Email: it@lgcstandards.com

Poland

LGC Standards Sp. z o.o.
ul. M. Konopnickiej 1
Dziekanów Lesny
05-092 Lomianki
Poland

Tel: +48 (0)22 751 31 40
Fax: +48 (0)22 751 58 45
Email: pl@lgcstandards.com

Romania

LGC Standards GmbH
Str. Muncitorilor nr.1/2
400424 Cluj-Napoca
Romania

Tel/Fax: +40 364 116890
Mobil: +40 749 255121
Email: ro@lgcstandards.com

Spain

LGC Standards S.L.U
C/Salvador Espriu 59, 2º planta
08005 Barcelona
España

Tel: +34 93 308 4181
Fax: +34 93 307 3612
Email: es@lgcstandards.com

Sweden

LGC Standards AB
Albanoliden 5, 4 trappor
Box 1737
SE-501 17 Borås
Sweden

Tel: +46 (0)33 20 90 60
Fax: +46 (0)33 20 90 79
Email: se@lgcstandards.com

United Kingdom

LGC Standards
Queens Road
Teddington
Middlesex
TW11 0LY
United Kingdom

Tel: +44 (0)20 8943 8489
Fax: +44 (0)20 8943 8405
Email: uksales@lgcstandards.com

USA

LGC Standards
5853 Post Road, Suite 201B
East Greenwich, RI 02818, USA

Tel: +1 401 398 7892
Fax: +1 401 398 7893
Email: David.Krost@lgcstandards.com



www.lgcstandards.com

Contents

Introduction

About LGC Standards	2
Sourcing of reference materials	2
Producers of reference materials	3
The use of reference materials	4
About proficiency testing	5
Ordering Information	6
Relevant training courses	6

Iron and steel products

Unalloyed steels	7
Alloy steels	11
Highly alloyed steels	16
Pure Iron	21
Special alloys	22
Cast irons	23
Ferro alloys	25
Ores, iron oxide	27
Ceramic materials and minerals	31
Slags	35
Steel with oxygen and nitrogen content	37
Setting up sample for spectrometric analysis of low alloyed steels	38

Non ferrous metals and alloys

Aluminium base standards	38
Copper base standards	40
Oxygen in copper	44
Unalloyed zinc	45
Zinc alloys	45
Lead base standards	47
Gold base standards	47
Nickel base standards	48
Titanium base standards	48
Zirconium base standards	48
Tungsten base standards	48

Bauxite and Fluorspar

	49
--	----

Ores

	50
--	----

Refractories and carbides

Refractories and carbides	52
Cement	54

Special materials

	54
High tech ceramics	54
Multielement glass for XRF analysis - type A - type B	55

Coals

	55
--	----

Introduction

About LGC Standards

LGC Standards is Europe's most comprehensive source of reference materials. We work closely with the leading manufacturers to offer laboratories worldwide improved access to reference materials, all designed to cover an increasingly large range of parameters.

Our network of offices in Europe, India and the United States, combined with our extensive experience in selling reference materials and our technical expertise, allow LGC Standards to work in partnership with our customers to provide fast delivery together with sound technical advice as well as dealing with complex import and export regulations.

LGC, whose Research and Technology Division acts as the designated UK National Measurement Institute (NMI) for Chemical and Biochemical measurement, has a long history in the development and validation of analytical methods and the production of reference materials.

Many of the analytical methods, which cover the food, environment, industrial, clinical and pharmaceutical sectors, are accredited to ISO/IEC 17025 (Requirements for the competency of testing and calibration laboratories). Full details of the accreditation schedule can be found on the UKAS website www.ukas.com.

LGC's Research and Technology Division is accredited to ISO Guide 34 (General Requirements for the competence of reference materials producers) for the production of reference materials. The measurement capabilities used to produce certified reference materials in its capacity as the UK NMI are verified through participation on key comparison studies organised by the Consultative Committee for the Amount of Substance (CCQM) of the International Weights and Measures Organisation (BIPM). The certificates are recognised by other National Measurement Institutes, under the Mutual Recognition Arrangement.

In its role as the NMI, LGC serves on the International Organisation for Standardisation (ISO) Committee for Reference Materials (REMCO), which aims to carry out and encourage a broad international effort for the harmonization, production and application of certified reference materials (CRMs).

LGC has a long history in the development and validation of analytical methods and the reference materials production facility at LGC, complements this expertise. LGC also offers a range of chemical and bio-analytical laboratory services across a variety of industries including: food and agriculture, environment, life sciences, chemicals and forensic science. LGC also plays a pivotal role in a number of programmes with both government and industry to improve analytical standards and laboratory performance.

This catalogue

The LGC Standards' catalogue for Metals and related materials presents a listing of reference materials grouped according to their matrix type and further subdivided by their chemistry. Listings give information on each product such as composition, weight of sample and format. The Reference Materials listed in this catalogue represent some of the most comprehensive ranges suitable for chemical and spectrometric techniques such as XRF/XRD, OES, AA, ICP and Combustion analysis.

Sourcing of reference materials

While this catalogue includes some of the more widely used materials, it is not possible to provide a comprehensive list of all the reference materials which can be provided by LGC Standards. Through our network of offices and extensive experience in sourcing reference materials, the list of materials we can provide is impractical to include in one catalogue.

Please use this catalogue as a first point of reference, and to locate your local LGC Standards office. Contact your nearest office to place an order and to request any products not included in this catalogue.

Producers of reference materials

Federal Institute for Materials Research and Testing (BAM)

The Federal Institute for Materials Research and Testing (BAM) has a long tradition in the production of Certified Reference Materials. Starting in 1912 with a "Normal Steel" for the determination of carbon, the development of new CRMs has increased continuously. One year later eight steel samples with different carbon contents were available. The development continued with the participation of regional German material research and testing institutes as well as industry (1957). In 1968 within the framework of EURONORM, the first European CRMs in the field of iron and steel were issued.

Today a large range of ferrous and non ferrous CRMs together with environmental CRMs and CRMs for engineering materials are offered.

IRMM

BCR® and IRMM reference materials (BCR® is a registered trademark of JRC-IRMM) are the products of both research funding and direct action programmes of the European Commission, in which new or improved measurement or testing methods are developed. These programmes are aimed at improving, harmonising or standardising measurements and testing in the European Union. As an authorised distributor of BCR® reference materials LGC Standards currently holds stock of more than 5000 units of certified BCR® and IRMM reference materials under carefully controlled and monitored conditions.

European Reference Materials (ERM®)

The ERM® range of reference materials was launched in May 2004. It is the result of collaboration between three major reference material producers, LGC in the UK, the Institute for Reference Materials and Measurement (IRMM) in Belgium and Bundesanstalt für Materialforschung und Prüfung (BAM) in Germany. The partners are committed to using the most advanced principles for the production of certified reference materials. The certified values have clearly defined and stated traceability and are internationally recognized through participation of the partners in key comparisons organized by the Bureau International des Poids et Mesures (BIPM). All ERM® materials are subject to rigorous homogeneity and stability testing guaranteeing the certified values for every unit over its entire shelf life.

The National Institute of Standards and Technology (NIST)

NIST produces standard reference materials (SRMs®). Based in the United States, NIST has provided reference materials to industry and commerce for nearly 100 years. NIST collaborates with companies to provide academia and industry with SRMs for expanding areas such as air and water pollution, which are international issues.

Bureau of Analysed Samples (BAS)

Bureau of Analysed Samples has a long history in producing certified reference materials (CRMs) for chemical analysis (BCS-CRMs and EURONORM-CRMs) including a wide range of plain carbon and alloy steels, cast iron, ferro alloys, non-ferrous, alloys, ores and ceramic materials for the iron and steel industry.

National Analysis Centre for Iron and Steel (NACIS)

National Analysis Centre for Iron and Steel is a research and development centre of analysis and test technology for iron and steel. Based in the Central Iron and Steel Research Institute in Beijing, China, NACIS has over fifty years experience in research of a wide variety of reference materials and production of metal alloy, ore, ferrous alloy certified reference materials (CRMs).

EURONORM certified reference materials

EURONORM certified reference materials for the chemical analysis of iron and steel products

EURONORM certified reference materials are prepared under the auspices of the European Committee for Iron and Steel Standardization (ECISS) in a collaboration between the producing organizations in: France: Institut de Recherches de la Sidérurgie (IRSID), Centre de Dévelopement des Industries de Mise en Forme des Matériaux (CTIF), the Federal Republic of Germany: Iron and Steel CRM Working Group comprising Bundesanstalt für Materialforschung und -prüfung (BAM), Max-Planck-Institut für Eisenforschung, Verein Deutscher Eisenhüttenleute (VDEh), the United Kingdom: Bureau of Analysed Samples Limited, Sweden/Finland: Jernkontoret, Swedish Institute for Metals Research. Starting in 1968 EURONORM-CRMs have been analysed by laboratories in the European Community (EC) and further European countries. A number of former national CRMs are also listed which after examination by laboratories in the EC they have been accepted as EURONORM-CRMs.

Approximately 20 laboratories take part in the analysis. Each laboratory is requested to analyse the elements to be determined four times. A statistical evaluation of the laboratory mean values is carried out with respect to their normal distribution and the identification of any outlying values.

The finely divided EURONORM-CRMs are supplied in glass bottles containing 100 g. Some EURONORM-CRMs are also available in solid form (discs). Samples in the form of chips, pins and balls with certified oxygen and nitrogen content are also available.

Types of material

The following types of material are available as EURONORM-CRMs:

Unalloyed steels (0), alloyed steels (1), highly alloyed steels (2), special alloys (3), cast iron (4), ferro-alloys (5), ores (6), ceramics (7) and slags (8).

The system of numbering of the samples allows an easy reference to the type of material. The first digit of the sample number shows the type of material (0 - unalloyed steel, 1 - low alloyed steel, 2 - highly alloyed steel etc.). The second and third digit characterises the single sample. Another digit, separated by a hyphen gives the number of editions of the material.

Content of the Euronorm certificate

On the head of the certificate the EURONORM-number and the type of material is given. The mean values of the laboratories involved in the certification campaign are given in a table together with indicative values.

The mean values of the accepted data sets, their standard deviations and the standard deviations of the laboratories are also given in the table. The sign "-" in the table stands for an outlier pointed out by statistical tests. The certified values are given in a second table together with their uncertainties (95%-level) or standard deviations. Additionally the following information is given: The owner of the material, a characterisation of the sample (e.g. grain size, dimensions of compact samples), the laboratories involved in the certification campaign, the analytical methods used for element determination, sources for additional information published by ECISS/EGKS.

The use of reference materials

Reference materials are instrumental in ensuring the reliability of analytical measurements and so ensuring that decisions are based on reliable data. When choosing a matrix reference material for a particular application the analyst should consider the following factors before selecting a material:

When choosing a matrix reference material for a particular application the analyst should consider the following factors before selecting a material:

- Matrix match and potential interferences
- Analytes
- Measurement range
- Measurement uncertainties
- Certification procedures used by the producer
- Documentation supplied with the material (e.g. certificate or report).

About proficiency testing (PT)

Proficiency testing (PT) is a powerful quality assurance tool for laboratories undertaking analytical measurements. A PT scheme provider distributes test materials on a regular basis to participating laboratories for independent testing. The results are returned to the organiser of the scheme who makes an analysis of the results and provides a report to all the participants.

There are a number of benefits of taking part in a PT scheme:

- Provides laboratories with a mechanism to compare their measurements with others
- Enables laboratories to demonstrate the quality of their results to third parties e.g. customers, regulators and accreditation bodies
- Facilitates the monitoring of trends, over time, in the quality of measurements
- Assists in the evaluation of methods and instrumentation
- Provides support and information to laboratory staff and their customers

The LGC Proficiency Testing Group (incorporating Quality Management Ltd and Aquacheck Ltd) is a major international provider of proficiency testing services. It has over twenty years experience in all aspects of proficiency testing services for laboratories undertaking chemical and microbiological analysis. Schemes of relevance to a wide range of sectors are available and most of them are run on an international basis. All LGC schemes are accredited by the United Kingdom Accreditation Service (UKAS).

Current schemes provided include:

Environmental:

- AQUACHECK - chemical analysis of clean waters, waste waters, sludges, sediments and soils
- QWAS - microbiological assessment of waters, effluents and sludges
- CONTEST - analysis of contaminated soils for a wide range of contaminants

Food:

- QMS - microbiological examination of food and food ingredients
- QDCS - composition and safety testing in the dairy sector
- QMAS - chemical analysis of meat
- QGS - analysis of gelatine samples
- QCS - microbiological testing of chocolate
- QFCS - chemical testing of food products

Beverage:

- BAPS - chemical, microbiological and sensory analysis of a range of beers.
- DAPS - analysis of a wide range of alcoholic beverages (except beer)
- QBS - microbiological and chemical testing of soft drinks
- SODAS - chemical analysis of carbonated and still soft drinks
- MAPS - analysis of malt and barley used by the malting, brewing and distilling industries
- SUPS - analysis of raw sugar used in soft drink production.

Others:

- QMIS - identification of micro organisms
- PHARMASSURE - measurement of a range of analytes in pharmaceutical products used in hospitals
- DIPS - dissolution proficiency scheme
- QUARTZ - toxicological analysis of drugs in post-mortem blood
- TOYTEST - toy safety testing to the European Standard EN71 and American Standard ASTM F963
- PACQS - particle analysis and characterisation quality scheme

In addition to these regular schemes LGC is able to offer customised or closed proficiency testing schemes tailored to a specific organisation's requirements.

LGC schemes are now available either through LGC Standards or through your local distributor.

Please contact your local LGC Standards office to find out more.

Relevant literature

Applications of Reference Materials in Analytical Chemistry VAM
V. Barwick, S. Burke, R. Lawn, P. Roper and R. Walker RSC

Reference Materials in Analytical Chemistry – A Guide for Selection and Use
A. Zchunke published by Springer-Verlag.

Reference Materials for Chemical Analysis – Certification, Availability and Proper Usage
M. Stoeppeler, W.R. Wolf, P.J. Jenks published by Wiley-VCH.

Proficiency Testing in Analytical Chemistry, R.Lawn, M.Thompson and R.Walker published by The Royal Society of Chemistry.

Relevant training courses

A range of training courses is available from LGC to help laboratory managers and analysts demonstrate competence in, and keep abreast of, quality assurance issues and practices. LGC's analytical quality training programme includes:

- Achieving traceability in chemical testing
- Using proficiency testing in the analytical laboratory
- Method validation
- Principles and practice of measurement uncertainty in chemical testing laboratories
- Quality systems in testing laboratories
- Statistics for analytical chemists
- Further statistical tools for analytical chemists
- Evaluating measurement uncertainty for chemical testing laboratories

The majority of the courses are run in Teddington, South West London, UK. In addition, LGC can provide training for groups of staff at your own site, where the courses can be customised to meet your exact needs. For further information, please contact:

Lorraine Didinal
LGC Training Centre
Queens Road
Teddington, Middlesex
TW11 0LY, UK
Tel: +44 (0)20 8943 7631
Fax: +44 (0)20 8943 7314
Email: training@lgc.co.uk
Web: www.lgc.co.uk

General Ordering Information

Prices and delivery procedures are shown in the price list that accompanies this catalogue, or are available from your local LGC Standards sales office. For products requiring special delivery procedures (cooled shipping, dangerous goods, etc) additional charges will be applied.

Please check with your local LGC Standards sales office for detailed procedures and transport charges where applicable.

Unless otherwise agreed in advance and in writing, orders are accepted only against LGC Standards' standard terms and conditions of sale.

Once delivered to the customer, reference substances are not returnable. For this reason it is very important for users to be certain that the product ordered meets their needs.

LGC Standards' technical staff are available to advise on the use and suitability of a particular product. Customers requiring assistance with the use or application of a particular reference substance should contact their local LGC Standards office, contact details are provided at the beginning of this catalogue.

For further information please contact your local office or visit www.lgcstandards.com

Iron and steel products

Unalloyed steels

Code	Product	Unit	
ECRM-D 030-4	Steel type 60-2 - chips	100 g	
	Year of issue: 1973		
	Certified values		
	C 0.456 ± 0.004 %	S 0.021 ± 0.002 %	As 0.012 ± 0.002 %
	Si 0.318 ± 0.007 %	Cr 0.117 ± 0.004 %	Cu 0.061 ± 0.002 %
	Mn 0.603 ± 0.004 %	Ni 0.042 ± 0.002 %	N 0.0051 ± 0.0003 %
	P 0.018 ± 0.002 %	Al _{total} 0.042 ± 0.006 %	Sn 0.0055 ± 0.0007 %
ECRM-D 031-3	Unalloyed steel - chips	100 g	
	Year of issue: 1972		
	Certified values		
	C 0.055 ± 0.002 %	P 0.014 ± 0.001 %	As 0.013 ± 0.002 %
	Si 0.037 ± 0.004 %	S 0.021 ± 0.001 %	Cu 0.020 ± 0.002 %
	Mn 0.329 ± 0.007 %	Al _{total} 0.054 ± 0.002 %	N 0.0050 ± 0.0004 %
ECRM-D 032-2	Steel with 0.27% carbon - 0.56% manganese steel - chips	100 g	
	Year of issue: 1968		
	Certified values		
	C 0.271 ± 0.007 %	S 0.0254 ± 0.0010 %	Cu 0.085 ± 0.002 %
	Si 0.282 ± 0.007 %	Cr (0.088) %	N 0.0044 ± 0.0009 %
	Mn 0.556 ± 0.008 %	Ni (0.040) %	Sn (0.006) %
	P 0.0129 ± 0.0007 %	As 0.020 ± 0.002 %	
	(Values in parenthesis are indicative values)		
ECRM-D 035-2	1.3% Carbon steel - powder	100 g	
	Year of issue: 1998		
	Certified values		
	C 1.277 ± 0.005 %	Cr 0.0104 ± 0.0003 %	As 0.0017 ± 0.0001 %
	Si 0.216 ± 0.004 %	Mo 0.0056 ± 0.0002 %	Cu 0.0085 ± 0.0002 %
	Mn 0.305 ± 0.002 %	Ni 0.0190 ± 0.0004 %	N 0.0230 ± 0.0004 %
	P 0.0038 ± 0.0003 %	Al _{total} 0.0193 ± 0.0005 %	Ti 0.0030 ± 0.0001 %
	S 0.0111 ± 0.0003 %	Al _{acid-sol} 0.0177 ± 0.0004 %	
ECRM-D 035-2-D	1.3% Carbon steel - disc (38 mm x 20 mm)	disc	
ECRM-D 036-1	Tool steel C85 W2, no. 1.1630 - chips	100 g	
	Year of issue: 1968		
	Certified values		
	C 0.858 ± 0.008 %	Cr (0.091) %	N 0.0100 ± 0.0008 %
	Si 0.194 ± 0.005 %	Ni (0.058) %	Sn (0.006) %
	Mn 0.327 ± 0.010 %	Al _{total} (0.015) %	V (0.019) %
	P 0.0074 ± 0.0009 %	As 0.0233 ± 0.0007 %	
	S 0.0095 ± 0.0009 %	Cu 0.065 ± 0.005 %	
	(Values in parenthesis are indicative values)		
ECRM-D 039-2	Free cutting steel with 0.2% lead- chips	100 g	
	Year of issue: 1971		
	Certified values		
	C 0.107 ± 0.003 %	S 0.310 ± 0.005 %	Cu 0.117 ± 0.006 %
	Si 0.011 ± 0.002 %	Cr 0.048 ± 0.003 %	N 0.0113 ± 0.0004 %
	Mn 1.274 ± 0.014 %	Ni 0.051 ± 0.003 %	Pb 0.207 ± 0.005 %
	P 0.083 ± 0.004 %	As 0.018 ± 0.001 %	Sn 0.016 ± 0.001 %
ECRM-D 042-1	0.05% Niobium steel- chips	100 g	
	Year of issue: 1972		
	Certified values		
	C 0.108 ± 0.003 %	S 0.024 ± 0.024 %	Cu 0.041 ± 0.002 %
	Si 0.037 ± 0.005 %	Cr 0.016 ± 0.004 %	N 0.0078 ± 0.0007 %
	Mn 0.666 ± 0.010 %	Ni 0.029 ± 0.002 %	Nb 0.054 ± 0.005 %
	P 0.0057R ± 0.0004 %	Al 0.010 ± 0.001 %	
	R: revised value		
ECRM-D 077-2	Steel X 60 with 0.06% vanadium- chips	100 g	
	Certified values		
	C 0.151 ± 0.004 %	Cr (0.016) %	Cu (0.029) %
	Si 0.293 ± 0.008 %	Mo (0.003) %	N 0.0054 ± 0.0005 %
	Mn 1.28 ± 0.02 %	Ni (0.021) %	Sn (0.003) %
	P 0.022 ± 0.001 %	Al 0.034 ± 0.002 %	V 0.058 ± 0.003 %
	S 0.014 ± 0.001 %	As 0.007 ± 0.001 %	
	(Values in parenthesis are indicative values)		

Iron and steel products

Code	Product	Unit
ECRM-D 079-2	Free cutting steel - chips Certified values C 0.596 ± 0.006 % Cr 0.0382 ± 0.0023 % N 0.0074 ± 0.0005 % Si 0.247 ± 0.006 % Ni 0.0219 ± 0.0010 % Sn 0.0037 ± 0.0008 % Mn 0.743 ± 0.013 % Al 0.0209 ± 0.0017 % Ti (0.0021) % P 0.0234 ± 0.0012 % As 0.0040 ± 0.0007 % S 0.192 ± 0.006 % Cu 0.0462 ± 0.0010 % (Values in parenthesis are indicative values)	100 g
ECRM-D 082-1	Free cutting steel - chips Year of issue: 1976 Certified values C 0.415 ± 0.003 % Cr 0.018 ± 0.001 % N (0.0047) % Si 0.235 ± 0.005 % Ni 0.027 ± 0.001 % Pb 0.149 ± 0.004 % Mn 0.769 ± 0.008 % Al 0.032 ± 0.002 % Te 0.030 ± 0.001 % P 0.013 ± 0.001 % As (0.029) % S 0.030 ± 0.001 % Cu 0.025 ± 0.001 % (Values in parenthesis are indicative values)	100 g
ECRM-D 083-1	Unalloyed steel - chips Year of issue: 1978 Certified values C 0.0262R ± 0.0003 %* Cr (0.0129) % Cu 0.016 ± 0.001 % Mn 0.289 ± 0.004 % Ni 0.014 ± 0.001 % N 0.0022 ± 0.0003 % P 0.0077 ± 0.0009 % Al (0.0044) % S 0.0100 ± 0.0005 % As (0.0043) % (Values in parenthesis are indicative values) * 95%-confidence interval	100 g
ECRM-B 051-1	0.1% Sulfur carbon steel - chips Year of issue: 1970 Certified values C 0.181 % P (0.025) % Ni (0.14) % Si (0.11) % S 0.126 % Cu (0.15) % Mn 1.18 % Cr (0.05) % (Values in parenthesis are indicative values)	100 g
ECRM-B 054-1	0.09% Phosphorus carbon steel - chips Year of issue: 1975 Certified values C 0.22 % P 0.092 % Ni (0.14) % Si (0.05) % S (0.10) % Cu (0.21) % Mn (0.88) % Cr (0.17) % (Values in parenthesis are indicative values)	100 g
ECRM-B 055-1	0.5% Carbon steel - chips Year of issue: 1972 Certified values C 0.51 % P 0.016 % Ni (0.12) % Si 0.24 % S 0.036 % Cu (0.15) % Mn 0.77 % Cr (0.16) % (Values in parenthesis are indicative values)	100 g
ECRM-B 056-1	0.8% Carbon steel - chips Year of issue: 1967 Certified values C 0.79 % P 0.043 % Ni (0.13) % Si 0.30 % S 0.030 % Cu (0.12) % Mn 1.02 % Cr (0.11) % (Values in parenthesis are indicative values)	100 g
ECRM-B 057-2	0.05% Carbon steel (BS/DIN EN 10130) - chips Year of issue: 1999 Certified values C 0.0507 ± 0.0009* % S 0.0127 ± 0.0003* % Al _{acid sol} (0.056) % Si (0.003) % Cr 0.0114 ± 0.0003* % Cu 0.0146 ± 0.0002* % Mn 0.246 ± 0.002* % Ni 0.0096 ± 0.0004* % N 0.0023 ± 0.0001* % P 0.0120 ± 0.0003* % Al _{total} 0.059 ± 0.001* % (Values in parenthesis are indicative values) * 95%-confidence interval	100 g
ECRM-B 057-2-D	0,05% Carbon steel (BS/DIN EN 10130) - disc (35-39 mm x 20-35 mm)	disc

Iron and steel products

Code	Product	Unit
ECRM-B 058-2	0.15% Sulfur steel (BS 1429) - chips Year of issue: 2002 Certified values C 0.424 ± 0.002 % S 0.1712 ± 0.0013* % As 0.0095 ± 0.0005* % Si 0.1080 ± 0.0012* % Cr 0.1211 ± 0.0010* % Cu 0.261 ± 0.002 % Mn 1.186 ± 0.005* % Mo 0.0589 ± 0.0007* % N 0.0107 ± 0.0002* % P 0.0098 ± 0.0003* % Ni 0.199 ± 0.002 % * 95%-confidence interval	100 g
ECRM-B 058-2-D	0,15% Sulfur steel (BS 1429) - disc (35-39 mm x 20-35 mm)	disc
ECRM-B 059-2	0.7% Carbon steel (BS EN 10277) - chips Year of issue: 2002 Certified values C 0.721 ± 0.003* % S 0.0084 ± 0.0003* % Al _{total} 0.00044 ± 0.00008* % Si 0.188 ± 0.002* % Cr 0.0090 ± 0.0003* % Al _{acid sol.} 0.00020 ± 0.00006* % Mn 0.495 ± 0.003* % Mo 0.0018 ± 0.0002* % Cu 0.0074 ± 0.0001* % P 0.0046 ± 0.0002* % Ni 0.0198 ± 0.0006* % N 0.0051 ± 0.0002* % * 95%-confidence interval	100 g
ECRM-B 059-2-D	0.7% Carbon steel (BS EN 10277) - disc (35-39 mm x 20-35 mm)	disc
ECRM-B 060-1	0.1% Carbon steel - chips Year of issue: 1975 Certified values C 0.122 % S (0.031) % Cu (0.060) % Si (0.17) % Cr (0.028) % N 0.0040 % Mn 0.45 % Ni (0.039) % P (0.024) % Al _{total} (0.004) % (Values in parenthesis are indicative values)	100 g
ECRM-B 061-1	0.2% Carbon steel - chips Year of issue: 1968 Certified values C 0.210 % P (0.019) % Cu (0.10) % Si 0.12 % S (0.034) % Mn (0.61) % Ni (0.21) % (Values in parenthesis are indicative values)	100 g
ECRM-B 063-1	1.2% Carbon steel - chips Year of issue: 1972 Certified values C 1.26 % P 0.019 % Ni (0.10) % Si 0.24 % S 0.022 % Cu (0.09) % Mn 0.30 % Cr (0.16) % N (0.005) % (Values in parenthesis are indicative values)	100 g
ECRM-B 064-1	Nb/Ti Interstitial free steel (EN 10130 and EN 10142) - chips Year of issue: 2002 Certified values C 0.0026 ± 0.0002 %* Ni 0.0115 ± 0.0003 %* Nb 0.0146 ± 0.0003 %* Si 0.0065 ± 0.0006 %* Al _{total} 0.0330 ± 0.0006 %* Pb 0.00018 ± 0.00002 %* Mn 0.1641 ± 0.0011 %* Al _{acid sol.} 0.0302 ± 0.0004 %* Sn 0.00051 ± 0.00005 %* P 0.0091 ± 0.0003 %* As 0.0036 ± 0.0002 %* Ti 0.0189 ± 0.0003 %* S 0.0104 ± 0.0003 %* Co 0.0027 ± 0.0001 %* V 0.00015 ± 0.00006 %* Cr 0.0184 ± 0.0004 %* Cu 0.0077 ± 0.0002 %* W (0.00022) % Mo 0.00077 ± 0.00007 %* N 0.0026 ± 0.0001 %* (Values in parenthesis are indicative values) * 95%-confidence interval	100 g
ECRM-B 064-1-D	Nb/Ti Interstitial free steel (EN 10130 and EN 10142) - disc (35-39 mm x 20-35 mm)	disc
ECRM-B 084-1	Unalloyed steel - chips Year of issue: 1979 Certified values C 0.391 ± 0.005 % P 0.018 ± 0.001 % Ni 0.154 ± 0.004 % Si 0.265 ± 0.004 % S 0.029 ± 0.001 % Cu 0.267 ± 0.004 % Mn 0.860 ± 0.011 % Mo 0.033 ± 0.002 % Sn 0.023 ± 0.003 % (Values in parenthesis are indicative values)	100 g
ECRM-B 085-1	Unalloyed steel - chips Year of issue: 1977 Certified values C 0.067 ± 0.002 % S 0.336 ± 0.008 % V 0.0021 ± 0.0007 % Si 0.008 ± 0.001 % Co 0.019 ± 0.001 % Sb 0.0073 ± 0.0005 % Mn 0.977 ± 0.015 % Cu 0.291 ± 0.008 % Zn 0.0025 ± 0.0006 % P 0.062 ± 0.002 % Pb 0.0010 ± 0.0001 %	100 g
ECRM-B 085-1-D	Unalloyed steel - disc (35-39 mm x 20-35 mm)	disc

Iron and steel products

Code	Product	Unit
ECRM-B 086-1	Unalloyed steel - chips Year of issue: 1978 Certified values C 0.297 ± 0.005 % S 0.037 ± 0.001 % Cu 0.320 ± 0.007 % Si 0.206 ± 0.006 % Cr 0.150 ± 0.006 % Sn 0.026 ± 0.002 % Mn 0.879 ± 0.010 % Ni 0.168 ± 0.003 % P 0.024 ± 0.001 % As 0.023 ± 0.002 %	100 g
ECRM-B 086-1-D	Unalloyed steel - disc (35-39 mm x 20-35 mm)	disc
ECRM-B 087-1	Unalloyed steel - chips Year of issue: 1980 Certified values C 0.174 ± 0.004 % Cr 0.078 ± 0.003 % Cu 0.171 ± 0.004 % Si 0.263 ± 0.004 % Mo 0.021 ± 0.002 % Sn 0.017 ± 0.001 % Mn 0.671 ± 0.006 % Ni 0.118 ± 0.003 % Sb 0.0046 ± 0.0004 % P 0.010 ± 0.001 % As 0.024 ± 0.002 % S 0.046 ± 0.002 % Co 0.015 ± 0.001 %	100 g
ECRM-B 090-1	1% Carbon steel - chips Year of issue: 1983 Certified values C 1.054 ± 0.005 %* Ni 0.053 ± 0.002 %* Ga 0.00228 ± 0.00020 %* Si 0.281 ± 0.003 %* N 0.0146 ± 0.0002 %* Hg (< 0.000010) % Mn 0.226 ± 0.002 %* Nb 0.00043 ± 0.00004 %* Sb 0.00090 ± 0.00008 %* P 0.0128 ± 0.0003 %* Pb 0.00239 ± 0.00006 %* Se (£ 0.0002) % S 0.0095 ± 0.0003 %* V 0.204 ± 0.006 %* Te < 0.0002 % Cr 0.121 ± 0.003 %* Bi < 0.00002 % Tl < 0.0001 % Mo 0.0089 ± 0.0004 %* Cd < 0.00002 % Zn 0.00209 ± 0.00009 %* (Values in parenthesis are indicative values) * 95%-confidence interval	100 g
ECRM-B 090-1-D	1% Carbon steel - disc (35-39 mm x 20-35 mm)	disc
ECRM-B 091-1	0.5% Carbon steel - chips Year of issue: 1983 Certified values C 0.518 ± 0.004 % Mo 0.098 ± 0.003 % N 0.0111 ± 0.0005 % Cr 0.312 ± 0.009 % Ni 0.310 ± 0.007 %	100 g
ECRM-B 096-1	Low sulfur low calcium steel - chips Year of issue: 1984 Certified values C 0.113 ± 0.002 % Cr 0.019 ± 0.001 % Cu 0.022 ± 0.001 % Si 0.263 ± 0.006 % Mo 0.003 ± 0.001 % Nb 0.029 ± 0.002 % Mn 1.35 ± 0.01 % Ni 0.020 ± 0.002 % V 0.0036 ± 0.0004 % P 0.019 ± 0.001 % Al _{total} 0.044 ± 0.002 % Ca 0.0021 ± 0.0003 % S 0.0009 ± 0.0002 % Al _{acid sol.} 0.040 ± 0.002 %	100 g
ECRM-B 096-2	Low sulfur low calcium steel (BS EN 10025) - chips Year of issue: 1999 Certified values C 0.105 ± 0.0007 %* S 0.0016 ± 0.0001 %* Al _{total} 0.0460 ± 0.0005 %* Si 0.262 ± 0.002 %* Cr 0.0243 ± 0.0004 %* Cu 0.0170 ± 0.0002 %* Mn 1.320 ± 0.005 %* Mo 0.0020 ± 0.0002 %* Nb 0.0252 ± 0.0005 %* P 0.0128 ± 0.0003 %* Ni 0.0253 ± 0.0004 %* Ca 0.0020 ± 0.0001 %*	100 g
ECRM-B 096-2-D	Low sulfur low calcium steel (BS EN 10025) - disc (35-39 mm x 20-35 mm)	disc
ECRM-F 012-1	Unalloyed steel - chips Year of issue: 1958 Certified values C 0.082 % P 0.083 % S 0.255 %	100 g
ECRM-F 013-1	Unalloyed steel - chips Year of issue: 1959 Certified values P 0.053 S 0.032	100 g
ECRM-F 017-1	Unalloyed steel - chips Year of issue: 1968 Certified values C 0.261 % P 0.0132R % Ni 0.0852 % Si 0.266 % S 0.0225 % Cu 0.0624 % Mn 0.7255 % Cr 0.0444 % N 0.0091 % R: revised value	100 g

Iron and steel products

Code	Product	Unit
ECRM-F 021-1	Unalloyed steel - chips Year of issue: 1959 Certified values C 0.243 % P 0.0121R % Cu 0.167 % Si 0.271 % Cr 0.125 % Mn 1.27 % Ni 0.255 % R: revised value	100 g
ECRM-F 022-1	Unalloyed steel - chips Year of issue: 1964 Certified values C 0.115 % P 0.057 % Mn 0.797 % S 0.300 %	100 g
ECRM-F 023-1	Unalloyed steel with 0.3% lead - chips Year of issue: 1964 Certified values C 0.331 % Mn 0.667 % S 0.0156 % Si 0.264 % P 0.021 % Pb 0.280 %	100 g
ECRM-F 024-1	Unalloyed steel with 0.3% lead - chips Year of issue: 1964 Certified values C 0.104 % Mn 0.726 % S 0.235 % Si 0.139 % P 0.0155R % Pb 0.287 % R: revised value	100 g
ECRM-F 080-1	Unalloyed steel - chips Year of issue: 1973 Certified values C 0.452 ± 0.003 % P 0.028 ± 0.001 % N 0.0073 ± 0.0005 % Si 0.317 ± 0.006 % S 0.024 ± 0.001 % Mn 1.116 ± 0.012 % Cr 0.025 ± 0.002 %	100 g
ECRM-F 081-1	Unalloyed steel - chips Year of issue: 1974 Certified values C 0.099 ± 0.002 % P 0.0129R± 0.0006 % Al _{total} 0.023 ± 0.001 % Si 0.105 ± 0.003 % S 0.014 ± 0.001 % Co 0.017 ± 0.004 % Mn 0.605 ± 0.006 % Ni 0.042 ± 0.002 % Cu 0.026 ± 0.001 % R: revised value	100 g

Alloy steels

ECRM-D 126-1	Tool steel 90MnCrV8, No. 1.2842 - chips Year of issue: 1963 Certified values C 0.841 ± 0.008 % S 0.0050 ± 0.0007 % V 0.143 ± 0.004 % Si (0.241) % Cr 0.317 ± 0.009 % Te (0.0002) % Mn 1.817 ± 0.009 % Ni (0.038) % P 0.0092 ± 0.0011 % Cu (0.098) % (Values in parenthesis are indicative values)	100 g
ECRM-D 128-1	Structural steel with 0.9% Ti, No. 1.5310 - chips Year of issue: 1972 Certified values C 0.085 ± 0.003 % S 0.007 ± 0.001 % Cu 0.055 ± 0.003 % Si 0.949 ± 0.010 % Cr 0.108 ± 0.003 % N (0.0024) % Mn 0.839 ± 0.010 % Ni 0.046 ± 0.006 % Ti 0.890 ± 0.013 % P 0.007 ± 0.001 % Al 0.286 ± 0.010 % V (0.008) % (Values in parenthesis are indicative values)	100 g
ECRM-D 130-1	Rail steel - chips Year of issue: 1968 Certified values C 0.546 ± 0.005 % Cr (0.032) % Cu 0.072 ± 0.003 % Si 0.313 ± 0.006 % Ni (0.031) % N 0.0093 ± 0.0008 % Mn 1.593 ± 0.009 % Al 0.0037 ± 0.0005 % Sn (0.006) % P 0.0209 ± 0.0017 % Al _{acid sol.} 0.0019 ± 0.0006 % V (0.003) % S 0.0158 ± 0.0011 % As 0.0167 ± 0.0011 % (Values in parenthesis are indicative values)	100 g

Iron and steel products

Code	Product	Unit
ECRM-D 179-2	Tool steel 60 WCrV 7, No. 1.2550 - chips Year of issue:1990 Certified values C 0.598 ± 0.009 % Cu 0.111 ± 0.004 % Ga 0.00129 ± 0.00012 % Si 0.579 ± 0.011 % N 0.0068 ± 0.0005 % Hg (< 0.00001) % Mn 0.539 ± 0.010 % Nb 0.00144 ± 0.00013 % Sb 0.00175 ± 0.00010 % P 0.0267 ± 0.0024 % Pb 0.00013 ± 0.00002 % Se (< 0.00020) % S (0.0006) % Ti (0.0014) % Te < 0.00020 % Cr 1.08 ± 0.03 % V 0.188 ± 0.007 % Tl (< 0.000035) % Mo 0.070 ± 0.006 % W 1.87 ± 0.05 % Zn 0.00023 ± 0.00004 % Ni 0.078 ± 0.007 % Bi < 0.00003 % Co (0.015) % Cd < 0.00003 % (Values in parenthesis are indicative values)	100 g
ECRM-D 179-2-D	Tool steel 60 WCrV 7, No. 1.2550 - disc (35-39 mm x 20-35 mm)	disc
ECRM-D 180-1	Low alloyed steel - chips Year of issue:1973 Certified values C 0.197 ± 0.005 % S 0.0249 ± 0.0010 % Cu 0.115 ± 0.004 % Si 0.362 ± 0.007 % Cr 1.250 ± 0.018 % N 0.0068 ± 0.0009 % Mn 1.286 ± 0.015 % Ni 0.096 ± 0.008 % P 0.0174 ± 0.0010 % As 0.030 ± 0.002 %	100 g
ECRM-D 181-1	Low alloyed steel - chips Year of issue:1973 Certified values C 0.590 ± 0.005 % Cr 0.126 ± 0.004 % N 0.0068 ± 0.0005 % Si 1.054 ± 0.015 % Ni 0.070 ± 0.004 % Sn (0.015) % Mn 1.047 ± 0.008 % Al 0.022 ± 0.004 % Sb (0.004) % P 0.018 ± 0.001 % As (0.026) % S 0.035 ± 0.001 % Cu 0.174 ± 0.005 % (Values in parenthesis are indicative values)	100 g
ECRM-D 182-1	Tool steel 80 CrV 2, No. 1.2235 - chips Year of issue:1974 Certified values C 0.790 ± 0.008 % Ni 0.152 ± 0.005 % Sn (0.0135) % Si 0.368 ± 0.014 % Al 0.020 ± 0.003 % V 0.177 ± 0.010 % Mn 0.389 ± 0.007 % As (0.0202) % Mg (0.0005) % P 0.0076R ± 0.0005 % Cu 0.141 ± 0.004 % Sb 0.0042 ± 0.0005 % S 0.011 ± 0.001 % N 0.0102 ± 0.0004 % Zn 0.0015 ± 0.0002 % Cr 0.591 ± 0.012 % Pb 0.0039 ± 0.0003 % (Values in parenthesis are indicative values) R: revised value	100 g
ECRM-D 183-1	High Copper steel, No. 1.8962 - chips Year of issue:1973 Certified values C 0.083 ± 0.002 % S 0.031 ± 0.001 % As (0.013) % Si 0.421 ± 0.006 % Cr 0.670 ± 0.013 % Cu 0.445 ± 0.010 % Mn 0.354 ± 0.004 % Ni 0.073 ± 0.004 % N 0.0064 ± 0.0006 % P 0.089 ± 0.002 % Al 0.027 ± 0.002 % (Values in parenthesis are indicative values)	100 g
ECRM-D 184-1	Low alloyed steel 35 NiCrMoV 14 6 - chips Year of issue:1978 Certified values C 0.333 ± 0.003 % Mo 0.457 ± 0.009 % N 0.0051 ± 0.0004 % Si 0.218 ± 0.005 % Ni 3.318 ± 0.015 % Sn 0.0044 ± 0.0004 % Mn 0.528 ± 0.006 % Al 0.0052 ± 0.0007 % V 0.108 ± 0.006 % P 0.0047R ± 0.0003 % As 0.0180 ± 0.0011 % Sb (0.0015) % S 0.0032 ± 0.0003 % Co 0.0560 ± 0.0019 % Cr 1.287 ± 0.011 % Cu 0.060 ± 0.002 % (Values in parenthesis are indicative values) R: revised value	100 g
ECRM-D 187-1	Low alloyed boron steel - chips Year of issue:1982 Certified values C 0.195 ± 0.003 % Cr 1.186 ± 0.015 % B 0.0004 ± 0.0002 % Si 0.026 ± 0.002 % Mo 0.035 ± 0.002 % Co 0.014 ± 0.001 % Mn 1.354 ± 0.011 % Ni 0.096 ± 0.003 % Cu 0.161 ± 0.003 % P 0.014 ± 0.001 % Al 0.046 ± 0.002 % N 0.014 ± 0.001 % S 0.025 ± 0.001 % As 0.018 ± 0.002 % Sn 0.011 ± 0.001 %	100 g

Code	Product	Unit
ECRM-D 191-1	3% Si-Stahl - disc (35-39 mm x 20-35 mm) Year of issue:1986 Certified values C 0.013 ± 0.002 % Cr 0.025 ± 0.002 % Cu 0.0080 ± 0.0006 % Si 3.140 ± 0.022 % Mo (0.0021) % N 0.0026 ± 0.0003 % Mn 0.025 ± 0.002 % Ni 0.018 ± 0.002 % Ti 0.009 ± 0.002 % P 0.011 ± 0.001 % Al 0.397 ± 0.015 % S 0.0017 ± 0.0003 % As 0.0031 ± 0.0006 % (Values in parenthesis are indicative values)	disc
ECRM-D 191-2	Dynamo steel, material number: 1.0812 - chips Year of issue:2006 Certified values C 0.0043 ± 0.0002 %* Cr 0.0314 ± 0.0006 %* Cu 0.0165 ± 0.0003 %* Si 3.267 ± 0.012 %* Mo 0.0020 ± 0.0002 %* N 0.00105 ± 0.00009 %* Mn 0.1334 ± 0.0019 %* Ni 0.0224 ± 0.0004 %* Sn 0.0050 ± 0.0005 %* P 0.0087 ± 0.0004 %* Al 0.985 ± 0.006 %* Ti 0.0024 ± 0.0002 %* S 0.0029 ± 0.0002 %* As 0.0018 ± 0.0003 %* Sb (0.0007) %	100 g
	* 95%-confidence interval	
ECRM-D 192-1	Reactor steel, No. 1.6310 - chips Year of issue:1995 Certified values C 0.1875 ± 0.0009 %* Mo 0.482 ± 0.004 %* Co 0.0055 ± 0.0002 %* Si 0.219 ± 0.004 %* Ni 0.755 ± 0.004 %* Cu 0.0453 ± 0.0008 %* Mn 1.377 ± 0.006 %* Al 0.0308 ± 0.0008 %* N 0.0118 ± 0.0002 %* P 0.0029 ± 0.0002 %* Al _{acid sol.} 0.0285 ± 0.0008 %* Sn (0.0030) % S 0.0010 ± 0.0001 %* As (0.003) % Ti (0.0009) % Cr 0.0717 ± 0.0018 %* B (0.00016) % V (0.003) % (Values in parenthesis are indicative values)	100 g
	* 95%-confidence interval	
ECRM-D 192-1-D	Reactor steel, No. 1.6310 - disc (35-39 mm x 20-35 mm)	disc
ECRM-D 193-1	Reactor steel - chips Year of issue:1990 Certified values C 0.139 ± 0.004 % Mo 0.347 ± 0.011 % Cu 0.598 ± 0.009 % Si 0.404 ± 0.006 % Ni 1.178 ± 0.019 % N 0.0108 ± 0.0004 % Mn 0.972 ± 0.017 % Al 0.0257 ± 0.0015 % Nb 0.0232 ± 0.0019 % P 0.0063 ± 0.0006 % As 0.0062 ± 0.0007 % Pb (0.0002) % S 0.0086 ± 0.0006 % B (0.0002) % Ti (0.0013) % Cr 0.182 ± 0.006 % Co 0.0073 ± 0.0007 % V (0.0019) % (Values in parenthesis are indicative values)	100 g
ECRM-D 193-1-D	Reactor steel - disc (40 mm x 35 mm)	disc
ECRM-D 194-1	Low sulfur steel - chips Year of issue:1993 Certified values C 0.1532 ± 0.0011 %* Cr 0.733 ± 0.006 %* B 0.0020 ± 0.0002 %* Si 0.431 ± 0.004 %* Mo 0.2857 ± 0.0026 %* Cu 0.0751 ± 0.0011 %* Mn 1.188 ± 0.004 %* Ni 0.3417 ± 0.0027 %* N 0.0115 ± 0.0002 %* P 0.0097 ± 0.0006 %* Al 0.0837 ± 0.0020 %* V 0.0243 ± 0.0009 %* S 0.00059R ± 0.00005 %* As 0.0042 ± 0.0004 %* Ca 0.0026 ± 0.0002 %*	100 g
	R: revised value * 95%-confidence interval	
ECRM-D 194-1-D	Low sulfur steel - disc (35 mm x 35 mm)	disc
ECRM-B 151-1	Low-carbon high silicon steel - chips Year of issue: 1966 Certified values C 0.028 % Mn 0.085 % S (0.025) % Si 3.49 % P 0.015 %	100 g
	(Values in parenthesis are indicative values)	
ECRM-B 152-1	Manganese-molybdenum steel - chips Year of issue: 1968 Certified values C 0.39 % S 0.043 % Cu 0.21 % Si 0.18 % Cr 0.09 % V (<0.01) % Mn 1.61 % Mo 0.26 % P 0.032 % Ni 0.15 %	100 g
	(Values in parenthesis are indicative values)	

Iron and steel products

Code	Product	Unit
ECRM-B 153-1	Nickel-chromium-molybdenum steel - chips Year of issue: 1977 Certified values C 0.315 % S 0.027 % Al (0.003) % Si 0.079 % Cr 0.66 % Cu 0.088 % Mn 0.81 % Mo 0.58 % Sn 0.011 % P 0.011 % Ni 2.55 % (Values in parenthesis are indicative values)	100 g
ECRM-B 154-1	3% Nickel steel - chips Year of issue: 1964 Certified values C (0.31) % P (0.020) % Mo (0.029) % Si (0.22) % S (0.009) % Ni 3.54 % Mn (0.62) % Cr (0.05) % Cu (0.14) % (Values in parenthesis are indicative values)	100 g
ECRM-B 155-1	Nickel-chromium-molybdenum steel - chips Year of issue: 1963 Certified values C 0.40 % S 0.012 % As (0.04) % Si 0.23 % Cr 1.08 % Cu (0.17) % Mn 0.56 % Mo 0.34 % V (<0.01) % P 0.019 % Ni 1.43 % (Values in parenthesis are indicative values)	100 g
ECRM-B 186-1	Alloy steel - chips Year of issue: 1981 Certified values C 0.610 ± 0.004 % S 0.035 ± 0.002 % Al 0.014 ± 0.001 % Si 1.72 ± 0.02 % Cr 0.218 ± 0.010 % Cu 0.281 ± 0.009 % Mn 0.870 ± 0.008 % Mo 0.048 ± 0.003 % P 0.022 ± 0.001 % Ni 0.190 ± 0.003 % (Values in parenthesis are indicative values)	100 g
ECRM-B 186-1-D	Alloy steel - disc (35-39 mm x 20-35 mm)	disc
ECRM-B 195-1	Chromium-molybdenum-nickel steel - chips Year of issue: 1992 Certified values C 0.756 ± 0.003 % Mo 0.768 ± 0.004 %* V 0.312 ± 0.004 %* Si 0.466 ± 0.004 %* Ni 0.327 ± 0.004 %* Ca 0.0017 ± 0.0001 %* Mn 0.571 ± 0.003 %* Cu 0.0355 ± 0.0005 %* Sb (0.0008) % P 0.0160 ± 0.0006 %* N 0.0100 ± 0.0002 %* Te (< 0.0002) % S 0.0121 ± 0.0003 %* Pb 0.0010 ± 0.0001 %* Zn 0.0046 ± 0.0002 %* Cr 1.566 ± 0.009 %* Sn (0.002) % (Values in parenthesis are indicative values) * 95%-confidence interval	100 g
ECRM-B 195-1-D	Chromium-molybdenum-nickel steel - disc (35-39 mm x 20-35 mm)	disc
ECRM-F 102-1	Low alloyed steel - chips Year of issue: 1958 Certified values C 0.389 % P 0.012 % Mo 1.20 % Si 0.281 % S (0.006) % Ni 4.40 % Mn 0.367 % Cr 0.261 % Cu 0.169 % (Values in parenthesis are indicative values)	100 g
ECRM-F 106-2	Low alloyed steel - chips Year of issue: 1968 Certified values C 0.153 % Cr 1.030 % Cu 0.078 % Si 0.238 % Mo 0.054 % Mn 0.727 % Ni 1.430 %	100 g
ECRM-F 107-1	Low alloyed steel - chips Year of issue: 1960 Certified values C 0.407 % Cr 1.443 % Al 0.345 % Si 0.286 % Mo 0.323 % Cu 0.108 % Mn 0.611 % Ni 0.191 %	100 g

Iron and steel products

Code	Product	Unit
ECRM-F 108-1	Low alloyed steel - chips Year of issue: 1965 Certified values C 0.384 % P 0.0176R % Mo 0.538 % Si 0.340 % S 0.017 % Ni 0.215 % Mn 0.690 % Cr 2.92 % Cu 0.118 % R: revised value	100 g
ECRM-F 112-1	Low alloyed steel - chips Year of issue: 1963 Certified values C 0.348 % Cr 4.78 % Cu 0.115 % Si 1.00 % Mo 1.21 % V 0.604 % Mn 0.191 % Ni 0.234 % W 1.78 %	100 g
ECRM-F 113-1	Low alloyed steel - chips Year of issue: 1963 Certified values C 0.680 % P (0.013) % Mo 0.413 % Si 0.249 % S (0.006) % Ni 1.744 % Mn 0.586 % Cr 0.816 % Cu 0.146 % (Values in parenthesis are indicative values)	100 g
ECRM-F 114-1	Low alloyed steel - chips Year of issue: 1964 Certified values C 0.044 % Mn 0.0655 % S 0.0037 % Si 4.00 % P 0.030 %	100 g
ECRM-F 185-1	Low alloyed steel - chips Year of issue: 1981 Certified values C 1.172 ± 0.005 % P (0.02) % Ni (0.17) % Si (0.25) % S 0.0209 ± 0.0010 % Cu (0.17) % Mn (0.41) % Cr 0.990 ± 0.011 % Sn 0.0212 ± 0.0013 % (Values in parenthesis are indicative values)	100 g
ECRM-F 188-1	Low alloyed steel - chips Year of issue: 1981 Certified values C 1.094 ± 0.008 % Bi < 0.00002 % Se (< 0.0002) % S 0.0113 ± 0.0006 % Cd < 0.00005 % Te < 0.0002 % Cr 1.538 ± 0.016 % Ga 0.00251 ± 0.00017 %* Tl (< 0.000020) % Nb 0.00013 ± 0.00003 %* Hg (< 0.00001) % Zn (< 0.00031) % Pb < 0.0001 % Sb 0.00048 ± 0.00005 %*	100 g
	(Values in parenthesis are indicative values) * 95%-confidence interval	
ECRM-F 190-1	Low alloyed steel - chips Year of issue: 1986 Certified values C 0.395 ± 0.004 % S 0.0044 ± 0.0003 % Co 0.034 ± 0.001 % Si 0.278 ± 0.015 % Cr 2.18 ± 0.03 % N 0.0096 ± 0.0005 % Mn 1.28 ± 0.02 % Mo 0.410 ± 0.007 % P 0.0112 ± 0.0002 % Ni 0.934 ± 0.017 %	100 g
ECRM-F 190-1-D	Low alloyed steel - disc (35-39 mm x 20-35 mm)	disc
ECRM-S 196-1	2% Silicium steel - disc (38 mm x 25 mm) Year of issue: 2000 Certified values C 0.0039 ± 0.0003 %* P 0.0076 ± 0.0005 %* N 0.0020 ± 0.0001 %* Si 1.908 ± 0.014 %* S 0.0005 ± 0.0001 %* Ca 0.0005 ± 0.0001 %* Mn 0.365 ± 0.003 %* Al 0.201 ± 0.002 %*	disc
	* 95%-confidence interval	

Iron and steel products

Code	Product	Unit
ECRM-S 197-1	Low-alloyed steel, No. 1.7321 Year of issue: 2007 Certified values C 0.219 ± 0.002 %* Ni 0.148 ± 0.003 %* Pb (0.0003) % Si 0.275 ± 0.005 %* Al 0.0313 ± 0.0014 %* Sn 0.0097 ± 0.0005 %* Mn 0.792 ± 0.006 %* Al _{acid sol.} (0.0270) % Ti 0.0005 ± 0.0001 %* P 0.0073 ± 0.0003 %* As 0.0083 ± 0.0004 %* V (0.0050) % S 0.0232 ± 0.0007 %* Co 0.0135 ± 0.0003 %* Bi (0.00001) % Cr 0.451 ± 0.003 %* Cu 0.152 ± 0.002 %* Sb (0.0018) % Mo 0.402 ± 0.004 %* N 0.0114 ± 0.0003 %* (Values in parenthesis are indicative values) * 95%-confidence interval	100 g
ECRM-S 197-1-D	Low-alloyed steel, No. 1.7321 - disc (38 mm x 25 mm)	disc
Highly alloyed steels		
ECRM-D 226-1	Stainless steel, 13.7% chromium no. 1.4034 - chips Year of issue: 1967 Certified values C 0.416 ± 0.007 % Cr 13.67 ± 0.06 % N 0.0362 ± 0.0017 % Si 0.514 ± 0.007 % Mo 0.024 ± 0.006 % Sn (0.0068) % Mn 0.434 ± 0.013 % Ni 0.139 ± 0.014 % V 0.022 ± 0.003 % P 0.0207 ± 0.0012 % As (0.0256) % S 0.0094 ± 0.0014 % Co (0.0246) % (Values in parenthesis are indicative values)	100 g
ECRM-D 227-1	High-speed steel, 4% chromium, 2,6% molybdenum, 2,4% vanadium und 3% tungsten, no. 1.3333 - chips Year of issue: 1971 Certified values C 0.950 ± 0.013 % Cr 4.25 ± 0.02 % Sn 0.0251 ± 0.0024 % Si 0.272 ± 0.013 % Mo 2.64 ± 0.05 % V 2.44 ± 0.03 % Mn 0.236 ± 0.007 % Ni 0.114 ± 0.008 % W 3.03 ± 0.06 % P 0.016 ± 0.001 % Cu 0.124 ± 0.005 % Ag (0.000064) % S 0.022 ± 0.002 % N 0.040 ± 0.002 % Sb 0.0035 ± 0.0005 % (Values in parenthesis are indicative values)	100 g
ECRM-D 231-2	Stainless steel, 18% chromium and 10% nickel, no. 1.4306 - chips Year of issue: 2002 Certified values C 0.0140 ± 0.0003 %* Ni 10.105 ± 0.021 %* Pb (0.00007) % Si 0.368 ± 0.006 %* Al 0.0032 ± 0.0004 %* Sn 0.0043 ± 0.0003 %* Mn 1.263 ± 0.009 %* As 0.0048 ± 0.0003 %* Ti 0.0007 ± 0.0002 %* P 0.0179 ± 0.0007 %* B 0.0020 ± 0.0002 %* V 0.0708 ± 0.0008 %* S 0.0250 ± 0.0007 %* Co 0.0402 ± 0.0011 %* W 0.0141 ± 0.0010 %* Cr 18.071 ± 0.018 %* Cu 0.0941 ± 0.0009 %* Sb 0.0011 ± 0.0001 %* Mo 0.301 ± 0.004 %* N 0.0444 ± 0.0004 %* Ca 0.00074 ± 0.00014 %* (Values in parenthesis are indicative values) * 95%-confidence interval	100 g
ECRM-D 235-1	Austenitic manganese steel, 12.7% manganese, no. 1.3401 - chips Year of issue: 1972 Certified values C 0.912 ± 0.014 % S 0.0072 ± 0.0007 % Cu 0.073 ± 0.002 % Si 0.094 ± 0.010 % Cr 0.354 ± 0.014 % N 0.020 ± 0.0008 % Mn 12.73 ± 0.07 % Mo 0.032 ± 0.003 % V (0.012) % P 0.045 ± 0.002 % Ni (0.08) % (Values in parenthesis are indicative values)	100 g
ECRM-D 237-1	Stainless steel, 17% chromium and 10% nickel, no. 1.4550 - chips Year of issue: 1973 Certified values C 0.068 ± 0.002 % Cr 17.24 ± 0.04 % N 0.035 ± 0.002 % Si 0.482 ± 0.013 % Mo 0.306 ± 0.006 % Nb 0.660 ± 0.023 % Mn 1.443 ± 0.018 % Ni 10.32 ± 0.04 % V 0.057 ± 0.005 % P 0.032 ± 0.002 % Co 0.221 ± 0.006 % S 0.012 ± 0.001 % Cu 0.123 ± 0.005 %	100 g

Iron and steel products

Code	Product	Unit
ECRM-D 271-1	High-alloyed steel (toolsteel 1.2344) - chips Year of issue: 2006 Certified values C 0.3698 ± 0.0021 %* Si 0.923 ± 0.006 %* Mn 0.437 ± 0.004 %* P 0.0120 ± 0.0004 %* S 0.00045 ± 0.00008 %* Cr 5.002 ± 0.019 %* Mo 1.247 ± 0.006 %* Ni 0.1552 ± 0.0020 %* Al 0.0234 ± 0.0011 %* As 0.0057 ± 0.0004 %* B (0.0003) % Co 0.0139 ± 0.0005 %* Cu 0.1371 ± 0.0015 %* N 0.0137 ± 0.0003 %* Nb (0.0009) % Pb (0.0005) % Sn 0.0084 ± 0.0002 %* Ti 0.0020 ± 0.0002 %* (Values in parenthesis are indicative values)	100 g
	* 95%-confidence interval	
ECRM-D 271-1-D	High-alloyed steel (toolsteel 1.2344) - disc (37 mm x 25 mm)	disc
ECRM-D 278-1	Stainless steel, 18% chromium, 1% molybdenum and 0,9% carbon, no. 1.4112 - chips Year of issue: 1973 Certified values C 0.903 ± 0.019 % Si 0.336 ± 0.008 % Mn 0.405 ± 0.006 % P 0.0154 ± 0.0014 % S 0.0052 ± 0.0011 % Cr 18.11 ± 0.08 % Mo 1.040 ± 0.030 % Ni 0.236 ± 0.024 % Cu 0.077 ± 0.008 % V 0.077 ± 0.008 %	100 g
ECRM-D 283-1	Highly alloyed steel - chips Year of issue: 1985 Certified values C 1.219 ± 0.009 % Si 0.345 ± 0.017 % Mn 0.217 ± 0.010 % P 0.022 ± 0.002 % S 0.029 ± 0.002 % Cr 4.15 ± 0.06 % Mo 3.41 ± 0.09 % Al 0.0099 ± 0.0014 % As (0.0096) % B 0.0003 ± 0.0001 % Co 10.27 ± 0.17 % N 0.033 ± 0.002 % Pb (< 0.0005) % Sn (0.0065) % V 3.28 ± 0.03 % W 9.66 ± 0.12 %	100 g
	(Values in parenthesis are indicative values)	
ECRM-D 284-2	Stainless steel, 17% chromium, 11% nickel and 2% molybdenum, no. 1.4571 - chips Year of issue: 2000 Certified values C 0.0201 ± 0.0005 %* Si 0.537 ± 0.008 %* Mn 1.745 ± 0.009 %* P 0.0258 ± 0.0008 %* S 0.0237 ± 0.0005 %* Cr 16.811 ± 0.019 %* Mo 2.111 ± 0.010 %* Ni 10.72 ± 0.05 %* Al 0.0027 ± 0.0004 %* As 0.0063 ± 0.0003 %* B 0.0026 ± 0.0001 %* Co 0.0525 ± 0.0011 %* Cu 0.1831 ± 0.0014 %* N 0.0151 ± 0.0002 %* Nb (0.0028) % Sn 0.0047 ± 0.0002 %* Ti 0.191 ± 0.004 %* V 0.0425 ± 0.0016 %* W (0.0183) % Zr (0.0005) % O 0.0099 ± 0.0007 ¹ % Ta (0.0013) %	100 g
	(Values in parenthesis are indicative values)	
	* 95%-confidence interval	
	¹ Oxygen content only certified for the chips	
ECRM-D 284-2-D	Stainless steel, 17% chromium, 11% nickel and 2% molybdenum, no. 1.4571 - disc (35-39 mm x 20-35 mm)	disc
ECRM-D 286-1	Stainless steel, 18% chromium, 8.5% nickel, 0.3% molybdenum and 0.3% sulfur, no. 1.4305 - chips Year of issue: 1985 Certified values C 0.100 ± 0.005 Mn 1.92 ± 0.03 % P 0.026 ± 0.002 % S 0.280 ± 0.014 % Cr 18.13 ± 0.08 % Mo 0.329 ± 0.009 % Ni 8.54 ± 0.04 % Al (0.0023) % B (0.0003) % Co 0.150 ± 0.008 % N 0.043 ± 0.002 % Pb (0.0003) % Sn 0.0084 ± 0.0009 % O (0.0315) % Sb 0.0014 ± 0.0004 %	100 g
	(Values in parenthesis are indicative values)	
ECRM-D 288-1	Steel with 2% C and 12% Cr, no. 1.2436 - chips Year of issue: 1986 Certified values C 2.08 ± 0.02 % Si 0.260 ± 0.012 % Mn 0.292 ± 0.008 % P 0.024 ± 0.002 % S (0.0012) % Cr 12.00 ± 0.08 % Mo 0.103 ± 0.007 % Ni 0.298 ± 0.007 % Al 0.012 ± 0.002 % As (0.0065) % Co 0.018 ± 0.002 % Cu 0.060 ± 0.004 % N 0.0151 ± 0.0004 % Sn (0.0043) % Ti 0.020 ± 0.002 % V 0.055 ± 0.004 % W (0.682) % Sb (0.0014) %	100 g
	(Values in parenthesis are indicative values)	
ECRM-D 288-1-D	Steel with 2% C and 12% Cr, no. 1.2436 - disc (35-39 mm x 20-35 mm)	disc

Iron and steel products

Code	Product	Unit
ECRM-D 289-1	Heat-resisting steel - chips Year of issue: 1990 Certified values C 0.0489 ± 0.0022 % Mo 1.102 ± 0.015 % Pb (0.0008) % Si 0.531 ± 0.013 % Ni 24.68 ± 0.19 % Sn 0.111 ± 0.010 % Mn 1.016 ± 0.016 % Al 0.199 ± 0.011 % Ti 2.01 ± 0.05 % P 0.0114 ± 0.0010 % As (0.0056) % V 0.260 ± 0.015 % S 0.0027 ± 0.0004 % B 0.0044 ± 0.0004 % Sb (0.0013) % Cr 14.63 ± 0.11 % Co 0.065 ± 0.006 % (Values in parenthesis are indicative values)	100 g
ECRM-D 289-1-D	Heat-resisting steel - disc (35-39 mm x 20-35 mm)	disc
ECRM-D 290-1	High-speed steel, 4% chromium, 5% cobalt und 5% molybdenum, no. 1.3243 - powder Year of issue: 1990 Certified values C 0.911 ± 0.010 % Cr 4.18 ± 0.06 % N 0.0325 ± 0.0012 % Si 0.072 ± 0.007 % Mo 4.83 ± 0.09 % V 1.91 ± 0.04 % Mn 0.244 ± 0.010 % Ni 0.329 ± 0.018 % W 6.27 ± 0.14 % P 0.0160 ± 0.0005 % Co 5.12 ± 0.12 % S 0.0160 ± 0.0008 % Cu 0.081 ± 0.004 %	100 g
ECRM-D 290-1-D	High-speed steel, 4% chromium, 5% cobalt und 5% molybdenum, no. 1.3243 - disc (35-39 mm x 20-35 mm)	disc
ECRM-D 291-1	Highly alloyed steel, 17% chromium and 2% molybdenum, no. 1.4748 - powder Year of issue: 1990 Certified values C 0.903 ± 0.008 % Cr 17.10 ± 0.10 % Cu 0.0711 ± 0.0019 % Si 0.907 ± 0.018 % Mo 2.10 ± 0.06 % N 0.1142 ± 0.0038 % Mn 0.808 ± 0.011 % Ni 0.563 ± 0.011 % Nb (0.0057) % P 0.0168 ± 0.0016 % Al 0.0030 ± 0.0006 % V 0.388 ± 0.016 % S 0.0087 ± 0.0007 % Co 0.0233 ± 0.0022 % (Values in parenthesis are indicative values)	100 g
ECRM-D 291-1-D	Highly alloyed steel, 17% chromium and 2% molybdenum, no. 1.4748 - disc (35-39 mm x 20-35 mm)	disc
ECRM-D 294-1	Manganese steel, material-no. 1.3816 - chips Year of issue: 2005 Certified values C 0.0657 ± 0.0010 %* Al (0.0095) % Sn (0.0014) % Si 0.283 ± 0.005 %* As 0.00365 ± 0.00029 %* Ti (0.0008) % Mn 18.68 ± 0.04 %* B (<0.00005) % V 0.0694 ± 0.0021 %* P 0.0273 ± 0.0013 %* Co 0.0288 ± 0.0009 % W (0.00114) % S 0.00031 ± 0.00009 %* Cu 0.0242 ± 0.0007 %* Zr (0.0001) % Cr 17.98 ± 0.05 %* N 0.566 ± 0.011 %* Ca (0.00026) % Mo 0.0861 ± 0.0022 %* Nb (0.00117) % Sb (0.00053) % Ni 0.427 ± 0.006 %* Pb (0.000128) % Te (<0.00008) % (Values in parenthesis are indicative values) * 95%-confidence interval	100 g
ECRM-D 294-1-D	Manganese steel, material-no. 1.3816 - disc (40 mm x 20 mm)	disc
ECRM-D 297-1	1.2 % boron steel, material-no. 1.4696 - chips Year of issue: 2005 Certified values C 0.0223 ± 0.0004 %* Ni 12.33 ± 0.02 %* Nb (0.0089) % Si 0.344 ± 0.006 %* Al 0.0195 ± 0.0009 %* Ti 0.0072 ± 0.0004 %* Mn 0.897 ± 0.007 %* As 0.0040 ± 0.0005 %* V 0.0535 ± 0.0008 %* P 0.0135 ± 0.0004 %* B 1.1462) ± 0.009 %* W (0.0057) % S 0.0101 ± 0.0003 %* Co 0.0413 ± 0.0007 %* Zr (0.0002) % Cr 18.37 ± 0.03 %* Cu 0.204 ± 0.004 %* Ca (0.0002) % Mo 0.290 ± 0.005 %* N 0.0152 ± 0.0007 %* (Values in parenthesis are indicative values) * 95%-confidence interval	100 g
ECRM-D 297-1-D	1.2 % boron steel, material-no. 1.4696 - disc (40 mm x 25 mm)	disc
ECRM-B 251-1	High-speed steel, 5,4% chromium, 5,7% cobalt, 0,5% molybdenum, 1,6% vanadium und 20% tungsten - chips Year of issue: 1968 Certified values C 0.84 % Cr 5.35 % Sn 0.025 % Si 0.21 % Mo 0.53 % V 1.59 % Mn 0.27 % Ni 0.15 % W 19.9 P 0.024 % Co 5.70 % S 0.025 % Cu 0.08 %	100 g

Iron and steel products

Code	Product	Unit
ECRM-B 253-1	13% Manganese steel - chips Year of issue: 1973 Certified values C 1.15 % S 0.019 % Co (0.35) % Si 0.34 % Cr 0.16 % Cu (0.177) % Mn 12.5 % Mo (0.031) % V (0.02) % P 0.042 % Ni 0.29 % (Values in parenthesis are indicative values)	100 g
ECRM-B 254-1	High-speed steel, 5,1% chromium, 4,9% molybdenum, 1,9% vanadium and 7% tungsten - chips Year of issue: 1971 Certified values C 0.88 % Cr 5.12 % Sn 0.019 % Si 0.19 % Mo 4.92 % V 1.94 % Mn 0.30 % Ni 0.12 % W 6.97 % P 0.023 % Co 0.32 % S 0.029 % Cu 0.09 %	100 g
ECRM-B 254-1-D	High-speed steel, 5,1% chromium, 4,9% molybdenum, 1,9% Vanadium and 7% tungsten	disc
ECRM-B 272-1	12% Chromium steel - chips Year of issue: 2005 Certified values C 0.2815 ± 0.0012 %* Ni 0.244 ± 0.003 %* Nb 0.0028 ± 0.0002 %* Si 0.420 ± 0.004 %* Al 0.0046 ± 0.0005 %* Ti 0.00096 ± 0.00009 %* Mn 0.600 ± 0.004 %* As 0.0116 ± 0.0005 %* V 0.0167 ± 0.0005 %* P 0.0156 ± 0.0005 %* B 0.0018 ± 0.0002 %* Ca 0.00090 ± 0.00011 %* S 0.0196 ± 0.0004 %* Co 0.0145 ± 0.0003 %* Sb 0.000070 ± 0.00006 %* Cr 11.927 ± 0.022 %* Cu 0.0192 ± 0.0004 %* Zn 0.0031 ± 0.0003 %* Mo 0.0030 ± 0.0002 %* N 0.0508 ± 0.0004 %*	100 g
ECRM-B 272-1-D	12% Chromium steel - disc (35-39 mm x 20-35 mm)	disc
ECRM-B 276-2	5% Cr-Mo-V steel - chips Year of issue: 1993 Certified values C 0.399 ± 0.002 %* S 0.0189 ± 0.0003 %* Cu 0.183 ± 0.002 %* Si 1.034 ± 0.004 %* Cr 4.975 ± 0.012 %* N 0.0116 ± 0.0001 %* Mn 0.365 ± 0.002 %* Mo 1.134 ± 0.010 %* Sn 0.0133 ± 0.0002 %* P 0.0093 ± 0.0003 %* Ni 0.203 ± 0.002 %* V 0.296 ± 0.005 %*	100 g
ECRM-B 276-2-D	5% Cr-Mo-V steel	disc
ECRM-B 281-1	Stainless steel, 18% chromium and 9% nickel - chips Year of issue: 1980 Certified values C 0.048 ± 0.002 % Cr 18.17 ± 0.05 % Cu 0.076 ± 0.003 % Si 0.929 ± 0.008 % Ni 9.37 ± 0.05 % N 0.023 ± 0.001 % Mn 0.786 ± 0.007 % Al 0.015 ± 0.001 % Pb (0.0005) % P 0.012 ± 0.001 % B (0.0012) % Sn 0.009 ± 0.001 % S 0.016 ± 0.001 % Co 0.023 ± 0.002 % Ti 0.216 ± 0.007 % (Values in parenthesis are indicative values)	100 g
ECRM-B 285-2	Maraging steel - chips Year of issue: 1997 Certified values C 0.0018 ± 0.0002 %* Cr 0.0236 ± 0.0011 %* Co 7.76 ± 0.04 %* Si 0.0117 ± 0.0008 %* Mo 4.99 ± 0.02 %* Cu 0.0094 ± 0.0004 %* Mn 0.0168 ± 0.0004 %* Ni 18.07 ± 0.04 %* N 0.0007 ± 0.0001 %* P 0.0053 ± 0.0005 %* Al 0.1067 ± 0.0018 %* Ti 0.520 ± 0.004 %* S 0.0025 ± 0.0002 %* B 0.0009 ± 0.0001 %* Zr 0.0050 ± 0.0001 %*	100 g
ECRM-B 285-2-D	Maraging steel - disc (35-39 mm x 20-35 mm)	disc
ECRM-B 287-1	High boron stainless steel - chips Year of issue: 1986 Certified values C 0.016 ± 0.0005 %* S 0.0014 ± 0.0002 %* B 0.924R ± 0.003 %* Si 0.569 ± 0.009 %* Cr 18.61 ± 0.04 %* Co 0.148 ± 0.002 %* Mn 1.478 ± 0.009 %* Mo 0.247 ± 0.004 %* Cu 0.203 ± 0.002 %* P 0.0267 ± 0.0006 %* Ni 10.35 ± 0.04 %* N 0.0194 ± 0.0005 %*	100 g
ECRM-B 287-1-D	High boron stainless steel - disc (35-39 mm x 20-35 mm)	disc

Iron and steel products

Code	Product	Unit
ECRM-B 292-1	Nb-stabilised stainless steel - chips Year of issue: 1990 Certified values C 0.0367 ± 0.0008 %* Mo..... 0.0464 ± 0.0011 %* Cu 0.0391 ± 0.0010 %* Si..... 0.402 ± 0.005 %* Ni..... 10.09 ± 0.02 %* N 0.0640 ± 0.0012 %* Mn 1.744 ± 0.006 %* Al (0.002) % Nb 0.571 ± 0.012 %* P 0.0175 ± 0.0007 %* As (0.008) % Ca (0.0006) % S 0.0055 ± 0.0002 %* B (0.0003) % Ta (0.001) % Cr 18.00 ± 0.02 %* Co 0.0255 ± 0.0011 %* (Values in parenthesis are indicative values)	100 g
	* 95%-confidence interval	
ECRM-B 292-1-D	Nb-stabilised stainless steel - disc (35-39 mm x 20-35 mm)	disc
ECRM-B 295-1	4% Mo-Cr-Ni steel - chips Year of issue: 1995 Certified values C 0.0166 ± 0.0004 %* Mo..... 3.996 ± 0.024 %* Cu 1.481 ± 0.012 %* Si..... 0.418 ± 0.008 %* Ni..... 24.40 ± 0.04 %* N 0.0615 ± 0.0008 %* Mn 1.758 ± 0.009 %* Al 0.0203 ± 0.0005 %* Sn 0.0025 ± 0.0003 %* P 0.0167 ± 0.0005 %* As 0.0041 ± 0.0002 %* V 0.0456 ± 0.0015 %* S 0.0003 ± 0.0001 %* B 0.0018 ± 0.0001 %* Fe 48.36 ± 0.11 %* Cr 19.51 ± 0.06 %* Co 0.0450 ± 0.0011 %* Sb 0.0007 ± 0.0001 %*	100 g
	* 95%-confidence interval	
ECRM-B 295-1-D	4% Mo-Cr-Ni steel - disc (35-39 mm x 20-35 mm)	disc
ECRM-B 296-1	Jethete steel - chips Year of issue: 1997 Certified values C 0.1166 ± 0.0008 %* Mo..... 1.700 ± 0.009 %* Cu 0.1498 ± 0.0013 %* Si..... 0.242 ± 0.003 %* Ni..... 2.790 ± 0.008 %* N 0.0214 ± 0.0005 %* Mn 0.676 ± 0.005 %* Al 0.0276 ± 0.0009 %* Pb 0.00016 ± 0.00004 %* P 0.0178 ± 0.0004 %* As 0.0139 ± 0.0006 %* Sn 0.0131 ± 0.0003 %* S 0.0026 ± 0.0001 %* B (0.0003) % V 0.363 ± 0.005 %* Cr 11.82 ± 0.04 %* Co 0.0218 ± 0.0006 %* Mg (0.0003) % (Values in parenthesis are indicative values)	100 g
	* 95%-confidence interval	
ECRM-B 296-1-D	Jethete steel - disc (35-39 mm x 20-35 mm)	disc
ECRM-F 201-1	Stainless steel, 12% chromium - chips Year of issue: 1960 Certified values C 0.291 % Ni 0.202 % Mg (0.0005) % Si 0.843 % Cu 0.099 % Sb (0.0033) % Mn 0.363 % Sn 0.0193 % Zn 0.0005 % Cr 12.33 % V (0.020) % Mo 0.050 % Ca 0.0018 %	100 g
	(Values in parenthesis are indicative values)	
ECRM-F 210-1	Stainless steel, 4% chromium, 8.2% molybdenum, 1.7% vanadium and 1.5% tungsten - chips Year of issue: 1967 Certified values C 0.762 % S (0.022) % V 1.65 % Si 0.200 % Cr 3.92 % W 1.54 % Mn 0.250 % Mo 8.15 % P (0.028) % Co 0.185 %	100 g
	(Values in parenthesis are indicative values)	
ECRM-F 279-2	Stainless steel, 15.6% chromium and 11% nickel - chips Year of issue: 1972 Certified values C 0.0885 ± 0.0010 % Cr 15.642 ± 0.030 % Si 0.5160 ± 0.0047 % Ni 1.603 ± 0.014 % Mn 0.2584 ± 0.0045 % Cu 0.1067 ± 0.0022 %	100 g
ECRM-F 282-1	Stainless steel, 17% chromium, 11% nickel, 2% molybdenum and 0.5% titanium - chips Year of issue: 1986 Certified values C 0.086 ± 0.003 % S 0.0042 ± 0.0005 % Cu 0.109 ± 0.003 % Si 0.734 ± 0.016 % Cr 16.72 ± 0.04 % Ti 0.488 ± 0.007 % Mn 1.64 ± 0.02 % Mo 2.19 ± 0.03 % P 0.019 ± 0.002 % Ni 10.86 ± 0.09 %	100 g

Iron and steel products

Code	Product	Unit																																																						
ECRM-F 273-1	<p>Stainless steel Z5CNU15.05 with 14.7 % Cr, 4.8 % Ni and 3 % Cu - chips</p> <p>Year of issue: 2005</p> <p>Certified values</p> <table> <tbody> <tr><td>C</td><td>0.0336 ± 0.0005 %*</td><td>Mo</td><td>0.246 ± 0.003 %*</td><td>N</td><td>0.0444 ± 0.0004 %*</td></tr> <tr><td>Si.....</td><td>0.378 ± 0.004 %*</td><td>Ni.....</td><td>4.852 ± 0.020 %*</td><td>Nb.....</td><td>0.221 ± 0.004 %*</td></tr> <tr><td>Mn.....</td><td>0.785 ± 0.004 %*</td><td>Al.....</td><td>(0.0062) %</td><td>Sn.....</td><td>0.0021 ± 0.0001 %*</td></tr> <tr><td>P</td><td>0.0131 ± 0.0006 %*</td><td>As</td><td>0.0030 ± 0.0002 %*</td><td>V</td><td>0.0512 ± 0.0008 %*</td></tr> <tr><td>S</td><td>0.00037 ± 0.00006 %*</td><td>Co</td><td>0.0391 ± 0.0005 %*</td><td>Ga</td><td>(0.0023) %</td></tr> <tr><td>Cr.....</td><td>14.747 ± 0.014 %*</td><td>Cu.....</td><td>3.047 ± 0.015 %*</td><td></td><td></td></tr> </tbody> </table> <p>(Values in parenthesis are indicative values)</p> <p>* 95%-confidence interval</p>	C	0.0336 ± 0.0005 %*	Mo	0.246 ± 0.003 %*	N	0.0444 ± 0.0004 %*	Si.....	0.378 ± 0.004 %*	Ni.....	4.852 ± 0.020 %*	Nb.....	0.221 ± 0.004 %*	Mn.....	0.785 ± 0.004 %*	Al.....	(0.0062) %	Sn.....	0.0021 ± 0.0001 %*	P	0.0131 ± 0.0006 %*	As	0.0030 ± 0.0002 %*	V	0.0512 ± 0.0008 %*	S	0.00037 ± 0.00006 %*	Co	0.0391 ± 0.0005 %*	Ga	(0.0023) %	Cr.....	14.747 ± 0.014 %*	Cu.....	3.047 ± 0.015 %*			100 g																		
C	0.0336 ± 0.0005 %*	Mo	0.246 ± 0.003 %*	N	0.0444 ± 0.0004 %*																																																			
Si.....	0.378 ± 0.004 %*	Ni.....	4.852 ± 0.020 %*	Nb.....	0.221 ± 0.004 %*																																																			
Mn.....	0.785 ± 0.004 %*	Al.....	(0.0062) %	Sn.....	0.0021 ± 0.0001 %*																																																			
P	0.0131 ± 0.0006 %*	As	0.0030 ± 0.0002 %*	V	0.0512 ± 0.0008 %*																																																			
S	0.00037 ± 0.00006 %*	Co	0.0391 ± 0.0005 %*	Ga	(0.0023) %																																																			
Cr.....	14.747 ± 0.014 %*	Cu.....	3.047 ± 0.015 %*																																																					
ECRM-F-273-1-D	Stainless steel Z5CNU15.05 with 14.7 % Cr, 4.8 % Ni and 3 % Cu - disc (40 mm x 20 mm)	disc																																																						
ECRM-S 274-1	<p>Vanadium steel - chips</p> <p>Year of issue: 2004</p> <p>Certified values</p> <table> <tbody> <tr><td>C</td><td>1.563 ± 0.006 %*</td><td>Ni</td><td>0.077 ± 0.002 %*</td><td>Pb</td><td>(0.0064) %</td></tr> <tr><td>Si.....</td><td>1.057 ± 0.008 %*</td><td>Al</td><td>(0.0025) %</td><td>Sn</td><td>(0.010) %</td></tr> <tr><td>Mn.....</td><td>0.397 ± 0.004 %*</td><td>As</td><td>(0.0013) %</td><td>Ti</td><td>(0.0011) %</td></tr> <tr><td>P</td><td>0.0148 ± 0.0005 %*</td><td>B</td><td>(0.0005) %</td><td>V</td><td>4.010 ± 0.018 %*</td></tr> <tr><td>S</td><td>0.0096 ± 0.0004 %*</td><td>Co</td><td>(0.0229) %</td><td>W</td><td>0.0087 ± 0.0007 %*</td></tr> <tr><td>Cr.....</td><td>8.036 ± 0.022 %*</td><td>Cu.....</td><td>0.0281 ± 0.0005 %*</td><td>O</td><td>0.0026 ± 0.0002+1) %</td></tr> <tr><td>Mo.....</td><td>1.455 ± 0.016 %*</td><td>N.....</td><td>0.0769 ± 0.0010 %*</td><td>Sb</td><td>(0.00019) %</td></tr> </tbody> </table> <p>(Values in parenthesis are indicative values)</p> <p>* 95%-confidence interval</p>	C	1.563 ± 0.006 %*	Ni	0.077 ± 0.002 %*	Pb	(0.0064) %	Si.....	1.057 ± 0.008 %*	Al	(0.0025) %	Sn	(0.010) %	Mn.....	0.397 ± 0.004 %*	As	(0.0013) %	Ti	(0.0011) %	P	0.0148 ± 0.0005 %*	B	(0.0005) %	V	4.010 ± 0.018 %*	S	0.0096 ± 0.0004 %*	Co	(0.0229) %	W	0.0087 ± 0.0007 %*	Cr.....	8.036 ± 0.022 %*	Cu.....	0.0281 ± 0.0005 %*	O	0.0026 ± 0.0002+1) %	Mo.....	1.455 ± 0.016 %*	N.....	0.0769 ± 0.0010 %*	Sb	(0.00019) %	100 g												
C	1.563 ± 0.006 %*	Ni	0.077 ± 0.002 %*	Pb	(0.0064) %																																																			
Si.....	1.057 ± 0.008 %*	Al	(0.0025) %	Sn	(0.010) %																																																			
Mn.....	0.397 ± 0.004 %*	As	(0.0013) %	Ti	(0.0011) %																																																			
P	0.0148 ± 0.0005 %*	B	(0.0005) %	V	4.010 ± 0.018 %*																																																			
S	0.0096 ± 0.0004 %*	Co	(0.0229) %	W	0.0087 ± 0.0007 %*																																																			
Cr.....	8.036 ± 0.022 %*	Cu.....	0.0281 ± 0.0005 %*	O	0.0026 ± 0.0002+1) %																																																			
Mo.....	1.455 ± 0.016 %*	N.....	0.0769 ± 0.0010 %*	Sb	(0.00019) %																																																			
ECRM-S 274-1-D	Vanadium steel - disc (38 mm x 25 mm)	disc																																																						
ECRM-S 298-1	<p>Duplex stainless steel (similar to EN 1.4410) - chips</p> <p>Year of issue: 2002</p> <p>Certified values</p> <table> <tbody> <tr><td>C</td><td>0.0146 ± 0.0003 %*</td><td>As</td><td>(0.0034) %</td><td>V</td><td>0.0607 ± 0.0011 %*</td></tr> <tr><td>Si.....</td><td>0.262 ± 0.005 %*</td><td>B</td><td>0.0021 ± 0.0001 %*</td><td>W</td><td>(0.0180) %</td></tr> <tr><td>Mn.....</td><td>0.398 ± 0.005 %*</td><td>Co</td><td>0.055 ± 0.001 %*</td><td>Ca</td><td>(0.0020) %</td></tr> <tr><td>P</td><td>0.0197 ± 0.0004 %*</td><td>Cu</td><td>0.201 ± 0.003 %*</td><td>Fe</td><td>63.38 ± 0.11 %*</td></tr> <tr><td>S</td><td>0.0006 ± 0.0001 %*</td><td>N</td><td>0.263 ± 0.0005 %*</td><td>Mg</td><td>(0.0008) %</td></tr> <tr><td>Cr.....</td><td>24.72 ± 0.04 %*</td><td>Nb</td><td>(0.0037) %</td><td>O</td><td>(0.0036) %</td></tr> <tr><td>Mo.....</td><td>3.799 ± 0.026 %*</td><td>Pb</td><td>0.00008 ± 0.00002 %*</td><td>Sb</td><td>(0.00072) %</td></tr> <tr><td>Ni.....</td><td>7.056 ± 0.018 %*</td><td>Sn</td><td>(0.0044) %</td><td>Ta</td><td>(0.00002) %</td></tr> <tr><td>Al.....</td><td>0.0285 ± 0.0011 %*</td><td>Ti</td><td>0.0014 ± 0.0002 %*</td><td>Zn</td><td>(0.0010) %</td></tr> </tbody> </table> <p>(Values in parenthesis are indicative values)</p> <p>* 95%-confidence interval</p>	C	0.0146 ± 0.0003 %*	As	(0.0034) %	V	0.0607 ± 0.0011 %*	Si.....	0.262 ± 0.005 %*	B	0.0021 ± 0.0001 %*	W	(0.0180) %	Mn.....	0.398 ± 0.005 %*	Co	0.055 ± 0.001 %*	Ca	(0.0020) %	P	0.0197 ± 0.0004 %*	Cu	0.201 ± 0.003 %*	Fe	63.38 ± 0.11 %*	S	0.0006 ± 0.0001 %*	N	0.263 ± 0.0005 %*	Mg	(0.0008) %	Cr.....	24.72 ± 0.04 %*	Nb	(0.0037) %	O	(0.0036) %	Mo.....	3.799 ± 0.026 %*	Pb	0.00008 ± 0.00002 %*	Sb	(0.00072) %	Ni.....	7.056 ± 0.018 %*	Sn	(0.0044) %	Ta	(0.00002) %	Al.....	0.0285 ± 0.0011 %*	Ti	0.0014 ± 0.0002 %*	Zn	(0.0010) %	100 g
C	0.0146 ± 0.0003 %*	As	(0.0034) %	V	0.0607 ± 0.0011 %*																																																			
Si.....	0.262 ± 0.005 %*	B	0.0021 ± 0.0001 %*	W	(0.0180) %																																																			
Mn.....	0.398 ± 0.005 %*	Co	0.055 ± 0.001 %*	Ca	(0.0020) %																																																			
P	0.0197 ± 0.0004 %*	Cu	0.201 ± 0.003 %*	Fe	63.38 ± 0.11 %*																																																			
S	0.0006 ± 0.0001 %*	N	0.263 ± 0.0005 %*	Mg	(0.0008) %																																																			
Cr.....	24.72 ± 0.04 %*	Nb	(0.0037) %	O	(0.0036) %																																																			
Mo.....	3.799 ± 0.026 %*	Pb	0.00008 ± 0.00002 %*	Sb	(0.00072) %																																																			
Ni.....	7.056 ± 0.018 %*	Sn	(0.0044) %	Ta	(0.00002) %																																																			
Al.....	0.0285 ± 0.0011 %*	Ti	0.0014 ± 0.0002 %*	Zn	(0.0010) %																																																			
ECRM-S 298-1-D	Duplex stainless steel (similar to EN 1.4410) - disc (38 mm x 25 mm)	disc																																																						
Pure iron																																																								
ECRM-D 098-1	<p>Pure iron - chips</p> <p>Year of issue: 1993</p> <p>Certified values</p> <table> <tbody> <tr><td>C</td><td>5.1 ± 1.3 µg/g*</td><td>P</td><td>(0.6) µg/g</td><td>Mo</td><td>8.5 ± 0.8 µg/g*</td></tr> <tr><td>Si.....</td><td>4.8 ± 1.1 µg/g*</td><td>S</td><td>3.1 ± 0.5 µg/g*</td><td>N</td><td>2.4 ± 0.7 µg/g*</td></tr> <tr><td>Mn.....</td><td>0.8 ± 0.4 µg/g*</td><td>Cr.....</td><td>57.1 ± 2.4 µg/g*</td><td></td><td></td></tr> </tbody> </table> <p>(Values in parenthesis are indicative values)</p> <p>* 95%-confidence interval</p>	C	5.1 ± 1.3 µg/g*	P	(0.6) µg/g	Mo	8.5 ± 0.8 µg/g*	Si.....	4.8 ± 1.1 µg/g*	S	3.1 ± 0.5 µg/g*	N	2.4 ± 0.7 µg/g*	Mn.....	0.8 ± 0.4 µg/g*	Cr.....	57.1 ± 2.4 µg/g*			100 g																																				
C	5.1 ± 1.3 µg/g*	P	(0.6) µg/g	Mo	8.5 ± 0.8 µg/g*																																																			
Si.....	4.8 ± 1.1 µg/g*	S	3.1 ± 0.5 µg/g*	N	2.4 ± 0.7 µg/g*																																																			
Mn.....	0.8 ± 0.4 µg/g*	Cr.....	57.1 ± 2.4 µg/g*																																																					
ECRM-B 088-2	<p>Pure iron - powder</p> <p>Year of issue: 2001</p> <p>Certified values</p> <table> <tbody> <tr><td>C</td><td>6 ± 2 µg/g*</td><td>S</td><td>70 ± 4 µg/g*</td><td>Cu</td><td>163 ± 3 µg/g*</td></tr> <tr><td>Si.....</td><td>52 ± 5 µg/g*</td><td>Cr.....</td><td>244 ± 8 µg/g*</td><td>V</td><td>2.9 ± 0.5 µg/g*</td></tr> <tr><td>Mn.....</td><td>809 ± 17 µg/g*</td><td>Ni</td><td>275 ± 7 µg/g*</td><td>Ca</td><td>7.2 ± 0.8 µg/g*</td></tr> <tr><td>P</td><td>48 ± 2 µg/g*</td><td>Co</td><td>61 ± 2 µg/g*</td><td></td><td></td></tr> </tbody> </table> <p>* 95%-confidence interval</p>	C	6 ± 2 µg/g*	S	70 ± 4 µg/g*	Cu	163 ± 3 µg/g*	Si.....	52 ± 5 µg/g*	Cr.....	244 ± 8 µg/g*	V	2.9 ± 0.5 µg/g*	Mn.....	809 ± 17 µg/g*	Ni	275 ± 7 µg/g*	Ca	7.2 ± 0.8 µg/g*	P	48 ± 2 µg/g*	Co	61 ± 2 µg/g*			100 g																														
C	6 ± 2 µg/g*	S	70 ± 4 µg/g*	Cu	163 ± 3 µg/g*																																																			
Si.....	52 ± 5 µg/g*	Cr.....	244 ± 8 µg/g*	V	2.9 ± 0.5 µg/g*																																																			
Mn.....	809 ± 17 µg/g*	Ni	275 ± 7 µg/g*	Ca	7.2 ± 0.8 µg/g*																																																			
P	48 ± 2 µg/g*	Co	61 ± 2 µg/g*																																																					
ECRM-B 097-1	<p>Pure iron - chips</p> <p>Year of issue: 1984</p> <p>Certified values</p> <table> <tbody> <tr><td>C</td><td>2.5 ± 0.7 µg/g</td><td>S</td><td>22 ± 2 µg/g</td><td>B</td><td>3 ± 1 µg/g</td></tr> <tr><td>Si.....</td><td>(<100) µg/g</td><td>Cr.....</td><td>16 ± 5 µg/g</td><td>Co</td><td>37 ± 4 µg/g</td></tr> <tr><td>Mn.....</td><td>64 ± 4 µg/g</td><td>Ni</td><td>25 ± 3 µg/g</td><td>Cu</td><td>20 ± 2 µg/g</td></tr> <tr><td>P</td><td>16R ± 3 µg/g</td><td>As</td><td>51 ± 5 µg/g</td><td>N</td><td>7 ± 1 µg/g</td></tr> </tbody> </table> <p>(Values in parenthesis are indicative values)</p> <p>R: revised value</p>	C	2.5 ± 0.7 µg/g	S	22 ± 2 µg/g	B	3 ± 1 µg/g	Si.....	(<100) µg/g	Cr.....	16 ± 5 µg/g	Co	37 ± 4 µg/g	Mn.....	64 ± 4 µg/g	Ni	25 ± 3 µg/g	Cu	20 ± 2 µg/g	P	16R ± 3 µg/g	As	51 ± 5 µg/g	N	7 ± 1 µg/g	100 g																														
C	2.5 ± 0.7 µg/g	S	22 ± 2 µg/g	B	3 ± 1 µg/g																																																			
Si.....	(<100) µg/g	Cr.....	16 ± 5 µg/g	Co	37 ± 4 µg/g																																																			
Mn.....	64 ± 4 µg/g	Ni	25 ± 3 µg/g	Cu	20 ± 2 µg/g																																																			
P	16R ± 3 µg/g	As	51 ± 5 µg/g	N	7 ± 1 µg/g																																																			
ECRM-B 097-1-D	Pure iron - disc (35-39 mm x 20-35 mm)	disc																																																						

Iron and steel products

Code	Product	Unit
Special alloys		
ECRM-D 326-1	Nickel-base alloy with 61,2% nickel content and 16,4% chromium content - chips	100 g
Year of issue: 1972		
Certified values		
C 0.092 ± 0.002 %	Cr 16.37 ± 0.05 %	Cu (0.027) %
Si 1.46 ± 0.025 %	Mo (0.025) %	N (0.0359) %
Mn 0.406 ± 0.008 %	Ni 61.16 ± 0.16 %	V (0.024) %
P 0.0093 ± 0.0009 %	Alttotal (0.79) %	Zr 0.129 ± 0.008 %
S 0.0028 ± 0.0006 %	Co 0.223 ± 0.011 %	
(Values in parenthesis are indicative values)		
ECRM-D 327-2	High-temperature steel with 24,4% chromium, 19,7% nickel, 2% silicon and 0,15% carbon - chips	100 g
Year of issue: 1972		
Certified values		
C 0.152 ± 0.003 %	Cr 24.35 ± 0.08 %	Cu 0.060 ± 0.003 %
Si 2.052 ± 0.028 %	Mo 0.174 ± 0.009 %	N 0.059 ± 0.0024 %
Mn 1.289 ± 0.018 %	Ni 19.72 ± 0.08 %	V 0.044 ± 0.004 %
P 0.0228 ± 0.0014 %	Al _{total} 0.070 ± 0.006 %	
S 0.0046 ± 0.0012 %	Co 0.159 ± 0.010 %	
(Values in parenthesis are indicative values)		
ECRM-D 328-1	High-temperature turbine steel with 20,5% chromium, 20,4% nickel and 41,7% cobalt - chips	100 g
Year of issue: 1973		
Certified values		
C 0.390 ± 0.005 %	Mo 4.41 ± 0.07 %	Nb 3.61 ± 0.22 %
Si 0.629 ± 0.014 %	Ni 20.38 ± 0.19 %	W 4.16 ± 0.04 %
Mn 1.395 ± 0.012 %	Al _{total} 0.070 ± 0.006 %	Fe 2.40 ± 0.06 %
P 0.005 ± 0.001 %	Co 41.65 ± 0.24 %	Ta 0.18 ± 0.02 %
S 0.003 ± 0.001 %	Cu 0.013 ± 0.003 %	
Cr 20.54 ± 0.07 %	N 0.027 ± 0.002 %	
ECRM-B 376-1	24% Co magnet alloy	100 g
Year of issue: 1980		
Certified values		
C 0.026 ± 0.002 %*	Ni 13.37 ± 0.08 %*	Nb 0.309 ± 0.012 %*
Si 0.313 ± 0.009 %*	Al 8.12 ± 0.10 %*	Ti 0.158R± 0.003 %*
Mn 0.046 ± 0.004 %*	Co 23.70 ± 0.19 %*	Ta (0.0155) %
S 0.004 ± 0.001 %*	Cu 2.94 ± 0.04 %*	
(Values in parenthesis are indicative values)		
R: revised value		
* 95%-confidence interval		
ECRM-F 377-1	Nickel alloy IN 625 - disc (40 mm x 20 mm)	disc
Year of issue: 1995		
Certified values		
C 0.0202 ± 0.0007 %*	Cr 21.72 ± 0.06 %*	Co 0.0348 ± 0.0009 %*
Si 0.077 ± 0.005 %*	Mo 8.94 ± 0.05 %*	Cu 0.0110 ± 0.0003 %*
Mn 0.0225 ± 0.0007 %*	Ni 61.45 ± 0.08 %*	Nb 3.50 ± 0.04 %*
P 0.0036 ± 0.0004 %*	Al 0.216 ± 0.007 %*	Ti 0.255 ± 0.008 %*
S 0.0006 ± 0.0001 %*	B (0.0006) %	Fe 3.77 ± 0.06 %*
* 95%-confidence interval		
ECRM-F 377-2	Nickel alloy IN 625 - disc (40 mm x 20 mm)	disc
Year of issue: 1995		
Certified values		
C 0.0202 ± 0.0007 %*	Cr 21.72 ± 0.06 %*	Co 0.0348 ± 0.0009 %*
Si 0.077 ± 0.005 %*	Mo 8.94 ± 0.05 %*	Cu 0.0104 ± 0.0003 %*
Mn 0.0225 ± 0.0007 %*	Ni 61.45 ± 0.08 %*	Nb 3.50 ± 0.04 %*
P 0.0036 ± 0.0004 %*	Al 0.232 ± 0.007 %*	Ti 0.264 ± 0.008 %*
S 0.0006 ± 0.0001 %*	B (0.0006) %	Fe 3.77 ± 0.06 %*
* 95%-confidence interval		
ECRM-F 378-1	Cobalt alloy, stellite grade 6 - powder	100 g
Year of issue: 1995		
Certified values		
C 1.181 ± 0.009 %*	S 0.0055 ± 0.0006 %*	Co 63.52 ± 0.19 %*
Si 1.172 ± 0.022 %*	Cr 28.22 ± 0.14 %*	W 4.43 ± 0.05 %*
Mn 0.0579 ± 0.0015 %*	Mo 0.0503 ± 0.0011 %*	Fe 0.606 ± 0.011 %*
P (0.0023) %	Ni 0.617 ± 0.011 %*	
(Values in parenthesis are indicative values)		
* 95%-confidence interval		
ECRM-F 378-1-D	Cobalt alloy, stellite grade 6 - disc (40 mm x 20 mm)	disc

Code	Product	Unit
Cast irons		
ECRM-D 428-2	Cast iron GG 25, powdered material produced by atomization of the melt - powder Year of issue: 1998 Certified values C _{total} 2.747 ± 0.009 %* S 0.1105 ± 0.0018 %* As 0.0156 ± 0.0005 %* Si 1.752 ± 0.007 %* Cr 0.0366 ± 0.0017 %* Cu 0.0996 ± 0.0014 %* Mn 0.750 ± 0.05 %* Mo (0.0014) % Ti 0.0311 ± 0.0005 %* P 0.0691 ± 0.0011 %* Ni 0.0358 ± 0.0005 %* V 0.0120 ± 0.0003 %* (Values in parenthesis are indicative values) * 95%-confidence interval	100 g
ECRM-D 476-3	Pig iron - powder Year of issue: 1996 Certified values C _{total} 3.390 ± 0.011 % S 0.0493 ± 0.0009 % Cu 0.2445 ± 0.0025 % Si 1.813 ± 0.005 % Cr 0.0648 ± 0.0012 % N 0.0038 ± 0.0001 % Mn 0.987 ± 0.008 % Ni 0.0549 ± 0.0014 % Ti 0.0222 ± 0.0005 % P 0.0908 ± 0.0023 % As 0.0145 ± 0.0007 % V 0.0115 ± 0.0002 % * 95%-confidence interval	100 g
ECRM-D 478-2	Hematite pig iron - powder Year of issue: 1996 Certified values C _{total} 4.003 ± 0.013 % Cr 0.251 ± 0.005 % N 0.0023 ± 0.0002 % Si 2.411 ± 0.021 % Ni 0.151 ± 0.007 % Ti 0.0328 ± 0.0007 % Mn 0.321 ± 0.005 % As (0.0018) % V 0.0113 ± 0.0003 % P 0.201 ± 0.006 % B 0.0006 ± 0.0001 % S 0.0460 ± 0.0015 % Cu 0.1276 ± 0.0019 % (Values in parenthesis are indicative values) * 95%-confidence interval	100 g
ECRM-D 479-1	Alloyed cast iron - powder Year of issue: 1978 Certified values C _{total} 2.86 ± 0.04 % P 0.076 ± 0.003 % Mo 0.196 ± 0.005 % Si 2.02 ± 0.02 % S 0.089 ± 0.003 % Ni 1.012 ± 0.015 % Mn 0.136 ± 0.008 % Cr 1.00 ± 0.02 % Al 0.014 ± 0.002 % (Values in parenthesis are indicative values)	100 g
ECRM-D 480-1	Nodular iron - powder Year of issue: 1979 Certified values C _{total} 3.03 ± 0.02 % S 0.0086 ± 0.0010 % Cu (0.0052) % Si 2.41 ± 0.02 % Cr (0.0164) % Mg 0.017 ± 0.001 % Mn 0.151 ± 0.005 % Ni 0.483 ± 0.007 % P 0.0021R ± 0.0005 % Al 0.016 ± 0.001 % (Values in parenthesis are indicative values) R: revised value	100 g
ECRM-B 451-2	Austenitic cast iron - chips Year of issue: 1999 Certified values C _{total} 2.059 ± 0.007 %* P 0.0593 ± 0.0011 %* Ni 14.01 ± 0.04 %* Si 2.092 ± 0.007 %* S 0.0315 ± 0.0008 %* Cu 6.26 ± 0.02 %* Mn 1.079 ± 0.08 %* Cr 1.097 ± 0.010 %* (Values in parenthesis are indicative values) * 95%-confidence interval	100 g
ECRM-B 453-1	High Si and P iron - chips Year of issue: 1972 Certified values C _{total} (2.44) % S 0.049 % Cu 0.10 % Si 3.17 % Cr 0.053 % Ti (0.040) % Mn 0.72 % Ni 0.068 % V 0.050 % P 1.63 % As 0.019 % (Values in parenthesis are indicative values)	100 g
ECRM-B 454-1	Hematite pig iron - chips Certified values Year of issue: 1968 C _{total} (2.53) % P 0.046 % As (0.025) % Si 2.00 % S 0.068 % Cu (0.07) % Mn 1.16 % Ni (0.21) % Ti 0.052 % (Values in parenthesis are indicative values)	100 g

Iron and steel products

Code	Product	Unit
ECRM-B 481-1	Nodular iron - chips Year of issue: 1976 Certified values C _{total}3.91 ± 0.02 % S.....0.004 ± 0.001 % Al.....0.023 ± 0.002 % Si.....2.29 ± 0.02 % Cr.....0.063 ± 0.006 % As.....0.010 ± 0.002 % Mn.....0.448 ± 0.014 % Mo.....0.011 ± 0.001 % Cu.....0.150 ± 0.004 % P.....0.019 ± 0.001 % Ni.....1.19 ± 0.02 % Mg.....0.051 ± 0.002 %	100 g
ECRM-B 482-2	Low alloyed cast iron - chips Year of issue: 1994 Certified values C _{total}2.599 ± 0.012 %* P.....0.0974±0.0015 %* Mo.....0.454 ± 0.004 %* Si.....1.815 ± 0.007 %* S.....0.0491±0.0015 %* Ni.....2.284 ± 0.009 %* Mn.....0.728 ± 0.005 %* Cr.....0.675 ± 0.006 %* Cu.....1.231 ± 0.008 %*	100 g
	* 95%-confidence interval	
ECRM-B 483-1	High duty iron - chips Year of issue: 1980 Certified values C _{total}2.46 ± 0.02 % Mn.....0.596 ± 0.011 % Cr.....0.039 ± 0.004 % Graphit(1.65) % P.....0.615 ± 0.018 % Sn.....0.130 ± 0.003 % Si.....1.75 ± 0.02 % S.....0.103 ± 0.003 %	100 g
	(Values in parenthesis are indicative values)	
ECRM-B 484-1	Whiteheart malleable iron - chips Year of issue: 1980 Certified values C _{total}3.20 ± 0.03 % Mn.....0.395 ± 0.005 % S.....0.230 ± 0.009 % Si.....0.717 ± 0.016 % P.....0.121 ± 0.005 % Cr.....0.155 ± 0.006 %	100 g
	(Values in parenthesis are indicative values)	
ECRM-B 486-1	Foundry iron - chips Year of issue: 1983 Certified values C _{total}2.21 ± 0.03 S.....0.023 ± 0.001 Sn.....0.074 ± 0.006 Si.....2.43 ± 0.02 Cr.....0.104 ± 0.006 V.....0.020 ± 0.003 Mn.....0.841 ± 0.017 Ni.....0.057 ± 0.004 P.....1.00 ± 0.04 Cu.....0.548 ± 0.013	100 g
ECRM-B 489-1	White cast iron - chips Year of issue: 1991 Certified values C _{total}2.860 ± 0.013 %* P.....0.815 ± 0.007 %* N(0.0056) % Si.....1.524 ± 0.006 %* S.....0.155 ± 0.002 %* Mn(0.510) Cu.....0.274 ± 0.003 %*	100 g
	(Values in parenthesis are indicative values) * 95%-confidence interval	
ECRM-F 485-2	Cast iron - chips Year of issue: 1994 Certified values C _{total}3.308 ± 0.015 % P0.610 ± 0.007 % Ni.....0.0742 ± 0.0010 % Si.....0.508 ± 0.005 % S.....0.165 ± 0.003 % Cu0.0303 ± 0.0006 % Mn0.190 ± 0.004 % Cr.....0.0845 ± 0.0009 % N0.0058 ± 0.0002 %	100 g
ECRM-F 485-3	Cast iron - chips Year of issue: 2006 Certified values C _{total}3.514 ± 0.009 %* S0.1488 ± 0.0024 %* N0.0081 ± 0.0002 %*	100 g
	* 95%-confidence interval	
ECRM-F 487-1	Cast iron - chips Year of issue: 1986 Certified values C _{total}3.27 ± 0.04 % S.....0.0007 ± 0.0002 % Co0.0088 ± 0.0007 % Si.....0.006 ± 0.001 % Cr.....0.063 ± 0.003 % N0.0042 ± 0.0004 % Mn0.094 ± 0.002 % Ni0.040 ± 0.002 % P.....0.0020 ± 0.0004 % As0.011 ± 0.001 %	100 g

Code	Product	Unit
ECRM-F 488-2	Low alloyed white cast iron - chips Year of issue: 1998 Certified values C _{total}3.965 ± 0.009 %* Si.....0.374 ± 0.005 %* Mn.....0.201 ± 0.002 %* P.....0.0111±0.0002 %* S.....0.1173±0.0011 %* Cr.....0.303 ± 0.004 %* Mo(0.0008) Ni.....0.1247±0.0009 %* Cu.....0.0256±0.0004 %* N.....0.0052±0.0003 %* Sn0.0013±0.0001 %* Ti.....0.0636±0.0009 %* V0.0545±0.0005 %* Te.....(0.0089) % Zn.....(0.0006) % (Values in parenthesis are indicative values) * 95%-confidence interval	100 g
ECRM-F 490-1	White cast iron (Spiegel) with 4.8% C and 10.8% Mn - chips Year of issue: 2003 Certified values C _{total}4.813 ± 0.011 %* Si.....(0.029) Mn.....10.83 ± 0.05 %* P.....0.0267 ± 0.0006 %* S.....0.0040 ± 0.0004 %* Cr.....0.0183 ± 0.0010 %* Ni.....(0.019) % Cu.....0.0088 ± 0.0002 %* N0.0030 ± 0.0002 %* Ti.....0.0035 ± 0.0003 %* V0.0152 ± 0.0007 %* (Values in parenthesis are indicative values) * 95%-confidence interval	100 g
ECRM-F 491-1	White cast iron, for C/S-determination - chips Year of issue: 2006 Certified values C _{total}3.616 ± 0.018 %* S.....0.0866 ± 0.0013 %* * 95%-confidence interval	100 g

Ferro alloys

ECRM-D 502-2	Ferromanganese, high-carbon - powder Year of issue: 2004 Certified values C6.94 ± 0.02 % Si.....(0.092) % Mn.....77.87 ± 0.11 % P0.148 ± 0.004 % S(0.0024) % Cr.....0.0265 ± 0.0006 % Ni.....0.0384 ± 0.0003 % B(0.0006) % Co.....(0.048) % Cu.....0.037 ± 0.0006 % N(0.017) % Ti.....0.0034 ± 0.0003 % Fe.....(14.6) % Pb0.0179 ± 0.0011 % (Values in parenthesis are indicative values)	100 g
ECRM-D 529-1	Ferro silicon, silicon content 91.1% - powder Year of issue: 1975 Certified values C0.10 ± 0.01 % Si.....91.11 ± 0.33 % Mn.....0.04 ± 0.005 % P0.013 ± 0.001 % Al.....0.86 ± 0.02 % Cu.....0.01 ± 0.001 % Ti.....0.09 ± 0.004 % Ca.....0.46 ± 0.04 % Fe.....6.15 ± 0.08 % Mg.....0.04 ± 0.006 %	100 g
ECRM-D 589-1	Ferro titanium with titanium content 68.4%, vanadium content 2.3%, aluminium content 5.3% and iron content 16.9% - powder Year of issue: 1991 Certified values C0.132 ± 0.008 % Si.....(0.41) % Mn.....0.151 ± 0.005 % P(0.0107) % S0.0152 ± 0.0011 % Cr.....0.506 ± 0.023 % Mo.....0.934 ± 0.017 % Ni.....0.663 ± 0.015 % Al.....5.34 ± 0.08 % Co.....0.115 ± 0.006 % Cu.....0.146 ± 0.006 % N.....0.64 ± 0.05 % Sn0.55 ± 0.03 % Ti.....68.4 ± 0.5 % V2.32 ± 0.07 % Zr.....0.866 ± 0.030 % Fe.....16.93 ± 0.17 % Zn.....(0.0103) % (Values in parenthesis are indicative values)	100 g
ECRM-D 591-1	Ferro vanadium with vanadium content 79.7%, aluminium content 3.2% and iron content 14.6% - powder Year of issue: 1996 Certified values C0.141 ± 0.004 % Si.....0.847 ± 0.012 % Mn.....0.307 ± 0.004 % P0.0299 ± 0.0017 % S0.0153 ± 0.0008 % Ni.....0.0141 ± 0.0014 % Al.....3.19 ± 0.05 % As0.0022 ± 0.0002 % B(0.0018) % Cu.....0.0596 ± 0.0016 % N.....(0.308) % Ti(0.044) % V79.72 ± 0.14 % Ca.....(0.0328) % Fe.....14.59 ± 0.10 % Mg.....(0.044) % O(0.516) % (Values in parenthesis are indicative values)	100 g

Iron and steel products

Code	Product	Unit
ECRM-B 554-1	Ferro chromium, low carbon with chromium content 74.6% - chips Year of issue: 1972 Certified values C 0.007 % P (0.018) % V.....(0.10) % Si.....0.53 % S.....(< 0.005) % Mn.....(0.09) % Cr.....74.6 % (Values in parenthesis are indicative values)	100 g
ECRM-B 555-1	Ferro tungsten, tungsten content 79.9% -powder Year of issue: 1975 Certified values C 0.025 % S.....(0.018) % W 79.9 % Si.....1.75 % Al 0.14 % P.....(0.02) % Sn 0.034 % (Values in parenthesis are indicative values)	100 g
ECRM-B 576-1	Ferro niobium, niobium content 40% - powder Year of issue: 1976 Certified values C 0.201 ± 0.004 % Nb 43.90 ± 0.22 % Ta.....0.306 ± 0.023 % Si.....1.79 ± 0.07 % Sn 0.195 ± 0.011 % Al 2.53 ± 0.05 % Ti 1.32 ± 0.04 %	100 g
ECRM-B 577-1	Ferro vanadium with vanadium content 50.2% and aluminium content 0.4% - powder Year of issue: 1978 Certified values C 0.089 ± 0.003 % P 0.035 ± 0.003 % Al.....0.414 ± 0.018 % Si.....1.79 ± 0.03 % S 0.034 ± 0.002 % Cu 0.054 ± 0.003 % Mn 0.158 ± 0.007 % Ni 0.053 ± 0.005 % V.....50.16 ± 0.13 %	100 g
ECRM-B 578-1	Ferro molybdenum, molybdenum content 72.2% - powder Year of issue: 1978 Certified values C 0.016 ± 0.002 % P 0.024 ± 0.002 % Mo 72.23 ± 0.16 % Si.....0.208 ± 0.010 % S 0.065 ± 0.003 % Cu 0.136 ± 0.010 %	100 g
ECRM-B 579-1	Ferro niobium, niobium content 60% - powder Year of issue: 1981 Certified values C 0.037 ± 0.002 % Al 1.86 ± 0.06 % Ti.....0.567 ± 0.053 % Si.....1.03 ± 0.05 % Co 0.0051 ± 0.0010 % Ta.....3.85 ± 0.15 % P.....0.064 ± 0.004 % Nb 62.87 ± 0.17 % S.....0.021 ± 0.002 % Sn 0.344 ± 0.016 %	100 g
ECRM-B 580-1	Ferro chromium, low carbon with chromium content 72,2% - chips Year of issue: 1981 Certified values C 0.019 ± 0.002 % Cr 72.18 ± 0.11 % V.....0.083 ± 0.006 % Si.....0.306 ± 0.011 % Co 0.047 ± 0.004 % P.....0.011 ± 0.001 % N 0.035 ± 0.002 %	100 g
ECRM-B 583-1	Ferro manganese, low carbon with manganese content 86,4% - powder Year of issue: 1982 Certified values C 0.333 ± 0.008 % Mn 86.42 ± 0.15 % S.....(0.007) % Si.....0.396 ± 0.012 % P 0.146 ± 0.006 % (Values in parenthesis are indicative values)	100 g
ECRM-B 584-1	Ferro titanium with titanium content 37.2% and aluminium content 7.2% - powder Year of issue: 1983 Certified values C 0.044 ± 0.004 % P 0.032 ± 0.003 % Ti.....37.17 ± 0.26 % Si.....1.80 ± 0.07 % S 0.030 ± 0.002 % Mn 1.13 ± 0.02 % Al 7.19 ± 0.23 % (Values in parenthesis are indicative values)	100 g
ECRM-B 585-2	High-Carbon Ferro-Chromium - powder	100 g
ECRM-B 587-1	Ferro boron, boron content 18.7% - powder Year of issue: 1986 Certified values C 0.738 ± 0.016 % Cr.....(0.104) % Ti.....(0.039) % Si.....(0.129) % Mo.....(0.005) % V.....(0.004) % Mn.....0.272 ± 0.008 % Al 0.047 ± 0.003 % Ca (0.048) % P.....(0.020) % B 18.7 ± 0.3 % S.....(0.0010) % Co (0.010) % (Values in parenthesis are indicative values)	100 g

Iron and steel products

Code	Product	Unit
ECRM-B 590-1	Ferro tungsten, tungsten content 79.5% - powder Year of issue: 1991 Certified values C 0.025 ± 0.007 % Mo 0.101 ± 0.003 % W 79.5 ± 0.11 % Si 1.05 ± 0.03 % Cu 0.048 ± 0.001 % Mn 0.136 ± 0.004 % Sn 0.045 ± 0.002 %	100 g
ECRM-F 503-1	Ferro manganese, low carbon with manganese content 80.8% - powder Year of issue: 1962 Certified values C 0.700 ± 0.002 % Mn 80.80 ± 0.10 % S (0.009) % Si 0.865 ± 0.006 % P 0.069 ± 0.002 % (Values in parenthesis are indicative values)	100 g
ECRM-F 507-1	Ferro chromium, high carbon with chromium content 70.3% - powder Year of issue: 1965 Certified values C 5.40 ± 0.10 % Mn 0.27 ± 0.03 % Cr 70.30 ± 0.10 % Si 1.20 ± 0.03 % P 0.017 ± 0.003 % N 0.049 ± 0.003 %	100 g
ECRM-F 509-1	Ferro chromium, low carbon with chromium content 72.9% - chips Year of issue: 1967 Certified values C 0.012 ± 0.001 % P (0.019) % N (0.032) % Si 0.230 ± 0.008 % Cr 72.85 ± 0.15 % (Values in parenthesis are indicative values)	100 g
ECRM-F 510-1	Ferro titanium with titanium content 26.9% - powder Year of issue: 1967 Certified values C 0.058 ± 0.003 % P (0.035) % Ti 26.95 ± 0.05 % Si 4.65 ± 0.08 % Al (4.9) % (Values in parenthesis are indicative values)	100 g
ECRM-F 511-1	Ferro vanadium with vanadium content 80.7% - powder Year of issue: 1967 Certified values C 0.049 ± 0.002 % P (0.016) % V 80.7 ± 0.2 % Si 0.341 ± 0.008 % S 0.018 ± 0.002 % (Values in parenthesis are indicative values)	100 g
ECRM-F 582-2	Ferro silicon, silicon content 75.2% - powder Year of issue: 1990 Certified values C 0.150 ± 0.013 % S (0.003) % Ca 0.405 ± 0.012 % Si 75.22 ± 0.20 % Cr 0.074 ± 0.004 % Fe 21.42 ± 0.22 % Mn 0.230 ± 0.006 % Al 1.154 ± 0.030 % P 0.018 ± 0.0014 % Ti 0.225 ± 0.011 % (Values in parenthesis are indicative values)	100 g
ECRM-F 586-1	Silico manganese with manganese content 62.5% and silicon content 34% - powder Year of issue: 1985 Certified values C 0.025 ± 0.003 % Cr 0.044 ± 0.002 % Ca 0.039 ± 0.003 % Si 34.0 ± 0.2 % Al 0.022 ± 0.003 % Fe 2.89 ± 0.05 % Mn 62.5 ± 0.2 % Co 0.007 ± 0.001 % P 0.041 ± 0.003 % V 0.041 ± 0.003 %	100 g

Code	Product	Unit
Ores, iron oxide		
ECRM-D 627-2	Minette CaO-rich with 31.8% iron content - powder Year of issue: 1966 Certified values Fe _{total} 31.77 ± 0.12 % MgO 1.57 ± 0.08 % S 0.020 ± 0.001 % SiO ₂ 9.24 ± 0.08 % Mn 0.250 ± 0.012 % Cr 0.018 ± 0.003 % Al ₂ O ₃ 4.49 ± 0.12 % P 0.661 ± 0.014 % Cu (0.002) % CaO 15.67 ± 0.21 % S 0.114 ± 0.009 % TiO ₂ 0.225 ± 0.014 % (Values in parenthesis are indicative values)	100 g

Iron and steel products

Code	Product	Unit
ECRM-D 629-1	Minette CaO-rich with 36.2% iron content - powder Year of issue:1966 Certified values Fe _{total} 36.21 ± 0.13 % MgO 1.64 ± 0.08 % As 0.023 ± 0.001 % SiO ₂ 19.25 ± 0.14 % Mn 0.390 ± 0.012 % Cr 0.016 ± 0.002 % Al ₂ O ₃ 4.07 ± 0.08 % P 0.696 ± 0.013 % Cu (0.001) % CaO 5.63 ± 0.08 % S 0.063 ± 0.006 % TiO ₂ 0.216 ± 0.013 % (Values in parenthesis are indicative values)	100 g
ECRM-D 630-1	Iron ore, Bomi hill concentrate with 65.6% iron content - powder Year of issue:1969 Certified values Fe _{total} 65.63 ± 0.17 % CaO 0.10 ± 0.017 % P 0.043 ± 0.003 % SiO ₂ 5.88 ± 0.07 % MgO 0.47 ± 0.046 % S 0.032 ± 0.004 % Al ₂ O ₃ 0.88 ± 0.038 % Mn 0.060 ± 0.005 % TiO ₂ 0.066 ± 0.013 %	100 g
ECRM-D 631-1	Venezuela iron ore with 61.1% iron content - powder Year of issue:1961 Certified values Fetotal 61.09 ± 0.09 % MgO 0.54 ± 0.059 % Na ₂ O (0.04) % SiO ₂ 3.20 ± 0.06 % Mn 0.044 ± 0.006 % K ₂ O (0.04) % Al ₂ O ₃ 1.06 ± 0.05 % P 0.114 ± 0.005 % TiO ₂ 0.109 ± 0.006 % CaO 0.75 ± 0.038 % S 0.033 ± 0.006 % (Values in parenthesis are indicative values)	100 g
ECRM-D 633-1	Manganese ore, Manganese Content 47.9% - powder Year of issue:1967 Certified values Fe _{total} 1.64 ± 0.04 % MgO 0.58 ± 0.10 % As (0.0040) % SiO ₂ 10.39 ± 0.15 % Mn 47.85 ± 0.21 % BaO 1.13 ± 0.08 % Al ₂ O ₃ 1.64 ± 0.12 % P 0.170 ± 0.007 % TiO ₂ 0.079 ± 0.009 % CaO 2.02 ± 0.12 % S 0.227 ± 0.009 % (Values in parenthesis are indicative values)	100 g
ECRM-D 678-1	Kiruna D Iron Ore with 60.8% iron content - powder Year of issue:1975 Certified values Fe _{total} 60.75 ± 0.07 % Mn 0.08 ± 0.01 % K 0.11 ± 0.01 % Si 1.73 ± 0.04 % P 1.61 ± 0.04 % K ₂ O 0.13 % Al 0.28 ± 0.03 % S 0.021 ± 0.002 % F 0.29 ± 0.02 % Ca 3.92 ± 0.09 % Na 0.11 ± 0.01 % Ti 0.13 ± 0.01 % Mg 0.57 ± 0.02 % Na ₂ O 0.15 % V 0.12 ± 0.01 %	100 g
ECRM-D 680-1	Iron ore, purple ore with 60% iron content - powder Year of issue:1977 Certified values Fe _{total} 59.98 ± 0.08 % MgO 0.23 ± 0.02 % Cu 0.063 ± 0.003 % Si 4.20 ± 0.02 % Mn 0.025 ± 0.002 % Ni 0.007 ± 0.001 % SiO ₂ 8.98 ± 0.04 % P 0.018 ± 0.002 % Pb 0.317 ± 0.008 % Al 0.66 ± 0.02 % S 0.544 ± 0.017 % Ti 0.045 ± 0.003 % Al ₂ O ₃ 1.23 ± 0.04 % Na 0.128 ± 0.004 % TiO ₂ 0.08 ± 0.005 % Ca 0.45 ± 0.02 % K 0.078 ± 0.003 % Zn 0.165 ± 0.004 % CaO 0.63 ± 0.03 % As 0.057 ± 0.003 % Mg 0.14 ± 0.01 % Cr 0.005 ± 0.001 % Sn 0.0025 ± 0.0002 %	100 g
ECRM-D 686-1	Iron oxide recovered from pickling bath (Ruthner) - powder Year of issue:2002 Certified values Fe _{total} 69.44 ± 0.11 % P 0.0078 ± 0.0001 % Ti 0.0014 ± 0.0001 % Si 0.0083 ± 0.0005 % Na 0.0058 ± 0.0005 % Zn 0.0004 ± 0.0001 % Al 0.0407 ± 0.0012 % K 0.0024 ± 0.0004 % Cl 0.095 ± 0.006 % Ca 0.0097 ± 0.0007 % Cr 0.0182 ± 0.0006 % Co 0.0019 ± 0.0001 % Mg 0.0027 ± 0.0002 % Cu 0.0038 ± 0.0003 % Mo 0.0007 ± 0.0001 % Mn 0.231 ± 0.004 % Ni 0.0127 ± 0.0004 % Sn 0.0025 ± 0.0002 %	100 g
ECRM-B 651-1	Lincolnshire iron ore - powder Year of issue: 1975 Certified values Fe _{total} 23.85 % Mg 1.04 % Na ₂ O 0.07 % Si 3.46 % MgO 1.73 % K 0.27 % SiO ₂ 7.40 % Mn 0.97 % K ₂ O 0.32 % Al 2.25 % MnO 1.25 % Ti 0.10 % Al ₂ O ₃ 4.26 % P 0.35 % TiO ₂ 0.16 % Ca 16.2 % S 0.40 % CaO 22.6 % Na 0.05 %	100 g

Iron and steel products

Code	Product	Unit																																																												
ECRM-B 676-1	Iron ore sinter - powder Year of issue: 1975 Certified values <table> <tr><td>Fe_{total}</td><td>39.76 ± 0.08 %</td><td>Mg</td><td>1.16 ± 0.04 %</td><td>Na₂O</td><td>0.128 %</td></tr> <tr><td>Si</td><td>6.40 ± 0.05 %</td><td>MgO</td><td>1.92 %</td><td>K</td><td>0.43 ± 0.02 %</td></tr> <tr><td>SiO₂</td><td>13.69 %</td><td>Mn</td><td>0.83 ± 0.01 %</td><td>K₂O</td><td>0.52 %</td></tr> <tr><td>Al</td><td>3.40 ± 0.07 %</td><td>P</td><td>0.59 ± 0.02 %</td><td>F</td><td>0.10 ± 0.01 %</td></tr> <tr><td>Al₂O₃</td><td>6.42 %</td><td>P₂O₅</td><td>1.35 %</td><td>Ti</td><td>0.19 ± 0.01 %</td></tr> <tr><td>Ca</td><td>12.78 ± 0.13 %</td><td>S</td><td>0.12 ± 0.006 %</td><td>TiO₂</td><td>0.32 %</td></tr> <tr><td>CaO</td><td>17.88 %</td><td>Na</td><td>0.095 ± 0.004 %</td><td>V</td><td>0.070 ± 0.004 %</td></tr> </table>	Fe _{total}	39.76 ± 0.08 %	Mg	1.16 ± 0.04 %	Na ₂ O	0.128 %	Si	6.40 ± 0.05 %	MgO	1.92 %	K	0.43 ± 0.02 %	SiO ₂	13.69 %	Mn	0.83 ± 0.01 %	K ₂ O	0.52 %	Al	3.40 ± 0.07 %	P	0.59 ± 0.02 %	F	0.10 ± 0.01 %	Al ₂ O ₃	6.42 %	P ₂ O ₅	1.35 %	Ti	0.19 ± 0.01 %	Ca	12.78 ± 0.13 %	S	0.12 ± 0.006 %	TiO ₂	0.32 %	CaO	17.88 %	Na	0.095 ± 0.004 %	V	0.070 ± 0.004 %	100 g																		
Fe _{total}	39.76 ± 0.08 %	Mg	1.16 ± 0.04 %	Na ₂ O	0.128 %																																																									
Si	6.40 ± 0.05 %	MgO	1.92 %	K	0.43 ± 0.02 %																																																									
SiO ₂	13.69 %	Mn	0.83 ± 0.01 %	K ₂ O	0.52 %																																																									
Al	3.40 ± 0.07 %	P	0.59 ± 0.02 %	F	0.10 ± 0.01 %																																																									
Al ₂ O ₃	6.42 %	P ₂ O ₅	1.35 %	Ti	0.19 ± 0.01 %																																																									
Ca	12.78 ± 0.13 %	S	0.12 ± 0.006 %	TiO ₂	0.32 %																																																									
CaO	17.88 %	Na	0.095 ± 0.004 %	V	0.070 ± 0.004 %																																																									
ECRM-B 681-1	Iron ore - powder Year of issue: 1976 Certified values <table> <tr><td>Fe_{total}</td><td>33.21 ± 0.13 %</td><td>MnO</td><td>0.29 %</td><td>Cr</td><td>0.041 ± 0.002 %</td></tr> <tr><td>Si</td><td>8.32 ± 0.08 %</td><td>P</td><td>0.88 ± 0.02 %</td><td>Cr₂O₃</td><td>0.060 %</td></tr> <tr><td>SiO₂</td><td>17.79 %</td><td>P₂O₅</td><td>2.01 %</td><td>F</td><td>0.19 ± 0.01 %</td></tr> <tr><td>Al</td><td>5.62 ± 0.10 %</td><td>S</td><td>0.103 ± 0.006 %</td><td>Ni</td><td>0.016 ± 0.001 %</td></tr> <tr><td>Al₂O₃</td><td>10.61 %</td><td>Na</td><td>0.068 ± 0.008 %</td><td>Pb</td><td>(0.0072) %</td></tr> <tr><td>Ca</td><td>2.80 ± 0.10 %</td><td>Na₂O</td><td>0.092 %</td><td>Ti</td><td>0.29 ± 0.01 %</td></tr> <tr><td>CaO</td><td>3.92 %</td><td>K</td><td>0.49 ± 0.01 %</td><td>TiO₂</td><td>0.48 %</td></tr> <tr><td>Mg</td><td>0.89 ± 0.02 %</td><td>K₂O</td><td>0.59 %</td><td>V</td><td>0.077 ± 0.004 %</td></tr> <tr><td>MgO</td><td>1.47 %</td><td>As</td><td>(0.0109) %</td><td>V₂O₅</td><td>0.138 %</td></tr> <tr><td>Mn</td><td>0.22 ± 0.01 %</td><td>C</td><td>1.80 ± 0.04 %</td><td></td><td></td></tr> </table>	Fe _{total}	33.21 ± 0.13 %	MnO	0.29 %	Cr	0.041 ± 0.002 %	Si	8.32 ± 0.08 %	P	0.88 ± 0.02 %	Cr ₂ O ₃	0.060 %	SiO ₂	17.79 %	P ₂ O ₅	2.01 %	F	0.19 ± 0.01 %	Al	5.62 ± 0.10 %	S	0.103 ± 0.006 %	Ni	0.016 ± 0.001 %	Al ₂ O ₃	10.61 %	Na	0.068 ± 0.008 %	Pb	(0.0072) %	Ca	2.80 ± 0.10 %	Na ₂ O	0.092 %	Ti	0.29 ± 0.01 %	CaO	3.92 %	K	0.49 ± 0.01 %	TiO ₂	0.48 %	Mg	0.89 ± 0.02 %	K ₂ O	0.59 %	V	0.077 ± 0.004 %	MgO	1.47 %	As	(0.0109) %	V ₂ O ₅	0.138 %	Mn	0.22 ± 0.01 %	C	1.80 ± 0.04 %			100 g
Fe _{total}	33.21 ± 0.13 %	MnO	0.29 %	Cr	0.041 ± 0.002 %																																																									
Si	8.32 ± 0.08 %	P	0.88 ± 0.02 %	Cr ₂ O ₃	0.060 %																																																									
SiO ₂	17.79 %	P ₂ O ₅	2.01 %	F	0.19 ± 0.01 %																																																									
Al	5.62 ± 0.10 %	S	0.103 ± 0.006 %	Ni	0.016 ± 0.001 %																																																									
Al ₂ O ₃	10.61 %	Na	0.068 ± 0.008 %	Pb	(0.0072) %																																																									
Ca	2.80 ± 0.10 %	Na ₂ O	0.092 %	Ti	0.29 ± 0.01 %																																																									
CaO	3.92 %	K	0.49 ± 0.01 %	TiO ₂	0.48 %																																																									
Mg	0.89 ± 0.02 %	K ₂ O	0.59 %	V	0.077 ± 0.004 %																																																									
MgO	1.47 %	As	(0.0109) %	V ₂ O ₅	0.138 %																																																									
Mn	0.22 ± 0.01 %	C	1.80 ± 0.04 %																																																											
	(Values in parenthesis are indicative values)																																																													
ECRM-B 683-1	Iron ore sinter - powder Year of issue: 1982 Certified values <table> <tr><td>Fe_{total}</td><td>56.06 ± 0.12 %</td><td>Mn</td><td>0.462 ± 0.007 %</td><td>Cr</td><td>0.018 ± 0.002 %</td></tr> <tr><td>Si</td><td>3.38 ± 0.03 %</td><td>MnO</td><td>0.596 %</td><td>Cr₂O₃</td><td>0.026 %</td></tr> <tr><td>SiO₂</td><td>7.23 %</td><td>P</td><td>0.148 ± 0.004 %</td><td>F</td><td>0.020 ± 0.002 %</td></tr> <tr><td>Al</td><td>1.30 ± 0.05 %</td><td>P₂O₅</td><td>0.339 %</td><td>Ti</td><td>0.097 ± 0.006 %</td></tr> <tr><td>Al₂O₃</td><td>2.46 %</td><td>S</td><td>(0.013) %</td><td>TiO₂</td><td>0.162 %</td></tr> <tr><td>Ca</td><td>5.70 ± 0.08 %</td><td>Na</td><td>0.045 ± 0.007 %</td><td>V</td><td>0.026 ± 0.004 %</td></tr> <tr><td>CaO</td><td>7.97 %</td><td>Na₂O</td><td>0.061 %</td><td>V₂O₅</td><td>0.046 %</td></tr> <tr><td>Mg</td><td>1.04 ± 0.02 %</td><td>K</td><td>0.148 ± 0.012 %</td><td>Zn</td><td>0.010 ± 0.001 %</td></tr> <tr><td>MgO</td><td>1.72 %</td><td>K₂O</td><td>0.178 %</td><td></td><td></td></tr> </table>	Fe _{total}	56.06 ± 0.12 %	Mn	0.462 ± 0.007 %	Cr	0.018 ± 0.002 %	Si	3.38 ± 0.03 %	MnO	0.596 %	Cr ₂ O ₃	0.026 %	SiO ₂	7.23 %	P	0.148 ± 0.004 %	F	0.020 ± 0.002 %	Al	1.30 ± 0.05 %	P ₂ O ₅	0.339 %	Ti	0.097 ± 0.006 %	Al ₂ O ₃	2.46 %	S	(0.013) %	TiO ₂	0.162 %	Ca	5.70 ± 0.08 %	Na	0.045 ± 0.007 %	V	0.026 ± 0.004 %	CaO	7.97 %	Na ₂ O	0.061 %	V ₂ O ₅	0.046 %	Mg	1.04 ± 0.02 %	K	0.148 ± 0.012 %	Zn	0.010 ± 0.001 %	MgO	1.72 %	K ₂ O	0.178 %			100 g						
Fe _{total}	56.06 ± 0.12 %	Mn	0.462 ± 0.007 %	Cr	0.018 ± 0.002 %																																																									
Si	3.38 ± 0.03 %	MnO	0.596 %	Cr ₂ O ₃	0.026 %																																																									
SiO ₂	7.23 %	P	0.148 ± 0.004 %	F	0.020 ± 0.002 %																																																									
Al	1.30 ± 0.05 %	P ₂ O ₅	0.339 %	Ti	0.097 ± 0.006 %																																																									
Al ₂ O ₃	2.46 %	S	(0.013) %	TiO ₂	0.162 %																																																									
Ca	5.70 ± 0.08 %	Na	0.045 ± 0.007 %	V	0.026 ± 0.004 %																																																									
CaO	7.97 %	Na ₂ O	0.061 %	V ₂ O ₅	0.046 %																																																									
Mg	1.04 ± 0.02 %	K	0.148 ± 0.012 %	Zn	0.010 ± 0.001 %																																																									
MgO	1.72 %	K ₂ O	0.178 %																																																											
ECRM-F 601-1	Iron ore 36.8% iron content - powder Year of issue: 1968 Certified values <table> <tr><td>Fe_{total}</td><td>36.76 ± 0.17 %</td><td>Ca</td><td>4.05 ± 0.15 %</td><td>S</td><td>0.065 ± 0.008 %</td></tr> <tr><td>Fe(II)</td><td>(8.8) %</td><td>Mg</td><td>1.21 ± 0.08 %</td><td>Ti</td><td>0.114 ± 0.005 %</td></tr> <tr><td>Si</td><td>8.95 ± 0.07 %</td><td>Mn</td><td>0.370 ± 0.013 %</td><td></td><td></td></tr> <tr><td>Al</td><td>2.33 ± 0.25 %</td><td>P</td><td>0.590 ± 0.016 %</td><td></td><td></td></tr> </table>	Fe _{total}	36.76 ± 0.17 %	Ca	4.05 ± 0.15 %	S	0.065 ± 0.008 %	Fe(II)	(8.8) %	Mg	1.21 ± 0.08 %	Ti	0.114 ± 0.005 %	Si	8.95 ± 0.07 %	Mn	0.370 ± 0.013 %			Al	2.33 ± 0.25 %	P	0.590 ± 0.016 %			100 g																																				
Fe _{total}	36.76 ± 0.17 %	Ca	4.05 ± 0.15 %	S	0.065 ± 0.008 %																																																									
Fe(II)	(8.8) %	Mg	1.21 ± 0.08 %	Ti	0.114 ± 0.005 %																																																									
Si	8.95 ± 0.07 %	Mn	0.370 ± 0.013 %																																																											
Al	2.33 ± 0.25 %	P	0.590 ± 0.016 %																																																											
	(Values in parenthesis are indicative values)																																																													
ECRM-F 603-1	Iron ore 53.7% iron content - powder Year of issue: 1968 Certified values <table> <tr><td>Fe_{total}</td><td>53.65 ± 0.15 %</td><td>Ca</td><td>(0.91) %</td><td>S</td><td>0.097 ± 0.006 %</td></tr> <tr><td>Fe(II)</td><td>(0.3) %</td><td>Mg</td><td>(0.2) %</td><td>Ti</td><td>0.137 ± 0.008 %</td></tr> <tr><td>Si</td><td>1.28 ± 0.03 %</td><td>Mn</td><td>0.440 ± 0.013 %</td><td></td><td></td></tr> <tr><td>Al</td><td>4.20 ± 0.16 %</td><td>P</td><td>0.084 ± 0.007 %</td><td></td><td></td></tr> </table>	Fe _{total}	53.65 ± 0.15 %	Ca	(0.91) %	S	0.097 ± 0.006 %	Fe(II)	(0.3) %	Mg	(0.2) %	Ti	0.137 ± 0.008 %	Si	1.28 ± 0.03 %	Mn	0.440 ± 0.013 %			Al	4.20 ± 0.16 %	P	0.084 ± 0.007 %			100 g																																				
Fe _{total}	53.65 ± 0.15 %	Ca	(0.91) %	S	0.097 ± 0.006 %																																																									
Fe(II)	(0.3) %	Mg	(0.2) %	Ti	0.137 ± 0.008 %																																																									
Si	1.28 ± 0.03 %	Mn	0.440 ± 0.013 %																																																											
Al	4.20 ± 0.16 %	P	0.084 ± 0.007 %																																																											
	(Values in parenthesis are indicative values)																																																													
ECRM-F 604-1	Iron ore 65.7% iron content - powder Year of issue: 1968 Certified values <table> <tr><td>Fe_{total}</td><td>65.69 ± 0.17 %</td><td>Ca</td><td>(0.13) %</td><td>S</td><td>0.015 ± 0.002 %</td></tr> <tr><td>Fe(II)</td><td>(0.5) %</td><td>Mg</td><td>(0.06) %</td><td>Ti</td><td>0.060 ± 0.006 %</td></tr> <tr><td>Si</td><td>1.27 ± 0.04 %</td><td>Mn</td><td>0.092 ± 0.008 %</td><td></td><td></td></tr> <tr><td>Al</td><td>0.93 ± 0.07 %</td><td>P</td><td>0.053 ± 0.010 %</td><td></td><td></td></tr> </table>	Fe _{total}	65.69 ± 0.17 %	Ca	(0.13) %	S	0.015 ± 0.002 %	Fe(II)	(0.5) %	Mg	(0.06) %	Ti	0.060 ± 0.006 %	Si	1.27 ± 0.04 %	Mn	0.092 ± 0.008 %			Al	0.93 ± 0.07 %	P	0.053 ± 0.010 %			100 g																																				
Fe _{total}	65.69 ± 0.17 %	Ca	(0.13) %	S	0.015 ± 0.002 %																																																									
Fe(II)	(0.5) %	Mg	(0.06) %	Ti	0.060 ± 0.006 %																																																									
Si	1.27 ± 0.04 %	Mn	0.092 ± 0.008 %																																																											
Al	0.93 ± 0.07 %	P	0.053 ± 0.010 %																																																											
	(Values in parenthesis are indicative values)																																																													
ECRM-F 606-1	Iron ore 59.7% iron content - powder Year of issue: 1970 Certified values <table> <tr><td>Fe_{total}</td><td>59.66 ± 0.14 %</td><td>Ca</td><td>1.04 ± 0.12 %</td><td>S</td><td>0.033 ± 0.004 %</td></tr> <tr><td>Fe(II)</td><td>(0.1) %</td><td>Mg</td><td>0.32 ± 0.05 %</td><td>Ti</td><td>0.019 ± 0.004 %</td></tr> <tr><td>Si</td><td>1.04 ± 0.03 %</td><td>Mn</td><td>2.59 ± 0.06 %</td><td></td><td></td></tr> <tr><td>Al</td><td>0.34 ± 0.06 %</td><td>P</td><td>0.026 ± 0.006 %</td><td></td><td></td></tr> </table>	Fe _{total}	59.66 ± 0.14 %	Ca	1.04 ± 0.12 %	S	0.033 ± 0.004 %	Fe(II)	(0.1) %	Mg	0.32 ± 0.05 %	Ti	0.019 ± 0.004 %	Si	1.04 ± 0.03 %	Mn	2.59 ± 0.06 %			Al	0.34 ± 0.06 %	P	0.026 ± 0.006 %			100 g																																				
Fe _{total}	59.66 ± 0.14 %	Ca	1.04 ± 0.12 %	S	0.033 ± 0.004 %																																																									
Fe(II)	(0.1) %	Mg	0.32 ± 0.05 %	Ti	0.019 ± 0.004 %																																																									
Si	1.04 ± 0.03 %	Mn	2.59 ± 0.06 %																																																											
Al	0.34 ± 0.06 %	P	0.026 ± 0.006 %																																																											
	(Values in parenthesis are indicative values)																																																													

Iron and steel products

Code	Product	Unit
ECRM-F 607-1	Iron ore 30.9% iron content Year of issue: 1970 Certified values Fe _{total} 30.89 ± 0.17 % Ca 13.74 ± 0.17 % S 0.050 ± 0.007 % Fe(II) (5.95) % Mg 0.77 ± 0.08 % Ti 0.123 ± 0.007 % Si 3.07 ± 0.07 % Mn 0.254 ± 0.011 % Al 2.48 ± 0.11 % P 0.529 ± 0.018 % (Values in parenthesis are indicative values)	100 g
ECRM-F 608-1	Ferriferous marl - powder Year of issue: 1970 Certified values Fe _{total} 4.00 ± 0.15 % Ca 6.22 ± 0.18 % S 0.455 ± 0.029 % Fe(II) (1.85) % Mg 0.81 ± 0.09 % Ti 0.428 ± 0.018 % Si 28.23 ± 0.27 % Mn 0.044 ± 0.007 % Al 5.26 ± 0.12 % P 0.053 ± 0.006 % (Values in parenthesis are indicative values)	100 g
ECRM-F 609-1	Iron ore 30.5% iron content - powder Year of issue: 1970 Certified values Fe _{total} 30.52 ± 0.16 % Ca 6.87 ± 0.12 % S 1.000 ± 0.057 % Fe(II) (15.65) % Mg 2.00 ± 0.09 % Ti 0.118 ± 0.007 % Si 7.83 ± 0.08 % Mn 0.472 ± 0.031 % Al 2.26 ± 0.11 % P 0.608 ± 0.043 % (Values in parenthesis are indicative values)	100 g
ECRM-F 610-1	Laterite, 47.5% iron content - powder Year of issue: 1971 Certified values Fe _{total} 47.46 ± 0.35 % Mg 1.86 ± 0.11 % Co 0.075 ± 0.006 % Si 3.16 ± 0.04 % Mn 0.581 ± 0.016 % Cr 1.84 ± 0.09 % Al 1.96 ± 0.10 % P 0.007 ± 0.002 % Ni 1.48 ± 0.07 % Ca (0.1) % S 0.189 ± 0.024 % Ti 0.015 ± 0.005 % (Values in parenthesis are indicative values)	100 g
ECRM-F 611-1	Iron ore sinter - powder Year of issue: 1971 Certified values Fe _{total} 62.22 ± 0.16 % Ca 2.85 ± 0.14 % S (0.008) % Fe(II) (13.84) % Mg 0.32 ± 0.04 % Ti 0.035 ± 0.005 % Si 2.07 ± 0.04 % Mn 1.97 ± 0.06 % Al 0.69 ± 0.05 % P 0.030 ± 0.008 % (Values in parenthesis are indicative values)	100 g
ECRM-F 612-1	Iron ore sinter - powder Year of issue: 1971 Certified values Fe _{total} 42.43 ± 0.13 % Ca 12.06 ± 0.17 % S 0.053 ± 0.007 % Fe(II) (9.19) % Mg 1.20 ± 0.07 % Ti 0.151 ± 0.006 % Si 5.94 ± 0.06 % Mn 0.363 ± 0.014 % Al 3.00 ± 0.13 % P 0.885 ± 0.027 % (Values in parenthesis are indicative values)	100 g
ECRM-F 677-1	Iron ore 51.5% iron content - powder Year of issue: 1975 Certified values Fe _{total} 51.54 ± 0.07 % P 0.017R ± 0.0007 % Cu (0.0012) % Si 11.78 ± 0.06 % S (0.005) % Ni (0.0015) % Al 0.32 ± 0.02 % Na 0.007 ± 0.002 % Pb (0.003) % Ca 0.038 ± 0.006 % K 0.008 ± 0.002 % Ti 0.013 ± 0.001 % Mg 0.012 ± 0.002 % Co (0.0006) % Zn (0.002) % Mn 0.016 ± 0.002 % Cr (0.002) % (Values in parenthesis are indicative values)	100 g

Iron and steel products

Code	Product	Unit																																				
ECRM-F 679-1	Iron ore 24.2% iron content - powder Year of issue: 1977 Certified values <table> <tr><td>Fe_{total}</td><td>24.20 ± 0.08 %</td><td>MgO..... 1.15 %</td><td>C 5.80 ± 0.07 %</td></tr> <tr><td>Fe(II)</td><td>(5.2) %</td><td>Mn 0.295 ± 0.005 %</td><td>Cr 0.012 ± 0.001 %</td></tr> <tr><td>Si.....</td><td>3.43 ± 0.04 %</td><td>P 0.557 ± 0.009 %</td><td>Ni (0.0095) %</td></tr> <tr><td>SiO₂.....</td><td>7.34 %</td><td>P₂O₅..... 1.275 %</td><td>Pb (0.0062) %</td></tr> <tr><td>Al.....</td><td>1.99 ± 0.05 %</td><td>S..... 0.099 ± 0.005 %</td><td>Ti 0.106 ± 0.005 %</td></tr> <tr><td>Al₂O₃.....</td><td>3.76 %</td><td>Na..... 0.054 ± 0.004 %</td><td>TiO₂..... 0.177 %</td></tr> <tr><td>Ca.....</td><td>18.07 ± 0.14 %</td><td>Na₂O..... 0.072 %</td><td>V 0.035 ± 0.002 %</td></tr> <tr><td>CaO.....</td><td>25.30 %</td><td>K..... 0.157 ± 0.007 %</td><td>Zn..... 0.021 ± 0.002 %</td></tr> <tr><td>Mg.....</td><td>0.70 ± 0.02 %</td><td>K₂O..... 0.189 %</td><td></td></tr> </table> (Values in parenthesis are indicative values)	Fe _{total}	24.20 ± 0.08 %	MgO..... 1.15 %	C 5.80 ± 0.07 %	Fe(II)	(5.2) %	Mn 0.295 ± 0.005 %	Cr 0.012 ± 0.001 %	Si.....	3.43 ± 0.04 %	P 0.557 ± 0.009 %	Ni (0.0095) %	SiO ₂	7.34 %	P ₂ O ₅ 1.275 %	Pb (0.0062) %	Al.....	1.99 ± 0.05 %	S..... 0.099 ± 0.005 %	Ti 0.106 ± 0.005 %	Al ₂ O ₃	3.76 %	Na..... 0.054 ± 0.004 %	TiO ₂ 0.177 %	Ca.....	18.07 ± 0.14 %	Na ₂ O..... 0.072 %	V 0.035 ± 0.002 %	CaO.....	25.30 %	K..... 0.157 ± 0.007 %	Zn..... 0.021 ± 0.002 %	Mg.....	0.70 ± 0.02 %	K ₂ O..... 0.189 %		100 g
Fe _{total}	24.20 ± 0.08 %	MgO..... 1.15 %	C 5.80 ± 0.07 %																																			
Fe(II)	(5.2) %	Mn 0.295 ± 0.005 %	Cr 0.012 ± 0.001 %																																			
Si.....	3.43 ± 0.04 %	P 0.557 ± 0.009 %	Ni (0.0095) %																																			
SiO ₂	7.34 %	P ₂ O ₅ 1.275 %	Pb (0.0062) %																																			
Al.....	1.99 ± 0.05 %	S..... 0.099 ± 0.005 %	Ti 0.106 ± 0.005 %																																			
Al ₂ O ₃	3.76 %	Na..... 0.054 ± 0.004 %	TiO ₂ 0.177 %																																			
Ca.....	18.07 ± 0.14 %	Na ₂ O..... 0.072 %	V 0.035 ± 0.002 %																																			
CaO.....	25.30 %	K..... 0.157 ± 0.007 %	Zn..... 0.021 ± 0.002 %																																			
Mg.....	0.70 ± 0.02 %	K ₂ O..... 0.189 %																																				
ECRM-F 685-1	Prereduced iron ore - powder Year of issue: 1983 Certified values <table> <tr><td>Fe_{total}</td><td>91.10 ± 0.27 %</td><td>Mg 0.239 ± 0.006 %</td><td>As..... (0.0017) %</td></tr> <tr><td>Fe(II)</td><td>(7.72) %</td><td>MgO..... 0.397 %</td><td>C 1.49 ± 0.02 %</td></tr> <tr><td>Femet.....</td><td>80.80 ± 0.44 %</td><td>Mn 0.042 ± 0.003 %</td><td>Co 0.013 ± 0.001 %</td></tr> <tr><td>Si.....</td><td>0.795 ± 0.052 %</td><td>P 0.017 ± 0.001 %</td><td>Cr (0.0035) %</td></tr> <tr><td>SiO₂.....</td><td>1.701 %</td><td>S..... 0.0031 ± 0.0005 %</td><td>Cu (0.0016) %</td></tr> <tr><td>Al.....</td><td>0.320 ± 0.014 %</td><td>Na..... 0.077 ± 0.003 %</td><td>Ni 0.018 ± 0.002 %</td></tr> <tr><td>Al₂O₃.....</td><td>0.604 %</td><td>Na₂O..... 0.104 %</td><td>Pb (0.0020) %</td></tr> <tr><td>Ca.....</td><td>0.140 ± 0.013 %</td><td>K..... 0.042 ± 0.004 %</td><td>Ti 0.220 ± 0.010 %</td></tr> <tr><td>CaO.....</td><td>0.195 %</td><td>K₂O..... 0.050 %</td><td>V 0.144 ± 0.009 %</td></tr> </table> (Values in parenthesis are indicative values)	Fe _{total}	91.10 ± 0.27 %	Mg 0.239 ± 0.006 %	As..... (0.0017) %	Fe(II)	(7.72) %	MgO..... 0.397 %	C 1.49 ± 0.02 %	Femet.....	80.80 ± 0.44 %	Mn 0.042 ± 0.003 %	Co 0.013 ± 0.001 %	Si.....	0.795 ± 0.052 %	P 0.017 ± 0.001 %	Cr (0.0035) %	SiO ₂	1.701 %	S..... 0.0031 ± 0.0005 %	Cu (0.0016) %	Al.....	0.320 ± 0.014 %	Na..... 0.077 ± 0.003 %	Ni 0.018 ± 0.002 %	Al ₂ O ₃	0.604 %	Na ₂ O..... 0.104 %	Pb (0.0020) %	Ca.....	0.140 ± 0.013 %	K..... 0.042 ± 0.004 %	Ti 0.220 ± 0.010 %	CaO.....	0.195 %	K ₂ O..... 0.050 %	V 0.144 ± 0.009 %	100 g
Fe _{total}	91.10 ± 0.27 %	Mg 0.239 ± 0.006 %	As..... (0.0017) %																																			
Fe(II)	(7.72) %	MgO..... 0.397 %	C 1.49 ± 0.02 %																																			
Femet.....	80.80 ± 0.44 %	Mn 0.042 ± 0.003 %	Co 0.013 ± 0.001 %																																			
Si.....	0.795 ± 0.052 %	P 0.017 ± 0.001 %	Cr (0.0035) %																																			
SiO ₂	1.701 %	S..... 0.0031 ± 0.0005 %	Cu (0.0016) %																																			
Al.....	0.320 ± 0.014 %	Na..... 0.077 ± 0.003 %	Ni 0.018 ± 0.002 %																																			
Al ₂ O ₃	0.604 %	Na ₂ O..... 0.104 %	Pb (0.0020) %																																			
Ca.....	0.140 ± 0.013 %	K..... 0.042 ± 0.004 %	Ti 0.220 ± 0.010 %																																			
CaO.....	0.195 %	K ₂ O..... 0.050 %	V 0.144 ± 0.009 %																																			
NIST-690	Iron ore, Canada - Constituents Certified values <table> <tr><td>Fe_{total}</td><td>66.85 %</td><td>CaO 0.20 %</td><td>Na₂O 0.003 %</td></tr> <tr><td>P</td><td>0.011 %</td><td>K₂O 0.0030 %</td><td>SiO₂..... 3.71 %</td></tr> <tr><td>S</td><td>0.003 %</td><td>MgO 0.18 %</td><td>TiO₂..... 0.022 %</td></tr> <tr><td>Al₂O₃</td><td>0.18 %</td><td>MnO..... 0.23 %</td><td></td></tr> </table>	Fe _{total}	66.85 %	CaO 0.20 %	Na ₂ O 0.003 %	P	0.011 %	K ₂ O 0.0030 %	SiO ₂ 3.71 %	S	0.003 %	MgO 0.18 %	TiO ₂ 0.022 %	Al ₂ O ₃	0.18 %	MnO..... 0.23 %		100 g																				
Fe _{total}	66.85 %	CaO 0.20 %	Na ₂ O 0.003 %																																			
P	0.011 %	K ₂ O 0.0030 %	SiO ₂ 3.71 %																																			
S	0.003 %	MgO 0.18 %	TiO ₂ 0.022 %																																			
Al ₂ O ₃	0.18 %	MnO..... 0.23 %																																				
NIST-691	Iron ore, reduced - Constituents Certified values <table> <tr><td>C</td><td>0.12 %</td><td>S 0.005 %</td><td>Na₂O 0.186 %</td></tr> <tr><td>Co</td><td>0.030 %</td><td>Al₂O₃..... 1.22 %</td><td>SiO₂..... 3.7 %</td></tr> <tr><td>Cu</td><td>0.032 %</td><td>CaO 0.63 %</td><td>TiO₂..... 0.27 %</td></tr> <tr><td>Fe_{total}</td><td>90.8 %</td><td>MgO 0.52 %</td><td></td></tr> <tr><td>P</td><td>0.006 %</td><td>MnO..... 0.043 %</td><td></td></tr> </table>	C	0.12 %	S 0.005 %	Na ₂ O 0.186 %	Co	0.030 %	Al ₂ O ₃ 1.22 %	SiO ₂ 3.7 %	Cu	0.032 %	CaO 0.63 %	TiO ₂ 0.27 %	Fe _{total}	90.8 %	MgO 0.52 %		P	0.006 %	MnO..... 0.043 %		100 g																
C	0.12 %	S 0.005 %	Na ₂ O 0.186 %																																			
Co	0.030 %	Al ₂ O ₃ 1.22 %	SiO ₂ 3.7 %																																			
Cu	0.032 %	CaO 0.63 %	TiO ₂ 0.27 %																																			
Fe _{total}	90.8 %	MgO 0.52 %																																				
P	0.006 %	MnO..... 0.043 %																																				
NIST-692	Iron ore, Labrador - Constituents Certified values <table> <tr><td>Fe_{total}</td><td>59.58 %</td><td>CaO 0.023 %</td><td>Na₂O 0.008 %</td></tr> <tr><td>P</td><td>0.039 %</td><td>K₂O 0.039 %</td><td>SiO₂..... 3.7 %</td></tr> <tr><td>S</td><td>0.005 %</td><td>MgO 0.035 %</td><td>TiO₂..... 0.045 %</td></tr> <tr><td>Al₂O₃</td><td>1.41 %</td><td>MnO..... 0.46 %</td><td></td></tr> </table>	Fe _{total}	59.58 %	CaO 0.023 %	Na ₂ O 0.008 %	P	0.039 %	K ₂ O 0.039 %	SiO ₂ 3.7 %	S	0.005 %	MgO 0.035 %	TiO ₂ 0.045 %	Al ₂ O ₃	1.41 %	MnO..... 0.46 %		100 g																				
Fe _{total}	59.58 %	CaO 0.023 %	Na ₂ O 0.008 %																																			
P	0.039 %	K ₂ O 0.039 %	SiO ₂ 3.7 %																																			
S	0.005 %	MgO 0.035 %	TiO ₂ 0.045 %																																			
Al ₂ O ₃	1.41 %	MnO..... 0.46 %																																				
NIST-693	Iron ore, Nimba - Constituents Certified values <table> <tr><td>Fe_{total}</td><td>65.11 %</td><td>CaO 0.016 %</td><td>Na₂O 0.0028 %</td></tr> <tr><td>P</td><td>0.056 %</td><td>K₂O 0.0028 %</td><td>SiO₂..... 3.87 %</td></tr> <tr><td>S</td><td>0.005 %</td><td>MgO 0.013 %</td><td>TiO₂..... 0.035 %</td></tr> <tr><td>Al₂O₃</td><td>1.02 %</td><td>MnO..... 0.091 %</td><td></td></tr> </table>	Fe _{total}	65.11 %	CaO 0.016 %	Na ₂ O 0.0028 %	P	0.056 %	K ₂ O 0.0028 %	SiO ₂ 3.87 %	S	0.005 %	MgO 0.013 %	TiO ₂ 0.035 %	Al ₂ O ₃	1.02 %	MnO..... 0.091 %		100 g																				
Fe _{total}	65.11 %	CaO 0.016 %	Na ₂ O 0.0028 %																																			
P	0.056 %	K ₂ O 0.0028 %	SiO ₂ 3.87 %																																			
S	0.005 %	MgO 0.013 %	TiO ₂ 0.035 %																																			
Al ₂ O ₃	1.02 %	MnO..... 0.091 %																																				

Ceramic materials and minerals

ECRM-D 777-1	Silica brick - powder Year of issue: 1983 Certified values <table> <tr><td>Si.....</td><td>44.44 ± 0.15 %</td><td>MgO..... 0.071 ± 0.012 %</td><td>K 0.13 ± 0.02 %</td></tr> <tr><td>SiO₂.....</td><td>95.06 ± 0.32 %</td><td>Al 0.42 ± 0.02 %</td><td>K₂O 0.15 ± 0.02 %</td></tr> <tr><td>Ca.....</td><td>2.02 ± 0.08 %</td><td>Al₂O₃..... 0.80 ± 0.04 %</td><td>Na 0.02 %</td></tr> <tr><td>CaO.....</td><td>2.83 ± 0.10 %</td><td>Fe 0.23 ± 0.03 %</td><td>Ti 0.27 ± 0.02 %</td></tr> <tr><td>Mg.....</td><td>0.043 ± 0.007 %</td><td>Fe₂O₃..... 0.33 ± 0.04 %</td><td></td></tr> </table>	Si.....	44.44 ± 0.15 %	MgO..... 0.071 ± 0.012 %	K 0.13 ± 0.02 %	SiO ₂	95.06 ± 0.32 %	Al 0.42 ± 0.02 %	K ₂ O 0.15 ± 0.02 %	Ca.....	2.02 ± 0.08 %	Al ₂ O ₃ 0.80 ± 0.04 %	Na 0.02 %	CaO.....	2.83 ± 0.10 %	Fe 0.23 ± 0.03 %	Ti 0.27 ± 0.02 %	Mg.....	0.043 ± 0.007 %	Fe ₂ O ₃ 0.33 ± 0.04 %		100 g
Si.....	44.44 ± 0.15 %	MgO..... 0.071 ± 0.012 %	K 0.13 ± 0.02 %																			
SiO ₂	95.06 ± 0.32 %	Al 0.42 ± 0.02 %	K ₂ O 0.15 ± 0.02 %																			
Ca.....	2.02 ± 0.08 %	Al ₂ O ₃ 0.80 ± 0.04 %	Na 0.02 %																			
CaO.....	2.83 ± 0.10 %	Fe 0.23 ± 0.03 %	Ti 0.27 ± 0.02 %																			
Mg.....	0.043 ± 0.007 %	Fe ₂ O ₃ 0.33 ± 0.04 %																				
ECRM-D 779-1	Magnesite - powder Year of issue: 1991 Certified values <table> <tr><td>Si.....</td><td>0.182 ± 0.015 %</td><td>B 0.0116 ± 0.0012 %</td><td>Mn 0.503 ± 0.017 %</td></tr> <tr><td>Ca.....</td><td>1.691 ± 0.023 %</td><td>Cr (0.0030) %</td><td>Na (0.0058) %</td></tr> <tr><td>Mg.....</td><td>(54.57) %</td><td>Fe 3.73 ± 0.06 %</td><td>P 0.0267 ± 0.0026 %</td></tr> <tr><td>Al.....</td><td>0.105 ± 0.007 %</td><td>K (0.0020) %</td><td>Ti 0.0081 ± 0.0012 %</td></tr> </table> (Values in parenthesis are indicative values)	Si.....	0.182 ± 0.015 %	B 0.0116 ± 0.0012 %	Mn 0.503 ± 0.017 %	Ca.....	1.691 ± 0.023 %	Cr (0.0030) %	Na (0.0058) %	Mg.....	(54.57) %	Fe 3.73 ± 0.06 %	P 0.0267 ± 0.0026 %	Al.....	0.105 ± 0.007 %	K (0.0020) %	Ti 0.0081 ± 0.0012 %	100 g				
Si.....	0.182 ± 0.015 %	B 0.0116 ± 0.0012 %	Mn 0.503 ± 0.017 %																			
Ca.....	1.691 ± 0.023 %	Cr (0.0030) %	Na (0.0058) %																			
Mg.....	(54.57) %	Fe 3.73 ± 0.06 %	P 0.0267 ± 0.0026 %																			
Al.....	0.105 ± 0.007 %	K (0.0020) %	Ti 0.0081 ± 0.0012 %																			

Iron and steel products

Code	Product	Unit
ECRM-B 752-1	Limestone - powder Year of issue: 1984 Certified values SiO ₂0.70 % Fe ₂ O ₃0.045 % S.....0.007 % CaO.....55.4 % K ₂ O.....0.02 % SrO.....0.019 % MgO.....0.15 % MnO.....0.010 % TiO ₂0.009 % Al ₂ O ₃0.12 % Na ₂ O.....(0.02) % L.O.I*.....43.4 % BaO.....0.006 % P ₂ O ₅(0.0055) %	100 g
	* Loss On Ignition	
ECRM-B 776-1	Firebrick - powder Year of issue: 1983 Certified values Si.....29.34 ± 0.18 % BaO0.122 ± 0.011 % Na0.362 ± 0.015 % SiO ₂62.76 ± 0.39 % Cr.....0.015 ± 0.002 % Na ₂ O0.488 ± 0.020 % Ca0.221 ± 0.013 % Cr ₂ O ₃0.022 ± 0.003 % P0.027 ± 0.003 % CaO.....0.310 ± 0.018 % Fe0.999 ± 0.017 % P ₂ O ₅0.062 ± 0.007 % Mg0.287 ± 0.010 % Fe ₂ O ₃1.43 ± 0.02 % Ti0.969 ± 0.015 % MgO.....0.476 ± 0.016 % K2.42 ± 0.05 % TiO ₂1.62 ± 0.03 % Al15.50 ± 0.12 % K ₂ O2.92 ± 0.07 % Zr(0.030) % Al ₂ O ₃29.28 ± 0.22 % Li0.009 ± 0.001 % L.O.I.*(0.3) % Ba.....0.109 ± 0.010 % Li ₂ O0.019 ± 0.002 %	100 g
	(Values in parenthesis are indicative values)	
	* Loss On Ignition	
ECRM-B 781-1	Silicon carbide refractory with 35.6% silicon and 48.3% carbon - powder Year of issue: 1993 Certified values Si.....35.56 ± 0.28 % Cfrei(37.22) % Na(0.0308) % Ca(0.0433) % Fe(0.806) % Ni(0.0210) % Mg(0.0421) % K(0.3765) % P(0.0117) % Al4.39 ± 0.04 % Mn(0.0274) % Ti(0.0320) % B(0.0149) % Mo(0.0264) % V(0.0216) % C48.25 ± 0.11 % N(0.0282) %	100 g
	(Values in parenthesis are indicative values)	
ECRM-B 782-1	Dolomite - powder Year of issue: 1996 Certified values Si.....0.124 ± 0.003 %* B(0.0012) % P0.0056 ± 0.0003 %* SiO ₂0.266 ± 0.007 %* Ba(0.0007) % P ₂ O ₅0.0128 ± 0.0007 %* Ca21.68 ± 0.05 %* Fe0.314 ± 0.005 %* Pb0.0027 ± 0.0002 %* CaO.....30.34 ± 0.07 %* K0.0216 ± 0.0013 %* S(0.0159) % Mg12.84 ± 0.09 %* K ₂ O0.0260 ± 0.0016 %* Ti0.0025 ± 0.0002 %* MgO.....21.29 ± 0.15 %* Mn0.063 ± 0.002 %* TiO ₂0.0042 ± 0.0004 %* Al0.055 ± 0.003 %* MnO0.081 ± 0.003 %* Zn0.0066 ± 0.0003 %* Al ₂ O ₃0.104 ± 0.006 %* Ni(0.0004) % L.O.I.*47.25 ± 0.12 %*	100 g
	(Values in parenthesis are indicative values)	
	* 95%-confidence interval	
	* Loss On Ignition	
ECRM-B 783-1	Tungsten carbide - powder Year of issue: 2005 Certified values C6.188 ± 0.011+ % Fe0.0022 ± 0.0002+ % C _{free}(0.042) % O(0.016) %	100 g
	(Values in parenthesis are indicative values)	
ECRM-F 701-1	Calcite - powder Year of issue: 1970 Certified values Si0.93 ± 0.03 % Fe0.73 ± 0.05 % P0.022 ± 0.004 % Ca37.66 ± 0.22 % K(0.090) % Pb(< 0.005) % Mg0.36 ± 0.04 % Mn0.022 ± 0.004 % S0.040 ± 0.010 % Al0.29 ± 0.04 % Mo(< 0.003) % Ti0.018 ± 0.003 % C(11.5) % Na(0.020) % Zn(< 0.005) % Cr(< 0.002) % Ni(0.003) % L.O.I.*42.4 ± 0.22 %	100 g
	(Values in parenthesis are indicative values)	
	* Loss On Ignition	
ECRM-F 702-1	Dolomite - powder Year of issue: 1971 Certified values Si1.04 ± 0.03 % Al0.21 ± 0.04 % P0.024 ± 0.003 % Ca21.48 ± 0.27 % Fe0.440 ± 0.025 % S0.027 ± 0.010 % Mg12.37 ± 0.33 % Mn0.098 ± 0.006 % Ti0.013 ± 0.005 %	100 g

Iron and steel products

Code	Product	Unit
ECRM-F 778-1	High carbon magnesia - powder Year of issue: 1986 Certified values Si 0.489 ± 0.019 % B ₂ O ₃ 0.0039 % Ni (0.007) % SiO ₂ 1.046 % C 14.00 ± 0.22 % P (0.004) % Ca 0.883 ± 0.024 % Cr 0.102 ± 0.004 % Pb (< 0.001) % CaO 1.236 % Cr ₂ O ₃ 0.149 % TiO ₂ (0.008) % Mg 48.87 ± 0.27 % Fe 0.67 ± 0.05 % Zn (< 0.002) % MgO 81.03 % K (0.020) % Zr (< 0.002) % Al 0.297 ± 0.008 % Mn 0.011 ± 0.001 % L.O.I.* 15.38 ± 0.14 % Al ₂ O ₃ 0.561 % MnO 0.014 % B 0.0012 ± 0.0001 % Na (0.023) % (Values in parenthesis are indicative values) * Loss On Ignition	100 g
ECRM-F 780-1	Silicon carbide refractory with 63.5% silicon and 26.4% carbon - powder Year of issue: 1992 Certified values Si 63.5 ± 0.5 % Cr (0.010) % Ni (0.015) % Ca 0.84 ± 0.03 % Fe 1.30 ± 0.06 % O (5.24) % Mg 0.051 ± 0.011 % K (0.011) % P (< 0.010) % Al 1.86 ± 0.06 % Mn 0.029 ± 0.004 % Ti (0.050) % C 26.38 ± 0.23 % N 0.32 ± 0.05 % V (0.025) % Cfrei (0.563) % Na (0.050) % Zr (0.01) % (Values in parenthesis are indicative values)	100 g
BCS CRM309	Sillimanite Certified values SiO ₂ 34.1 % CaO 0.22 % BaO (0.006) % Al ₂ O ₃ 61.1 % MgO 0.17 % Li ₂ O (0.01) % TiO ₂ 1.92 % Na ₂ O 0.34 % SrO (0.003) % Fe ₂ O ₃ 1.51 % K ₂ O 0.46 % MnO (0.03) % L.O.I. (0.1) % (Values in parenthesis are indicative values)	100 g
BCS CRM313/1	High purity silica Certified values SiO ₂ 99.78 % CaO 0.006 % L.O.I. (0.1) % Al ₂ O ₃ 0.036 % MgO 0.0013 % Cr ₂ O ₃ (<0.0002) % TiO ₂ 0.017 % Na ₂ O 0.003 % Li ₂ O (0.0005) % Fe ₂ O ₃ 0.012 % K ₂ O 0.005 % MnO 0.00013 % ZrO ₂ (0.002) % (Values in parenthesis are indicative values)	100 g
BCS CRM319/1	Magnesite Certified values SiO ₂ 1.093 % CaO 3.00 % BaO (0.0038) % Al ₂ O ₃ 0.109 % MgO 95.38 % Cr ₂ O ₃ 0.0035 % TiO ₂ 0.007 % P ₂ O ₅ (0.007) % SrO (0.006) % Fe ₂ O ₃ 0.291 % ZrO ₂ (0.0008) % Y ₂ O ₃ (0.0014) % MnO 0.108 % B ₂ O ₃ (0.002) % Ni (0.0075) % (Values in parenthesis are indicative values)	100 g
BCS CRM348	Ball clay Certified values SiO ₂ 51.1 % MgO 0.3 % S (0.1) % Al ₂ O ₃ 31.6 % Na ₂ O 0.34 % L.O.I. 11.8 % TiO ₂ 1.08 % K ₂ O 2.23 % BaO (0.04) % Fe ₂ O ₃ 1.04 % P ₂ O ₅ 0.071 % Cr ₂ O ₃ 0.016 % CaO 0.17 % ZrO ₂ (0.03) % C (1.64) % (Values in parenthesis are indicative values)	100 g
BCS CRM358	Zirconia Certified values SiO ₂ 0.2 % MgO 3.42 % BaO 0.1 % Al ₂ O ₃ 0.08 % Na ₂ O (<0.01) % HfO ₂ 1.63 % TiO ₂ 0.2 % K ₂ O (<0.01) % SrO 0.07 % Fe ₂ O ₃ 0.064 % ZrO ₂ 92.7 % ThO ₂ (0.0007) % CaO 1.5 % L.O.I. 0.08 % U ₃ O ₈ (0.08) % (Values in parenthesis are indicative values)	100 g

Iron and steel products

Code	Product	Unit																																										
BCS CRM362	Mine tailings sample Certified values	100 g																																										
	<table> <tbody> <tr><td>SiO₂.....</td><td>9.03 %</td><td>Na₂O.....</td><td>0.084 %</td><td>BaO.....</td><td>(2.02) %</td></tr> <tr><td>Al₂O₃.....</td><td>0.667 %</td><td>K₂O.....</td><td>0.14 %</td><td>Cr₂O₃.....</td><td>(0.003) %</td></tr> <tr><td>TiO₂.....</td><td>0.047 %</td><td>PbO.....</td><td>2.63 %</td><td>SrO.....</td><td>0.034 %</td></tr> <tr><td>Fe₂O₃.....</td><td>0.483 %</td><td>ZnO.....</td><td>2.59 %</td><td>C.....</td><td>(9.9) %</td></tr> <tr><td>Mn₃O₄.....</td><td>0.829 %</td><td>P₂O₅.....</td><td>(0.014) %</td><td>Cd.....</td><td>0.02 %</td></tr> <tr><td>CaO.....</td><td>44.21 %</td><td>S.....</td><td>1.48 %</td><td>Ni.....</td><td>(0.001) %</td></tr> <tr><td>MgO.....</td><td>0.068 %</td><td>L.O.I.....</td><td>32.81 %</td><td></td><td></td></tr> </tbody> </table> <p>(Values in parenthesis are indicative values)</p>	SiO ₂	9.03 %	Na ₂ O.....	0.084 %	BaO.....	(2.02) %	Al ₂ O ₃	0.667 %	K ₂ O.....	0.14 %	Cr ₂ O ₃	(0.003) %	TiO ₂	0.047 %	PbO.....	2.63 %	SrO.....	0.034 %	Fe ₂ O ₃	0.483 %	ZnO.....	2.59 %	C.....	(9.9) %	Mn ₃ O ₄	0.829 %	P ₂ O ₅	(0.014) %	Cd.....	0.02 %	CaO.....	44.21 %	S.....	1.48 %	Ni.....	(0.001) %	MgO.....	0.068 %	L.O.I.....	32.81 %			
SiO ₂	9.03 %	Na ₂ O.....	0.084 %	BaO.....	(2.02) %																																							
Al ₂ O ₃	0.667 %	K ₂ O.....	0.14 %	Cr ₂ O ₃	(0.003) %																																							
TiO ₂	0.047 %	PbO.....	2.63 %	SrO.....	0.034 %																																							
Fe ₂ O ₃	0.483 %	ZnO.....	2.59 %	C.....	(9.9) %																																							
Mn ₃ O ₄	0.829 %	P ₂ O ₅	(0.014) %	Cd.....	0.02 %																																							
CaO.....	44.21 %	S.....	1.48 %	Ni.....	(0.001) %																																							
MgO.....	0.068 %	L.O.I.....	32.81 %																																									
BCS CRM369	Magnesite-Chrome Certified values	100 g																																										
	<table> <tbody> <tr><td>SiO₂.....</td><td>2.59 %</td><td>CaO.....</td><td>1.17 %</td><td>Cr₂O₃.....</td><td>17.2 %</td></tr> <tr><td>Al₂O₃.....</td><td>14.7 %</td><td>MgO.....</td><td>53.5 %</td><td>Li₂O.....</td><td>0.03 %</td></tr> <tr><td>TiO₂.....</td><td>0.14 %</td><td>Na₂O.....</td><td>0.05 %</td><td>SrO.....</td><td><0.01 %</td></tr> <tr><td>Fe₂O₃.....</td><td>10.3 %</td><td>K₂O.....</td><td>0.03 %</td><td>Ni.....</td><td>(0.15) %</td></tr> <tr><td>MnO.....</td><td>0.11 %</td><td>BaO.....</td><td>(<0.01) %</td><td></td><td></td></tr> </tbody> </table> <p>(Values in parenthesis are indicative values)</p>	SiO ₂	2.59 %	CaO.....	1.17 %	Cr ₂ O ₃	17.2 %	Al ₂ O ₃	14.7 %	MgO.....	53.5 %	Li ₂ O.....	0.03 %	TiO ₂	0.14 %	Na ₂ O.....	0.05 %	SrO.....	<0.01 %	Fe ₂ O ₃	10.3 %	K ₂ O.....	0.03 %	Ni.....	(0.15) %	MnO.....	0.11 %	BaO.....	(<0.01) %															
SiO ₂	2.59 %	CaO.....	1.17 %	Cr ₂ O ₃	17.2 %																																							
Al ₂ O ₃	14.7 %	MgO.....	53.5 %	Li ₂ O.....	0.03 %																																							
TiO ₂	0.14 %	Na ₂ O.....	0.05 %	SrO.....	<0.01 %																																							
Fe ₂ O ₃	10.3 %	K ₂ O.....	0.03 %	Ni.....	(0.15) %																																							
MnO.....	0.11 %	BaO.....	(<0.01) %																																									
BCS CRM370	Magnesite-Chrome Certified values	100 g																																										
	<table> <tbody> <tr><td>SiO₂.....</td><td>3.01 %</td><td>CaO.....</td><td>1.54 %</td><td>Cr₂O₃.....</td><td>13.4 %</td></tr> <tr><td>Al₂O₃.....</td><td>12.3 %</td><td>MgO.....</td><td>61.8 %</td><td>Li₂O.....</td><td>0.03 %</td></tr> <tr><td>TiO₂.....</td><td>0.13 %</td><td>Na₂O.....</td><td>0.06 %</td><td>SrO.....</td><td><0.01 %</td></tr> <tr><td>Fe₂O₃.....</td><td>7.23 %</td><td>K₂O.....</td><td>0.03 %</td><td>Ni.....</td><td>0.08 %</td></tr> <tr><td>MnO.....</td><td>0.11 %</td><td>BaO.....</td><td><0.01 %</td><td></td><td></td></tr> </tbody> </table>	SiO ₂	3.01 %	CaO.....	1.54 %	Cr ₂ O ₃	13.4 %	Al ₂ O ₃	12.3 %	MgO.....	61.8 %	Li ₂ O.....	0.03 %	TiO ₂	0.13 %	Na ₂ O.....	0.06 %	SrO.....	<0.01 %	Fe ₂ O ₃	7.23 %	K ₂ O.....	0.03 %	Ni.....	0.08 %	MnO.....	0.11 %	BaO.....	<0.01 %															
SiO ₂	3.01 %	CaO.....	1.54 %	Cr ₂ O ₃	13.4 %																																							
Al ₂ O ₃	12.3 %	MgO.....	61.8 %	Li ₂ O.....	0.03 %																																							
TiO ₂	0.13 %	Na ₂ O.....	0.06 %	SrO.....	<0.01 %																																							
Fe ₂ O ₃	7.23 %	K ₂ O.....	0.03 %	Ni.....	0.08 %																																							
MnO.....	0.11 %	BaO.....	<0.01 %																																									
BCS CRM375/1	Soda feldspar Certified values	100 g																																										
	<table> <tbody> <tr><td>SiO₂.....</td><td>69.26 %</td><td>K₂O.....</td><td>1.47 %</td><td>Cr₂O₃.....</td><td>(0.0018) %</td></tr> <tr><td>Al₂O₃.....</td><td>17.89 %</td><td>PbO.....</td><td>(0.0004) %</td><td>HfO₂.....</td><td>(0.0004) %</td></tr> <tr><td>TiO₂.....</td><td>0.313 %</td><td>ZnO.....</td><td>(0.0005) %</td><td>SrO.....</td><td>(0.012) %</td></tr> <tr><td>Fe₂O₃.....</td><td>0.291 %</td><td>P₂O₅.....</td><td>0.226 %</td><td>ThO₂.....</td><td>(0.0011) %</td></tr> <tr><td>CaO.....</td><td>0.78 %</td><td>ZrO₂.....</td><td>(0.0107) %</td><td>U₃O₈.....</td><td>(0.0002) %</td></tr> <tr><td>MgO.....</td><td>0.18 %</td><td>L.O.I.....</td><td>0.72 %</td><td>Y₂O₃.....</td><td>(0.0023) %</td></tr> <tr><td>Na₂O.....</td><td>8.89 %</td><td>BaO.....</td><td>(0.0106) %</td><td></td><td></td></tr> </tbody> </table> <p>(Values in parenthesis are indicative values)</p>	SiO ₂	69.26 %	K ₂ O.....	1.47 %	Cr ₂ O ₃	(0.0018) %	Al ₂ O ₃	17.89 %	PbO.....	(0.0004) %	HfO ₂	(0.0004) %	TiO ₂	0.313 %	ZnO.....	(0.0005) %	SrO.....	(0.012) %	Fe ₂ O ₃	0.291 %	P ₂ O ₅	0.226 %	ThO ₂	(0.0011) %	CaO.....	0.78 %	ZrO ₂	(0.0107) %	U ₃ O ₈	(0.0002) %	MgO.....	0.18 %	L.O.I.....	0.72 %	Y ₂ O ₃	(0.0023) %	Na ₂ O.....	8.89 %	BaO.....	(0.0106) %			
SiO ₂	69.26 %	K ₂ O.....	1.47 %	Cr ₂ O ₃	(0.0018) %																																							
Al ₂ O ₃	17.89 %	PbO.....	(0.0004) %	HfO ₂	(0.0004) %																																							
TiO ₂	0.313 %	ZnO.....	(0.0005) %	SrO.....	(0.012) %																																							
Fe ₂ O ₃	0.291 %	P ₂ O ₅	0.226 %	ThO ₂	(0.0011) %																																							
CaO.....	0.78 %	ZrO ₂	(0.0107) %	U ₃ O ₈	(0.0002) %																																							
MgO.....	0.18 %	L.O.I.....	0.72 %	Y ₂ O ₃	(0.0023) %																																							
Na ₂ O.....	8.89 %	BaO.....	(0.0106) %																																									
BCS CRM376/1	Potash feldspar Certified values	100 g																																										
	<table> <tbody> <tr><td>SiO₂.....</td><td>65.77 %</td><td>CaO.....</td><td>0.421 %</td><td>P₂O₅.....</td><td>(0.02) %</td></tr> <tr><td>Al₂O₃.....</td><td>18.63 %</td><td>MgO.....</td><td>(0.03) %</td><td>ZrO₂.....</td><td>(<0.01) %</td></tr> <tr><td>TiO₂.....</td><td>(<0.01) %</td><td>Na₂O.....</td><td>3.00 %</td><td>L.O.I.....</td><td>0.203 %</td></tr> <tr><td>Fe₂O₃.....</td><td>0.085 %</td><td>K₂O.....</td><td>11.59 %</td><td>BaO.....</td><td>0.021 %</td></tr> <tr><td>Mn₃O₄.....</td><td>(0.004) %</td><td>PbO.....</td><td>0.009 %</td><td>Cr₂O₃.....</td><td>0.001 %</td></tr> </tbody> </table> <p>(Values in parenthesis are indicative values)</p>	SiO ₂	65.77 %	CaO.....	0.421 %	P ₂ O ₅	(0.02) %	Al ₂ O ₃	18.63 %	MgO.....	(0.03) %	ZrO ₂	(<0.01) %	TiO ₂	(<0.01) %	Na ₂ O.....	3.00 %	L.O.I.....	0.203 %	Fe ₂ O ₃	0.085 %	K ₂ O.....	11.59 %	BaO.....	0.021 %	Mn ₃ O ₄	(0.004) %	PbO.....	0.009 %	Cr ₂ O ₃	0.001 %													
SiO ₂	65.77 %	CaO.....	0.421 %	P ₂ O ₅	(0.02) %																																							
Al ₂ O ₃	18.63 %	MgO.....	(0.03) %	ZrO ₂	(<0.01) %																																							
TiO ₂	(<0.01) %	Na ₂ O.....	3.00 %	L.O.I.....	0.203 %																																							
Fe ₂ O ₃	0.085 %	K ₂ O.....	11.59 %	BaO.....	0.021 %																																							
Mn ₃ O ₄	(0.004) %	PbO.....	0.009 %	Cr ₂ O ₃	0.001 %																																							
BCS CRM388	Zircon Certified values	100 g																																										
	<table> <tbody> <tr><td>SiO₂.....</td><td>32.7 %</td><td>MgO.....</td><td>(<0.05) %</td><td>L.O.I.....</td><td>(0.2) %</td></tr> <tr><td>Al₂O₃.....</td><td>0.291 %</td><td>Na₂O.....</td><td>(<0.02) %</td><td>HfO₂.....</td><td>1.3 %</td></tr> <tr><td>TiO₂.....</td><td>0.232 %</td><td>K₂O.....</td><td>(<0.03) %</td><td>ThO₂.....</td><td>0.018 %</td></tr> <tr><td>Fe₂O₃.....</td><td>0.049 %</td><td>P₂O₅.....</td><td>0.12 %</td><td>U₃O₈.....</td><td>0.034 %</td></tr> <tr><td>CaO.....</td><td>(0.04) %</td><td>ZrO₂.....</td><td>64.9 %</td><td>Y₂O₃.....</td><td>0.136 %</td></tr> </tbody> </table> <p>(Values in parenthesis are indicative values)</p>	SiO ₂	32.7 %	MgO.....	(<0.05) %	L.O.I.....	(0.2) %	Al ₂ O ₃	0.291 %	Na ₂ O.....	(<0.02) %	HfO ₂	1.3 %	TiO ₂	0.232 %	K ₂ O.....	(<0.03) %	ThO ₂	0.018 %	Fe ₂ O ₃	0.049 %	P ₂ O ₅	0.12 %	U ₃ O ₈	0.034 %	CaO.....	(0.04) %	ZrO ₂	64.9 %	Y ₂ O ₃	0.136 %													
SiO ₂	32.7 %	MgO.....	(<0.05) %	L.O.I.....	(0.2) %																																							
Al ₂ O ₃	0.291 %	Na ₂ O.....	(<0.02) %	HfO ₂	1.3 %																																							
TiO ₂	0.232 %	K ₂ O.....	(<0.03) %	ThO ₂	0.018 %																																							
Fe ₂ O ₃	0.049 %	P ₂ O ₅	0.12 %	U ₃ O ₈	0.034 %																																							
CaO.....	(0.04) %	ZrO ₂	64.9 %	Y ₂ O ₃	0.136 %																																							
BCS CRM389/1	High purity magnesia Certified values	100 g																																										
	<table> <tbody> <tr><td>SiO₂.....</td><td>0.274 %</td><td>MgO.....</td><td>97.89 %</td><td>Cr₂O₃.....</td><td>0.004 %</td></tr> <tr><td>Al₂O₃.....</td><td>0.104 %</td><td>ZnO.....</td><td>(0.0029) %</td><td>SrO.....</td><td>0.0007 %</td></tr> <tr><td>TiO₂.....</td><td>0.0051 %</td><td>P₂O₅.....</td><td>0.0295 %</td><td>Y₂O₃.....</td><td>0.0029 %</td></tr> <tr><td>Fe₂O₃.....</td><td>0.607 %</td><td>ZrO₂.....</td><td>(0.0008) %</td><td>Ni.....</td><td>0.0012 %</td></tr> <tr><td>MnO.....</td><td>0.1 %</td><td>B₂O₃.....</td><td>0.015 %</td><td></td><td></td></tr> <tr><td>CaO.....</td><td>0.88 %</td><td>BaO.....</td><td>0.0015 %</td><td></td><td></td></tr> </tbody> </table> <p>(Values in parenthesis are indicative values)</p>	SiO ₂	0.274 %	MgO.....	97.89 %	Cr ₂ O ₃	0.004 %	Al ₂ O ₃	0.104 %	ZnO.....	(0.0029) %	SrO.....	0.0007 %	TiO ₂	0.0051 %	P ₂ O ₅	0.0295 %	Y ₂ O ₃	0.0029 %	Fe ₂ O ₃	0.607 %	ZrO ₂	(0.0008) %	Ni.....	0.0012 %	MnO.....	0.1 %	B ₂ O ₃	0.015 %			CaO.....	0.88 %	BaO.....	0.0015 %									
SiO ₂	0.274 %	MgO.....	97.89 %	Cr ₂ O ₃	0.004 %																																							
Al ₂ O ₃	0.104 %	ZnO.....	(0.0029) %	SrO.....	0.0007 %																																							
TiO ₂	0.0051 %	P ₂ O ₅	0.0295 %	Y ₂ O ₃	0.0029 %																																							
Fe ₂ O ₃	0.607 %	ZrO ₂	(0.0008) %	Ni.....	0.0012 %																																							
MnO.....	0.1 %	B ₂ O ₃	0.015 %																																									
CaO.....	0.88 %	BaO.....	0.0015 %																																									
BCS CRM394	Calcined bauxite Certified values	100 g																																										
	<table> <tbody> <tr><td>SiO₂.....</td><td>4.98 %</td><td>MgO.....</td><td>0.12 %</td><td>L.O.I.....</td><td>(0.4) %</td></tr> <tr><td>Al₂O₃.....</td><td>88.8 %</td><td>Na₂O.....</td><td>0.02 %</td><td>Cr₂O₃.....</td><td>(0.08) %</td></tr> <tr><td>TiO₂.....</td><td>3.11 %</td><td>K₂O.....</td><td>0.02 %</td><td>Li₂O.....</td><td><0.01 %</td></tr> <tr><td>Fe₂O₃.....</td><td>1.9 %</td><td>P₂O₅.....</td><td>0.22 %</td><td></td><td></td></tr> <tr><td>CaO.....</td><td>0.08 %</td><td>ZrO₂.....</td><td>(0.159) %</td><td></td><td></td></tr> </tbody> </table> <p>(Values in parenthesis are indicative values)</p>	SiO ₂	4.98 %	MgO.....	0.12 %	L.O.I.....	(0.4) %	Al ₂ O ₃	88.8 %	Na ₂ O.....	0.02 %	Cr ₂ O ₃	(0.08) %	TiO ₂	3.11 %	K ₂ O.....	0.02 %	Li ₂ O.....	<0.01 %	Fe ₂ O ₃	1.9 %	P ₂ O ₅	0.22 %			CaO.....	0.08 %	ZrO ₂	(0.159) %															
SiO ₂	4.98 %	MgO.....	0.12 %	L.O.I.....	(0.4) %																																							
Al ₂ O ₃	88.8 %	Na ₂ O.....	0.02 %	Cr ₂ O ₃	(0.08) %																																							
TiO ₂	3.11 %	K ₂ O.....	0.02 %	Li ₂ O.....	<0.01 %																																							
Fe ₂ O ₃	1.9 %	P ₂ O ₅	0.22 %																																									
CaO.....	0.08 %	ZrO ₂	(0.159) %																																									

Code	Product	Unit
BCS CRM396	Low silica magnesite chrome Certified values	100 g
	SiO ₂1.37 % Al ₂ O ₃5.73 % TiO ₂0.26 % Fe ₂ O ₃10.9 % MnO.....0.17 % (Values in parenthesis are indicative values)	CaO.....1.12 % MgO.....64.6 % Na ₂ O.....(0.06) % K ₂ O.....(0.03) % L.O.I.....(0.04) % B ₂ O ₃0.09 % Cr ₂ O ₃15.6 % Li ₂ O.....(0.05) %
BCS CRM512	Dolomite Certified values	100 g
	SiO ₂0.379 % Al ₂ O ₃0.055 % TiO ₂0.002 % Fe ₂ O ₃0.03 % MnO.....0.0036 % CaO.....30.61 % MgO.....21.59 % Na ₂ O.....(0.1) % (Values in parenthesis are indicative values)	K ₂ O.....(<0.02) % PbO.....(<0.001) % ZnO.....(<0.01) % P ₂ O ₅(<0.02) % S.....(<0.05) % L.O.I.....46.8 % BaO.....(<0.02) % Cr ₂ O ₃(<0.001) % SrO.....0.024 % As.....(<0.003) % C.....(12.4) % Cd.....(<0.0003) % F.....(0.01) % Ni.....(<0.001) %
BCS CRM513	Limestone Certified values	100 g
	SiO ₂0.228 % Al ₂ O ₃0.108 % TiO ₂(0.004) % Fe ₂ O ₃0.0275 % MnO.....0.0095 % CaO.....55.59 % MgO.....0.182 % Na ₂ O.....(<0.3) % (Values in parenthesis are indicative values)	K ₂ O.....0.015 % PbO.....0.0009Pb % ZnO.....0.0014Zn % P ₂ O ₅(0.005) % S.....0.0097 % L.O.I.....43.61 % BaO.....0.01 % Cr ₂ O ₃0.0012 % SrO.....0.0176 % As.....<0.001 % C.....11.9 % Cd.....<0.001 % F.....0.002 % Ni.....<0.001 %

Slags

ECRM-D 826-1	Phosphate slag - powder Year of issue: 1976 Certified values	100 g
	SiO ₂8.96 ± 0.15 % Al.....0.696 ± 0.008 % CaO.....46.48 ± 0.54 % MgO.....(2.46) % P ₂ O ₅14.65 ± 0.15 % P ₂ O ₅ citric acid sol....10.73 ± 0.14 % B.....(0.0029) % (Values in parenthesis are indicative values)	Cr.....0.182 ± 0.005 % Cu.....(0.0019) % F.....(0.3667) % Fe _{total}(20.73) % K.....0.0278 ± 0.0017 % Mn _{total}(3.46) % Mo.....(0.0011) % Na.....0.375 ± 0.009 % Ni.....(0.0017) % Pb.....(0.0049) % V.....0.503 ± 0.008 % V ₂ O ₅(0.89) %
ECRM-D 827-1	Thomas phosphate - powder Year of issue: 1976 Certified values	100 g
	SiO ₂6.21 ± 0.15 % Al ₂ O ₃(0.57) % CaO.....47.38 ± 0.49 % MgO.....(3.70) % (Values in parenthesis are indicative values)	P ₂ O ₅20.70 ± 0.16 % P ₂ O ₅ citric acid sol....18.79 ± 0.22 % Cr ₂ O ₃(0.14) % Fe _{total}(15.72) % Mn _{total}(2.34) % V ₂ O ₅(1.15) %
ECRM-B 851-1	Basic slag - powder Year of issue: 1979 Certified values	100 g
	Si.....6.09 % Al.....2.0 % Ca.....28.66 % Mg.....2.25 % P.....1.33 % B.....(0.0144) % Cr.....0.55 % (Values in parenthesis are indicative values)	F.....0.10 % Fe.....19.9 % K.....0.0089 % Mn.....6.17 % Mo.....(0.0050) % Na.....0.061 % Ni.....(0.0037) % Pb.....(0.0028) % S.....0.37 % Ti.....0.25 % V.....0.13 % Zn.....0.0141 %
ECRM-B 877-1	Furnace dust - powder Year of issue: 1977 Certified values	100 g
	Si.....1.08 ± 0.05 % Al.....0.044 ± 0.003 % Ca.....3.23 ± 0.11 % Mg.....0.28 ± 0.01 % P.....0.18 ± 0.01 % As.....0.014 ± 0.002 % C.....0.83 ± 0.02 % (Values in parenthesis are indicative values)	Cr.....0.017 ± 0.001 % F.....0.78 ± 0.05 % Fe.....62.07 ± 0.35 % K.....0.058 ± 0.003 % Mn.....1.37 ± 0.02 % Na.....0.23 ± 0.02 % Ni.....0.010 ± 0.001 % Pb.....1.00 ± 0.03 % S.....0.18 ± 0.01 % Ti.....0.032 ± 0.003 % V.....0.029 ± 0.003 % Zn.....1.16 ± 0.04 %

Iron and steel products

Code	Product	Unit
ECRM-B 879-1	Basic slag - powder Year of issue: 1980 Certified values Si 4.12 ± 0.03 % SiO ₂ 8.82 % Al 0.425 ± 0.016 % Al ₂ O ₃ 0.803 % Ca 31.23 ± 0.26 % CaO 43.70 % Mg 1.32 ± 0.03 % MgO 2.19 % P 3.69 ± 0.03 % P citric acid sol. 3.31 ± 0.10 % P ₂ O ₅ 8.46 % P ₂ O ₅ citric acid sol. 7.59 % Cr 0.326 ± 0.013 % Cr ₂ O ₃ 0.477 % F 0.368 ± 0.035 % Fe 18.97 ± 0.19 % Mn 3.45 ± 0.04 % MnO 4.45 % S 0.102 ± 0.008 % Ti 0.321 ± 0.005 % TiO ₂ 0.535 % V 0.414 ± 0.012 % V ₂ O ₅ 0.738 %	100 g
	(Values in parenthesis are indicative values)	
BCS CRM381	Basic slag Certified values SiO ₂ 8.78 % TiO ₂ 0.35 % Al ₂ O ₃ 0.67 % Fe 13.3 % FeO 3.69 % CaO 49 % MgO 1.03 % Cr ₂ O ₃ 0.33 % MnO 3.16 % V ₂ O ₅ 0.94 % P ₂ O ₅ Cit Sol 15.2 % P ₂ O ₅ Total 15.7 % S 0.19 %	100 g
ECRM-F 802-1	Blast furnace slag - powder Year of issue: 1971 Certified values Si 15.16 ± 0.14 % Al 8.53 ± 0.31 % Ca 30.62 ± 0.27 % Mg 2.87 ± 0.17 % P 0.109 ± 0.011 % B 0.0245 ± 0.003 %* Cr 0.0053 ± 0.0004 %* F 0.243 ± 0.008 %* Fe 0.576 ± 0.021 % K 0.491 ± 0.006 %* Mn 0.460 ± 0.016 % Na 0.236 ± 0.004 %* Pb (0.0015) % S 0.714 ± 0.036 % Ti 0.366 ± 0.018 % V 0.028 ± 0.002 %* Zn 0.0025 ± 0.0005 %*	100 g
	(Values in parenthesis are indicative values) * 95%-confidence interval	
ECRM-F 803-1	Blast furnace slag - powder Year of issue: 1971 Certified values Si 17.01 ± 0.13 % Al 6.98 ± 0.28 % Ca 30.93 ± 0.25 % Mg 2.44 ± 0.19 % P 0.118 ± 0.011 % Fe 0.613 ± 0.025 % Mn 0.552 ± 0.004 % S 0.767 ± 0.044 % Ti 0.301 ± 0.017 %	100 g
ECRM-F 804-1	Basic slag - powder Year of issue: 1971 Certified values Si 2.59 ± 0.09 % Al (0.407) % Ca 36.88 ± 0.34 % Mg 0.88 ± 0.15 % P 7.67 ± 0.10 % Fe 11.92 ± 0.27 % Mn 1.48 ± 0.04 % S 0.127 ± 0.015 % Ti 0.152 ± 0.014 % V 0.460 ± 0.024 %	100 g
	(Values in parenthesis are indicative values)	
ECRM-F 805-1	Basic slag - powder Year of issue: 1972 Certified values Si 3.10 ± 0.04 % Al 0.326 ± 0.013 % Ca 34.96 ± 0.12 % Mg 1.12 ± 0.05 % P 7.07 ± 0.05 % Fe 14.87 ± 0.07 % Mn 1.59 ± 0.03 % S 0.092 ± 0.005 % Ti 0.205 ± 0.008 % V 0.514 ± 0.009 %	100 g
ECRM-F 806-1	Basic slag - powder Year of issue: 1972 Certified values Si 5.48 ± 0.06 % Al 0.477 ± 0.020 % Ca 32.97 ± 0.19 % Mg 1.82 ± 0.06 % P 0.982 ± 0.022 % Fe 17.89 ± 0.10 % Mn 4.60 ± 0.06 % S 0.110 ± 0.004 % Ti 0.302 ± 0.008 % V 0.288 ± 0.009 %	100 g
	(Values in parenthesis are indicative values)	

Code	Product	Unit	
ECRM-F 876-1	Electric furnace dust - powder	100 g	
	Year of issue: 1977		
	Certified values		
	Si.....1.72 ± 0.04 %	C.....0.26 ± 0.01 %	Na.....1.98 ± 0.07 %
	SiO ₂3.67 %	Cd.....0.13 ± 0.01 %	Na ₂ O.....2.67 %
	Al.....0.34 ± 0.03 %	Cl.....3.63 ± 0.08 %	Ni.....0.034 ± 0.003 %
	Al ₂ O ₃0.63 %	Cr.....0.17 ± 0.01 %	Pb.....7.82 ± 0.23 %
	Ca.....3.43 ± 0.13 %	F.....0.24 ± 0.03 %	S.....0.87 ± 0.07 %
	CaO.....4.80 %	Fe.....24.85 ± 0.17 %	Sn.....0.094 ± 0.004 %
	Mg.....1.31 ± 0.03 %	K.....1.63 ± 0.06 %	Ti.....0.048 ± 0.004 %
	MgO.....2.17 %	K ₂ O.....1.96 %	TiO ₂0.080 %
	P.....0.128 ± 0.007 %	Mn.....2.84 ± 0.04 %	Zn.....23.29 ± 0.32 %
	As.....0.023 ± 0.001 %	MnO.....3.67 %	
	(Values in parenthesis are indicative values)		
ECRM-F 878-1	Blast furnace slag - powder	100 g	
	Year of issue: 1977		
	Certified values		
	Si.....15.77 ± 0.08 %	MgO.....9.62 %	Mn.....0.99 ± 0.01 %
	SiO ₂33.73 %	P.....0.013 ± 0.001 %	MnO.....1.28 %
	Al.....8.45 ± 0.10 %	P ₂ O ₅0.029 %	Na.....0.35 ± 0.02 %
	Al ₂ O ₃15.97 %	F.....0.15 ± 0.01 %	Na ₂ O.....0.47 %
	Ca.....25.46 ± 0.13 %	Fe.....0.60 ± 0.03 %	S.....0.83 ± 0.02 %
	CaO.....35.63 %	K.....1.07 ± 0.03 %	Ti.....0.37 ± 0.01 %
	Mg.....5.80 ± 0.07 %	K ₂ O.....1.29 %	TiO ₂0.62 %
	(Values in parenthesis are indicative values)		

Steel with oxygen and nitrogen content

ECRM-D 026-1	Rolled wire for the determination of oxygen and nitrogen - rods	100 g
	Year of issue: 1969	
	Certified values	
	O.....0.0031 ± 0.0003 %	N.....0.0053 ± 0.0004 %
ECRM-D 026-1-S	Rolled wire for the determination of oxygen and nitrogen - set rods	5 x 100 g
ECRM-D 026-2	Rolled wire for the determination of oxygen and nitrogen - rod	100 g
	Year of issue: 1973	
	Certified values	
	O.....0.0025 ± 0.0004 %	N.....0.0042 ± 0.0003 %
ECRM-D 026-2-S	Rolled wire for the determination of oxygen and nitrogen - set rod	5 x 100 g
ECRM-D 027-1	Rolled wire for the determination of oxygen and nitrogen - rod	100 g
	Year of issue: 1970	
	Certified values	
	O.....0.0084 ± 0.0006 %	N.....0.0157 ± 0.0010 %
ECRM-D 027-1-S	Rolled wire for the determination of oxygen and nitrogen - set rod	5 x 100 g
ECRM-D 028-1	Rolled wire for the determination of oxygen and nitrogen - rod	100 g
	Year of issue: 1970	
	Certified values	
	O.....0.0113 ± 0.0007 %	N.....0.0029 ± 0.0005 %
ECRM-D 028-1-S	Rolled wire for the determination of oxygen and nitrogen - set rod	5 x 100 g
ECRM-D 029-1	Rolled wire for the determination of oxygen and nitrogen - rod	100 g
	Year of issue: 1970	
	Certified values	
	O.....0.0312 ± 0.0010 %	N.....0.0083 ± 0.0008 %
ECRM-D 029-1-S	Rolled wire for the determination of oxygen and nitrogen - set rod	5 x 100 g
ECRM-D 099-1	Ball bearing steel for the determination of oxygen and nitrogen	100 g
	Year of issue: 1987	
	Certified values	
	O.....0.0008 ± 0.0002 %	N.....0.0078 ± 0.0005 %

Non ferrous metals and alloys

Code	Product	Unit
ECRM-D 284-2	Stainless steel, 17% chromium, 11% nickel and 2% molybdenum, no. 1.4571 - chips Year of issue: 2000 Certified values C 0.0201 ± 0.0005 %* Al 0.0027 ± 0.0004 %* Ti 0.191 ± 0.004 %* Si 0.537 ± 0.008 %* As 0.0063 ± 0.0003 %* V 0.0425 ± 0.0016 %* Mn 1.745 ± 0.009 %* B 0.0026 ± 0.0001 %* W (0.0183) % P 0.0258 ± 0.0008 %* Co 0.0525 ± 0.0011 %* Zr (0.0005) % S 0.0237 ± 0.0005 %* Cu 0.1831 ± 0.0014 %* O 0.0099 ± 0.0007 ¹ % Cr 16.811 ± 0.019 %* N 0.0151 ± 0.0002 %* Ta (0.0013) % Mo 2.111 ± 0.010 %* Nb (0.0028) % Ni 10.72 ± 0.05 %* Sn 0.0047 ± 0.0002 %* <p>(Values in parenthesis are indicative values)</p> <p>* 95%-confidence interval</p> <p>¹ Oxygen content only certified for the chips</p>	100 g
ECRM-D 286-1	Stainless steel, 18% chromium, 8.5% nickel, 0.3% molybdenum and 0.3% sulfur, no. 1.4305 - chips Year of issue: 1985 Certified values C 0.100 ± 0.005 Mo 0.329 ± 0.009 % N 0.043 ± 0.002 % Mn 1.92 ± 0.03 % Ni 8.54 ± 0.04 % Pb (0.0003) % P 0.026 ± 0.002 % Al (0.0023) % Sn 0.0084 ± 0.0009 % S 0.280 ± 0.014 % B (0.0003) % O (0.0315) % Cr 18.13 ± 0.08 % Co 0.150 ± 0.008 % Sb 0.0014 ± 0.0004 % <p>(Values in parenthesis are indicative values)</p>	100 g
ECRM-S 274-1	Vanadium steel - chips Year of issue: 2004 Certified values C 1.563 ± 0.006 %* Ni 0.077 ± 0.002 %* Pb (0.0064) % Si 1.057 ± 0.008 %* Al (0.0025) % Sn (0.010) % Mn 0.397 ± 0.004 %* As (0.0013) % Ti (0.0011) % P 0.0148 ± 0.0005 %* B (0.0005) % V 4.010 ± 0.018 %* S 0.0096 ± 0.0004 %* Co (0.0229) % W 0.0087 ± 0.0007 %* Cr 8.036 ± 0.022 %* Cu 0.0281 ± 0.0005 %* O 0.0026 ± 0.0002+1) % Mo 1.455 ± 0.016 %* N 0.0769 ± 0.0010 %* Sb (0.00019) % <p>(Values in parenthesis are indicative values)</p> <p>* 95%-confidence interval</p>	100 g
ECRM-S 274-1-D	Vanadium steel - disc (38 mm x 25 mm)	disc

Setting-up sample for spectrometric analysis of low alloyed steels

BAM SUS-1 R	Setting-up sample for spectrometric analysis of low alloyed steels	cylinder
Uncertified mass fractions		
C	0.9 %	Cr 1.7 %
Si	0.8 %	Mo 0.9 %
Mn	1.1 %	Ni 2.9 %
P	0.02 %	V 0.5 %
S	0.017 %	W 0.7 %

Non ferrous metals and alloys

Aluminium base standards

Code	Product	Unit
BAM-M315	Aluminium alloy AlSi9Cu3	disc
	Year of issue: 2006	
	Certified values	
Si	9.18 ± 0.21 %	B (< 3) µg/g
Fe	0.59 ± 0.02 %	Be 5 ± 2 µg/g
Cu	2.51 ± 0.09 %	Bi 41 ± 7 µg/g
Mn	0.314 ± 0.007 %	Ca (~ 15) µg/g **)
Mg	0.422 ± 0.012 %	Cd 11 ± 4 µg/g
Cr	0.0311 ± 0.0007 %	Co (< 3) µg/g
Ni	0.096 ± 0.003 %	Ga 101 ± 5 µg/g
Zn	0.77 ± 0.02 %	Hg (33 ± 2) µg/g
Ti	0.143 ± 0.005 %	Li (~ 7) µg/g **)
	(Values in parenthesis are indicative values)	
	** the given values are average values. the exact value must be calculated for each single sample	

Code	Product	Unit
BAM-307	Aluminium alloy AlMg4,5Mn - disc (60 mm x 25 mm) Year of issue: 1990 Certified values Si.....0.155 ± 0.005 % Cr.....0.162 ± 0.003 % Cd.....0.00489 ± 0.00009 % Fe.....0.412 ± 0.004 % Zn.....0.0634 ± 0.0006 % Li.....0.00044 ± 0.00003 % Cu.....0.1043 ± 0.0012 % Ti.....0.1009 ± 0.0012 % Na.....0.00214 ± 0.00008 % Mn.....0.701 ± 0.004 % Be.....0.0011 ± 0.00003 % Mg.....4.576 ± 0.021 % Ca.....0.00053 ± 0.00005 %	disc
BAM-308	Aluminium alloy AlZnMgCu1,5 - disc (60 mm x 25 mm) Year of issue: 1990 Certified values Si.....0.0707 ± 0.0024 % Mg.....2.290 ± 0.013 % Ti.....0.0285 ± 0.0009 % Fe.....0.1634 ± 0.0027 % Cr.....0.1962 ± 0.0024 % Be.....0.00022 ± 0.00001 % Cu.....1.315 ± 0.011 % Ni.....0.0122 ± 0.0004 % Zr.....0.0078 ± 0.0004 % Mn.....0.0342 ± 0.0009 % Zn.....5.67 ± 0.04 %	disc
BAM-310	Aluminium alloy Al99,85Mg1 - disc (60 mm x 25 mm) Year of issue: 1993 Certified values Si.....0.0797 ± 0.0012 % Ti.....0.00301 ± 0.00011 % Na.....(0.0003) % Fe.....0.0705 ± 0.0012 % B.....(0.0006) % P.....(0.0003) % Cu.....0.00169 ± 0.00009 % Be.....0.000128 ± 0.000014 % Pb.....0.00347 ± 0.00025 % Mn.....0.00307 ± 0.00011 % Ca.....0.00073 ± 0.00004 % Sn.....0.00238 ± 0.00018 % Mg.....0.994 ± 0.015 % Cd.....0.00237 ± 0.00007 % V.....0.00444 ± 0.00023 % Cr.....0.00090 ± 0.00012 % Co.....(0.0009) % Zr.....0.00135 ± 0.00019 % Ni.....0.00244 ± 0.00014 % Ga.....0.01152 ± 0.00024 % Zn.....0.0086 ± 0.0004 % Li.....0.000366 ± 0.000012 %	disc
BAM-311	(Values in parenthesis are indicative values)	
BAM-311	Aluminium alloy AlCuMg2 - disc (60 mm x 25 mm) Year of issue: 1993 Certified values Si.....0.2040 ± 0.0029 % Zn.....0.2005 ± 0.0022 % Ga.....0.0159 ± 0.0005 % Fe.....0.310 ± 0.006 % Ti.....0.0562 ± 0.0006 % Li.....0.00053 ± 0.00005 % Cu.....4.653 ± 0.028 % Be.....0.00052 ± 0.00004 % Na.....(0.0018) % Mn.....0.694 ± 0.006 % Bi.....0.0500 ± 0.0030 % Pb.....0.0504 ± 0.0011 % Mg.....1.567 ± 0.014 % Ca.....(0.0006) % Sn.....0.0127 ± 0.0012 % Cr.....0.1037 ± 0.0014 % Cd.....0.00127 ± 0.00005 % V.....0.0240 ± 0.0008 % Ni.....0.0519 ± 0.0009 % Co.....0.00115 ± 0.00010 % Zr.....0.140 ± 0.005 %	disc
BAM-311	(Values in parenthesis are indicative values)	
BAM-312	Aluminium alloy AlMgSi0,5 - disc (60 mm x 25 mm) Year of issue: 1995 Certified values Si.....0.415 ± 0.006 % Ni.....0.00452 ± 0.00015 % Pb.....0.00439 ± 0.00025 % Fe.....0.185 ± 0.004 % Zn.....0.0290 ± 0.0004 % Sn.....(0.002) % Cu.....0.0419 ± 0.0008 % Ti.....0.0288 ± 0.0004 % Sr.....0.00082 ± 0.00010 % Mn.....0.0416 ± 0.0008 % Bi.....0.0023 ± 0.0004 % V.....0.00615 ± 0.00023 % Mg.....0.410 ± 0.005 % Cd.....0.00226 ± 0.00010 % Zr.....0.00101 ± 0.00005 % Cr.....0.0276 ± 0.0008 % Ga.....0.0115 ± 0.0004 %	disc
BAM-312	(Values in parenthesis are indicative values)	
BAM-314	Aluminium alloy AlSi11Cu2(Fe) - disc (60 mm x 25 mm) Year of issue: 1999 Certified values Si.....11.49 ± 0.10 % Zn.....1.195 ± 0.012 % Ga.....0.0154 ± 0.0006 % Fe.....0.757 ± 0.007 % Ti.....0.1638 ± 0.0025 % Pb.....0.221 ± 0.006 % Cu.....2.071 ± 0.0019 % As.....0.00279 ± 0.00020 % Sb.....0.0093 ± 0.0005 % Mn.....0.400 ± 0.003 % Be.....0.000396 ± 0.000021 % Sn.....0.199 ± 0.005 % Mg.....0.1805 ± 0.0029 % Bi.....0.0094 ± 0.0004 % V.....0.0192 ± 0.0005 % Cr.....0.0517 ± 0.0008 % Cd.....0.00130 ± 0.00005 % Zr.....0.00552 ± 0.00013 % Ni.....0.221 ± 0.003 % Co.....0.00532 ± 0.00021 %	disc
BAM-314	(Values in parenthesis are indicative values)	
ERM-EB313	Aluminium alloy AlMg3 - disc (60 mm x 25 mm) (BAM-313) Year of issue: 1997 Certified values Si.....0.363 ± 0.007 % As.....0.00072 ± 0.00007 % Na.....0.00370 ± 0.00024 % Fe.....0.391 ± 0.003 % Be.....0.000547 ± 0.00002 % Pb.....0.00433 ± 0.00028 % Cu.....0.0932 ± 0.0014 % Bi.....0.0095 ± 0.0008 % Sb.....0.00087 ± 0.00019 % Mn.....0.495 ± 0.003 % Ca.....0.00057 ± 0.00008 % Sn.....0.0197 ± 0.0006 % Mg.....3.40 ± 0.04 % Cd.....0.00074 ± 0.00004 % Tl.....0.00064 ± 0.00004 % Cr.....0.1224 ± 0.0012 % Ga.....0.0121 ± 0.00005 % V.....0.0299 ± 0.0006 % Ni.....0.0278 ± 0.0006 % Hg.....0.00041 ± 0.00004 % Zr.....0.0359 ± 0.0019 % Zn.....0.1579 ± 0.0015 % Li.....0.000604 ± 0.00001 % Ti.....0.0947 ± 0.0014 % Mo.....0.00053 ± 0.00012 %	disc
ERM-EB313	(Values in parenthesis are indicative values)	

Non ferrous metals and alloys

Code	Product	Unit
BAM 201	Aluminium alloy GAlSi12 - chips Year of issue: 1963 Certified values Si 13.20 % Fe 0.18 % Zn 0.038 % Mg 0.0024 % Mn 0.38 % Cu 0.009 % Ti 0.011 %	100 g
BAM 209	Aluminium alloy GAlSi10Mg - chips Year of issue: 1963 Certified values Si 9.65 % Fe 0.18 % Zn 0.021 % Mg 0.31 % Mn 0.36 % Cu 0.004 % Ti 0.023 %	100 g
BAM 300	Aluminium alloy AlMg3 - chips Year of issue: 1959 Certified values Si 0.14 % Mn 0.018 % Ti 0.011 % Mg 2.67 % Cr 0.23 % Zn 0.128 % Cu 0.046 % Pb 0.016 % Fe 0.203 % Sn (< 0.0005) %	100 g
BAM 301	(Values in parenthesis are indicative values) Aluminium 99.8 - chips Year of issue: 1961 Certified values Si 0.061 % Fe 0.054 % Ti 0.005 % Mg 0.0008 % Mn 0.001 % V 0.0018 % Cu 0.0016 % Sn (< 0.0005) % Zn 0.033 %	100 g
	(Values in parenthesis are indicative values)	

Copper base standards

BAM-M381	Pure copper - disc (40 mm x 30 mm) Year of issue: 2006 Certified values Al (< 1) µg/g Bi < 0.3 µg/g Sb < 1 µg/g Ni 0.7 ± 0.2 µg/g Cd < 0.4 µg/g Se (< 1) µg/g Fe 3.3 ± 0.2 µg/g Co < 0.3 µg/g Si (< 3) µg/g Mn 0.22 ± 0.03 µg/g Cr < 0.4 µg/g Sn 3.86 ± 0.25 µg/g Zn 5.3 ± 0.3 µg/g Mg < 0.6 µg/g Te (< 0.3) µg/g Ag < 1 µg/g Pb 0.59 ± 0.07 µg/g Ti (< 0.3) µg/g As < 0.5 µg/g S (3.2 ± 1.3) µg/g Zr < 6 µg/g	disc
BAM-M382	(Values in parenthesis are indicative values) Pure copper - disc (40 mm x 30 mm) Year of issue: 2006 Certified values Al < 2.5 µg/g Bi 0.53 ± 0.03 µg/g Sb 0.7 ± 0.2 µg/g Ni 1.7 ± 0.2 µg/g Cd 0.90 ± 0.09 µg/g Se 0.6 ± 0.1 µg/g Fe 6.0 ± 0.4 µg/g Co 0.73 ± 0.07 µg/g Si < 6 µg/g Mn 0.76 ± 0.06 µg/g Cr 0.56 ± 0.06 µg/g Sn 4.29 ± 0.21 µg/g Zn 6.0 ± 0.5 µg/g Mg (1.4 ± 0.3) µg/g Te 0.61 ± 0.06 µg/g Ag 1.8 ± 0.2 µg/g Pb 1.0 ± 0.2 µg/g Ti (0.6 ± 0.2) µg/g As (0.6 ± 0.2) µg/g S (3.2 ± 1.4) µg/g Zr < 3 µg/g	disc
BAM-M383	(Values in parenthesis are indicative values) Pure copper - disc (40 mm x 30 mm) Year of issue: 2005 Certified values Al (2.3 ± 0.6) µg/g Bi 1.02 ± 0.09 µg/g Sb 1.44 ± 0.17 µg/g Ni 3.59 ± 0.21 µg/g Cd 1.48 ± 0.15 µg/g Se (1.16 ± 0.19) µg/g Fe 10.9 ± 0.5 µg/g Co 1.37 ± 0.05 µg/g Si < 10 µg/g Mn 1.24 ± 0.05 µg/g Cr 1.03 ± 0.09 µg/g Sn 4.7 ± 0.6 µg/g Zn (7.8 ± 0.4) µg/g Mg 2.37 ± 0.29 µg/g Te 1.40 ± 0.16 µg/g Ag 4.70 ± 0.20 µg/g Pb 1.31 ± 0.20 µg/g Ti 1.56 ± 0.16 µg/g As 1.93 ± 0.15 µg/g S (2.8 ± 1.4) µg/g	disc

Code	Product	Unit
BAM-M384	Pure copper - disc (40 mm x 30 mm) Year of issue: 2005 Certified values Al 13.0 ± 0.8 µg/g Bi 3.34 ± 0.22 µg/g Sb 12.0 ± 0.4 µg/g Ni 5.7 ± 0.4 µg/g Cd 3.95 ± 0.09 µg/g Se 4.24 ± 0.19 µg/g Fe 32.8 ± 1.9 µg/g Co 3.88 ± 0.16 µg/g Si (5.0 ± 0.7) µg/g Mn 6.88 ± 0.15 µg/g Cr 6.53 ± 0.21 µg/g Sn (10.2 ± 0.9) µg/g Zn (12.7 ± 2.1) µg/g Mg 14.6 ± 0.5 µg/g Te 7.0 ± 0.5 µg/g Ag 10.3 ± 0.4 µg/g Pb 5.7 ± 0.5 µg/g Ti (2.10 ± 0.23) µg/g As 5.0 ± 0.4 µg/g S (4.1 ± 1.0) µg/g Zr < 9 µg/g (Values in parenthesis are indicative values)	disc
BAM-229	Copper-alloy CuZn37 - chips Year of issue: 1996 Certified values Cu 63.334 ± 0.007% Fe 106.1 ± 2.1 µg/g Sb 7.2 ± 0.7 µg/g Zn 36.63 ± 0.04% Ni 111.4 ± 0.9 µg/g Se 34 ± 4 µg/g Sn 48.5 ± 1.1 µg/g As 21.7 ± 0.8 µg/g Pb 192 ± 5 µg/g P (10.6 ± 1.6) µg/g (Values in parenthesis are indicative values)	100 g
BAM-365	Refined copper - chips Year of issue: 1996 Certified values Cu 99.937 ± 0.012% Mn (< 1) µg/g S (7.7 ± 1.4) µg/g Sn (< 5) µg/g Ag 102.7 ± 1.7 µg/g Sb 8.8 ± 0.3 µg/g Pb 28.8 ± 1.3 µg/g As 29.8 ± 1.0 µg/g Te 4.6 ± 0.6 µg/g Fe 22.3 ± 1.3 µg/g Bi 29.4 ± 1.4 µg/g Ni 175.3 ± 1.5 µg/g Co 23.6 ± 1.4 µg/g (Values in parenthesis are indicative values)	100 g
BAM-366	SF-copper - chips Year of issue: 1992 Certified values Zn 15.6 ± 1.2 % Ag 7.9 ± 0.8 % S 8.7 ± 0.6 % Sn 111 ± 3 % As 1.11 ± 0.08 % Sb 0.99 ± 0.10 % Pb 10.8 ± 0.5 % Bi (< 0.3) % Se (< 1.1) % Fe 23.4 ± 0.5 % Cd 0.27 ± 0.04 % Te (< 0.3) % Ni 3.2 ± 0.7 % P 263 ± 6 % (Values in parenthesis are indicative values)	100 g
BAM-366-D	SF-copper	disc
BAM-367	Copper-alloy CuNi10Fe1Mn - disc (40 mm x 30 mm) Year of issue: 1995 Certified values Cu 87.88 ± 0.04% Zn 715 ± 9 % P 124 ± 6 % Ni 9.72 ± 0.05% C 28.7 ± 0.6 % Pb 298 ± 6 % Fe 1.443 ± 0.012% Co 498 ± 3 % S 162 ± 9 % Mn 0.723 ± 0.005% Mg 347 ± 13 % Sn 105 ± 4 %	disc
BAM-368	Copper-alloy CuZn20Al2- disc (40 mm x 30 mm) Year of issue: 1993 Certified values Cu 77.049 ± 0.018% Mn 202.8 ± 2.4 µg/g Pb 131.3 ± 2.4 µg/g Al 1.972 ± 0.014% As 246 ± 9 µg/g S (18.5 ± 2.9) µg/g Ni 258 ± 4 µg/g Mg 62.1 ± 1.5 µg/g Si 130 ± 7 µg/g Fe 192.7 ± 2.9 µg/g P 89.9 ± 1.6 µg/g Sn 147 ± 4 µg/g (Values in parenthesis are indicative values)	disc
BAM-369	OF-copper- disc (40 mm x 30 mm) Year of issue: 1993 Certified values Zn 22.0 ± 0.6 % Co 10.42 ± 0.29 % Mg 3.60 ± 0.18 % Bi 9.7 ± 0.4 % Cr 9.2 ± 0.5 %	disc
BAM-370	OF-copper- disc (40 mm x 30 mm) Year of issue: 1993 Certified values Al 12.6 ± 0.8 % Pb 15.8 ± 1.1 % Si 18.7 ± 3.0 % P 11.7 ± 0.7 % Sb 15.6 ± 1.3 % Sn 16.8 ± 0.9 %	disc
BAM-371	OF-copper- disc (40 mm x 30 mm) Year of issue: 1995 Certified values Fe 18.3 ± 0.7 % Cd 1.63 ± 0.08 % Te 14.4 ± 0.6 % Be 11.5 ± 0.6 % S 12.1 ± 0.9 % Ti 12.9 ± 0.7 %	disc

Non ferrous metals and alloys

Code	Product	Unit
BAM-372	OF-copper- disc (40 mm x 30 mm) Year of issue: 1995 Certified values Ni.....11.66 ± 0.24 % Ag.....9.01 ± 0.29 % Se.....(8.4 ± 0.6) % Mn.....11.4 ± 0.4 % As.....10.3 ± 0.6 % Zr.....5.8 ± 0.4 % (Values in parenthesis are indicative values)	disc
BAM-373/1-3	E-copper with increasing P-contents - disc (50 mm x 30 mm) Set of BAM-373/1, BAM-373/2 and BAM-373/3 Certified values BAM-373/1 P.....33.8 ± 1.2 µg/g BAM-373/2 P.....226.5 ± 1.7 µg/g BAM-373/3 P.....455.7 ± 1.7 µg/g	set of discs
BAM-376	Pure copper - disc (40 mm x 30 mm) Year of issue: 1996 Certified values Al.....(181.5 ± 10) % Bi.....200 ± 5 % Sb.....202 ± 5 % Ni.....209 ± 6 % Cd.....186.1 ± 2.5 % Se.....210 ± 4 % Fe.....234.6 ± 2.7 % Co.....207.9 ± 1.8 % Sn.....247.3 ± 2.9 % Mn.....205.9 ± 2.5 % Cr.....(400 ± 9) % Te.....215 ± 7 % Zn.....217.3 ± 2.7 % Mg.....124 ± 4 % Ti.....(4.5 ± 1.7) % Ag.....163.0 ± 2.4 % P.....203 ± 5 % Zr.....42.2 ± 1.9 % As.....199.9 ± 2.5 % Pb.....236 ± 4 % Be.....40.6 ± 0.9 % S.....133 ± 5 % (Values in parenthesis are indicative values)	disc
ERM-EB374	Copper-alloy CuSn8 - cylinder (40 mm x 30 mm) (BAM-374) Year of issue: 1999 Certified values Cu92.22 ± 0.04% Bi(2.2 ± 1.3) µg/g Sb.....(6.3 ± 1.4) µg/g Al.....(< 1) µg/g Cd.....(< 1) µg/g Se.....(< 2) µg/g Ni.....32.7 ± 1.3 µg/g Co.....(< 1) µg/g Si.....(< 10) µg/g Fe.....40 ± 4 µg/g Cr.....(< 1) µg/g Sn.....7.60% ± 0.13% Mn.....4.3 ± 0.3 µg/g Mg.....(< 1) µg/g Te.....(< 1) µg/g Zn.....40.4 ± 1.9 µg/g P.....0.1697 ± 0.0023% Ti.....(< 1) µg/g Ag.....12.1 ± 1.3 µg/g Pb.....8.3 ± 0.9 µg/g Zr.....(< 1) µg/g As.....(4.3 ± 1.2) µg/g S.....(13 ± 5) µg/g (Values in parenthesis are indicative values)	cylinder
ERM-EB375	Copper-alloy CuZn39Pb3 - cylinder (40 mm x 30 mm) (BAM-375) Year of issue: 1999 Certified values Cu58.32 ± 0.05 % Ag166 ± 4 µg/g Pb.....2.90 ± 0.03 % Al.....270 ± 5 µg/g As231 ± 4 µg/g Sb.....122 ± 4 µg/g Ni.....0.1053 ± 0.0015 % Bi68.6 ± 2.5 µg/g Si.....211 ± 14 µg/g Fe.....0.207 ± 0.004 % Cd85.9 ± 2.1 µg/g Sn.....0.2090 ± 0.0024 % Mn.....222 ± 3 µg/g Co196.4 ± 2.8 µg/g Te.....53.8 ± 2.4 µg/g Zn.....38.02 ± 0.08 % P(8.6 ± 1.2) µg/g	cylinder
ERM-EB377	Copper-alloy CuSn6 - cylinder (40 mm x 30 mm) Year of issue: 1999 Certified values Cu94.04 ± 0.05 % Bi42.2 ± 1.5 µg/g Sb.....13.0 ± 1.3 µg/g Al.....45.1 ± 1.2 µg/g Cd.....(< 1) µg/g Se.....55 ± 4 µg/g Ni.....107.4 ± 1.5 µg/g Co.....(< 2) µg/g Si.....(134) µg/g Fe.....104.2 ± 2.7 µg/g Cr66.9 ± 2.1 µg/g Sn.....5.92 ± 0.13 % Mn.....92.1 ± 2.1 µg/g Mg.....(< 1) µg/g Te.....(< 1) µg/g Zn.....100.6 ± 3.0 µg/g P(< 10) µg/g Ti.....(< 1) µg/g Ag.....64.4 ± 1.1 µg/g Pb44.9 ± 2.3 µg/g As.....(< 10) µg/g S(6.8 ± 0.8) µg/g (Values in parenthesis are indicative values)	cylinder
ERM-EB378	Copper-alloy CuSn6 - cylinder (40 mm x 30 mm) (BAM-378) Year of issue: 2000 Certified values Cu94.13 ± 0.04 % Bi(< 1) µg/g Sb.....86.1 ± 2.6 µg/g Al.....(< 1) µg/g Cd.....100.7 ± 2.2 µg/g Se.....(< 2) µg/g Ni.....18.3 ± 0.9 µg/g Co89 ± 5 µg/g Si.....(< 10) µg/g Fe.....182 ± 7 µg/g Cr311 ± 5 µg/g Sn.....5.74 ± 0.21 % Mn.....(0.74 ± 0.24) µg/g Mg28.7 ± 0.8 µg/g Te.....85.0 ± 2.6 µg/g Zn.....(7.4 ± 1.0) µg/g P602 ± 23 µg/g Ti.....(29.4 ± 4) µg/g Ag.....26.6 ± 1.3 µg/g Pb4.2 ± 0.7 µg/g Zr.....(1.7 ± 0.09) µg/g As.....99.5 ± 2.5 µg/g S9.1 ± 1.9 µg/g (Values in parenthesis are indicative values)	cylinder

Code	Product	Unit
ERM-EB385	Pure copper - cylinder (40 mm x 30 mm) (BAM-M385) Year of issue: 2003 Certified values Al 28.6 ± 2.5 µg/g Ni 11.9 ± 0.8 µg/g Fe 45.4 ± 1.4 µg/g Mn 10.1 ± 0.2 µg/g Zn 57.9 ± 4.0 µg/g Ag 28.6 ± 0.8 µg/g As 11.4 ± 0.8 µg/g Bi 5.81 ± 0.17 µg/g Cd 5.8 ± 0.3 µg/g Co 6.93 ± 0.15 µg/g Cr 9.81 ± 0.20 µg/g Mg 29.1 ± 1.3 µg/g P 12.9 ± 1.0 µg/g Pb 11.3 ± 0.5 µg/g S 31.3 ± 1.5 µg/g Sb 19.9 ± 0.8 µg/g Se 7.2 ± 0.5 µg/g Si (7.2 ± 1.5) µg/g Sn 18.0 ± 0.9 µg/g Te 10.0 ± 0.4 µg/g Ti 3.83 ± 0.17 µg/g Zr (< 7) µg/g (Values in parenthesis are indicative values)	cylinder
ERM-EB386	Pure copper- cylinder (40 mm x 30 mm) (BAM-M386) Year of issue: 2003 Certified values Al 36.5 ± 2.5 µg/g Ni 25.0 ± 1.0 µg/g Fe 64.7 ± 1.8 µg/g Mn 13.3 ± 0.2 µg/g Zn 49.5 ± 1.6 µg/g Ag 47.4 ± 1.2 µg/g As 24.2 ± 1.0 µg/g Bi 9.6 ± 0.5 µg/g Cd 7.8 ± 0.4 µg/g Co 5.20 ± 0.14 µg/g Cr 12.4 ± 0.7 µg/g Mg 36.1 ± 1.2 µg/g P 7.2 ± 0.7 µg/g Pb 23.4 ± 1.2 µg/g S 21.9 ± 2.1 µg/g Sb 31.2 ± 1.1 µg/g Se 11.6 ± 0.3 µg/g Si (14.3 ± 4.3) µg/g Sn 28.3 ± 0.8 µg/g Te 38.3 ± 0.9 µg/g Ti 33.1 ± 1.3 µg/g Zr (8.9 ± 1.7) µg/g (Values in parenthesis are indicative values)	cylinder
ERM-EB387	Copper alloy CuZn20Ni5 - cylinder (40 mm x 30 mm) (BAM-M387) Year of issue: 2004 Certified values Cu 75.18 ± 0.04 % Ni 5.020 ± 0.025 % Fe 617 ± 10 µg/g Mn 796 ± 6 µg/g Zn 19.57 ± 0.06 % Pb 10.8 ± 0.8 µg/g Sn 30.1 ± 1.2 µg/g	cylinder
ERM-EB388	Copper alloy CuAl5Zn5Sn - cylinder (40 mm x 30 mm) (BAM-M388) Year of issue: 2004 Certified values Cu 89.27 ± 0.05 % Al 4.972 ± 0.024 % Ni 73.6 ± 2.0 µg/g Fe 303 ± 9 µg/g Mn 512 ± 6 µg/g Zn 4.81 ± 0.03 % Pb 9.69 ± 0.83 µg/g Sn 0.857 ± 0.011 %	cylinder
BAM 211	Copper-alloy G-SnBz10 Year of issue: 1974 Certified values Cu 87.71 ± 0.03 % Sn 10.60* ± 0.04 % Zn 0.56 ± 0.02 % Pb 0.74 ± 0.02 % Fe 0.110 ± 0.003 % Ni 0.122 ± 0.002 % Mn 0.0019 ± 0.0002 % Ag 0.059 ± 0.002 % As 0.0213 ± 0.0008 % Bi 0.0020 ± 0.0002 % Cd 0.00144 ± 0.00005 % P 0.0267 ± 0.0005 % S 0.0211 ± 0.0006 % Sb 0.033 ± 0.001 % Se 0.00114 ± 0.00005 %	100 g
BAM 223	Copper-alloy CuZn39Pb2 Year of issue: 1974 Certified values Cu 58.74 ± 0.02 % Sn 0.089 ± 0.004 % Zn 38.82 ± 0.09 % Pb 2.13 ± 0.02 % Fe 0.091 ± 0.002 % Ni 0.0214 ± 0.0005 % Mn (< 0.001) % Al (< 0.002) % As 0.0084 ± 0.0005 % Bi 0.0018 ± 0.0001 % P 0.0003 ± 0.00015 % S 0.0011 ± 0.0001 % Sb 0.0040 ± 0.0002 % Se (< 0.0001) % Si (< 0.003) % (Values in parenthesis are indicative values)	100 g
BAM 224	Copper-alloy CuZn40MnPb Year of issue: 1975 Certified values Cu 57.40 ± 0.02 % Sn 0.066 ± 0.003 % Zn 39.40 ± 0.04 % Pb 1.13 ± 0.04 % Fe 0.136 ± 0.002 % Ni 0.038 ± 0.001 % Mn 1.70 ± 0.03 % Al 0.0012 ± 0.0002 % As 0.0025 ± 0.0002 % Bi 0.0006 ± 0.0001 % P 0.0112 ± 0.0002 % S 0.0004 ± 0.0001 % Sb 0.0026 ± 0.0001 % Si (0.002) %	100 g
BAM 227	Copper-alloy Rg7 Year of issue: 1979 Certified values Cu 85.57 ± 0.03 % Sn 6.01 ± 0.07 % Zn 3.46 ± 0.03 % Pb 4.12 ± 0.04 % Fe 0.129 ± 0.002 % Ni 0.284 ± 0.003 % Al (< 0.0001) % As 0.081 ± 0.002 % Bi 0.0088 ± 0.0002 % P (0.0002) % S 0.122 ± 0.005 % Sb 0.160 ± 0.002 % Se 0.0028 ± 0.0002 % Si (< 0.01) % Te 0.0012 ± 0.0003 % (Values in parenthesis are indicative values)	100 g

Non ferrous metals and alloys

Code	Product	Unit
BAM 228	Copper-alloy Rg10 Year of issue: 1979 Certified values Cu 85.34 ± 0.03 % Ni 0.109 ± 0.005 % P 0.019 ± 0.001 % Sn 9.76 ± 0.05 % Mn (< 0.001) % S 0.036 ± 0.002 % Zn 3.32 ± 0.05 % Al (0.0001) % Sb 0.078 ± 0.001 % Pb 1.24 ± 0.03 % As 0.024 ± 0.001 % Se 0.0012 ± 0.0001 % Fe 0.036 ± 0.002 % Bi 0.0086 ± 0.0003 %	100 g
	(Values in parenthesis are indicative values)	
BCR-017A	Copper - Sulphur and phosphorous 42 mm diameter x 30 mm disc Certified values P 6.85 mg/kg S 10.4 mg/kg	370 g
BCR-017B	Copper - Sulphur and phosphorous Chips (bottle with 50 g)	50 g

Oxygen in copper

BAM-379/1	Oxygen in copper - disc (40 mm x 30 mm)	disc
	Certified value O 38 ± 4 µg/g	
BAM-379/2	Oxygen in copper - disc (40 mm x 30 mm)	disc
	Certified value O 212 ± 8 µg/g	
BAM-379/3	Oxygen in copper - disc (40 mm x 30 mm)	disc
	Certified value O 378 ± 12 µg/g	
BAM-379/1-3	Oxygen in copper - disc (40 mm x 30 mm) Set of BAM-379/1, BAM-379/2 and BAM-379/3 Certified values BAM-379/1 O 38 ± 4 µg/g BAM-379/2 O 212 ± 8 µg/g BAM-379/3 O 378 ± 12 µg/g	set of discs
BCR-022A	Copper (electrolytic tough pitch) - Oxygen 26 mm diameter x 9 mm disc Certified value O 138 mg/kg	42 g
BCR-022B	Copper (electrolytic tough pitch) - Oxygen 9 mm diameter x 50 mm rods	28 g
BCR-054R	Copper - Oxygen 7 mm diameter x 50 mm rod Certified value O 0.47 mg/kg	42 g
BCR-058	Copper (continuous cast) - Oxygen 7 mm diameter x 50 mm rods Certified value O 390 mg/kg	17 g

BCR-074A	Electrolytic copper - Trace elements Cylinder of 40 mm diameter and 30 mm height Certified values Ag 12.8 mg/kg Fe 1.14 mg/kg Se 0.37 mg/kg As 0.78 mg/kg Mn 1.27 mg/kg Sn <0.07 mg/kg Cd <0.02 mg/kg Ni 1.04 mg/kg Zn 0.46 mg/kg Co <0.05 mg/kg Pb 0.97 mg/kg Cr <0.1 mg/kg Sb 0.576 mg/kg	disc
BCR-074C	Indicative values for Bi, Te Electrolytic copper - Trace elements Small pieces (approx. 4 mm x 4 mm x 2 mm) in a bottle containing about 50 g	50 g

Code	Product	Unit
BCR-691	Copper alloys - Alloying elements Set of five discs (one of each composition) of 35 diameter and 2 mm thickness Certified values Substance Quaternary bronze Brass Arsenic-Copper Lead-bronze Tin-bronze As 1.94 g/kg 0.99 g/kg 46.0 g/kg 2.85 g/kg 1.94 g/kg Pb 79 g/kg 3.9 g/kg 1.75 g/kg 92 g/kg 2.04 g/kg Sn 71.6 g/kg 20.6 g/kg 2.02 g/kg 101 g/kg 70 g/kg Zn 60.2 g/kg 148 g/kg 0.55 g/kg 1.48 g/kg 1.57 g/kg	Set (5)

Unalloyed zinc

BAM-M601	Pure zinc- disc (45 mm x 30 mm) Year of issue:2005 Certified values Cd 0.55 ± 0.06 µg/g Tl 2.25 ± 0.09 µg/g In < 0.05 µg/g Fe 2.20 ± 0.09 µg/g Pb 15.7 ± 0.3 µg/g Cu 1.89 ± 0.11 µg/g Al < 0.5 µg/g (Values in parenthesis are indicative values)	disc
BCR-321	Unalloyed zinc - Trace elements Disc of 80 mm diameter and 20 mm thickness Certified values Al <0.7 mg/kg Pb 4.85 mg/kg Tl 0.78 mg/kg In <0.2 mg/kg Sn <0.5 mg/kg	0.7 kg
ERM-EB322	Unalloyed zinc - Trace elements Disc of 60 mm diameter and 30 mm thickness Certified values Cd 15.08 mg/kg Fe 19.1 mg/kg Sn 5.6 mg/kg Cu 5.89 mg/kg Pb 15.0 mg/kg Tl 5.28 mg/kg	0.7 kg
ERM-EB323	Unalloyed zinc - Trace elements Disc of 60 mm diameter and 30 mm thickness Certified values Cd 6.51 mg/kg Fe 11.3 mg/kg Sn 18.7 mg/kg Cu 18.9 mg/kg Pb 48.6 mg/kg Tl 10.8 mg/kg	0.7 kg
ERM-EB324	Unalloyed zinc - Trace elements Disc of 60 mm diameter and 30 mm thickness Certified values Cd 48.6 mg/kg Fe 58.5 mg/kg Sn 9.8 mg/kg Cu 9.87 mg/kg Pb 26.1 mg/kg Tl 19.9 mg/kg	0.7 kg
ERM-EB325	Unalloyed zinc - Trace elements Disc of 60 mm diameter and 30 mm thickness Certified values Cd 94.7 mg/kg Fe 56.1 mg/kg Sn 46.1 mg/kg Cu 47.5 mg/kg Pb 142 mg/kg Tl 36.8 mg/kg	0.7 kg
BCR-326	Unalloyed zinc - Trace elements Discs of 80 mm diameter and 20 mm thickness Certified values Cd 203.0 mg/kg Fe 264.8 mg/kg Cu 104.8 mg/kg Pb 307.0 mg/kg Indicative value for Al	0.7 kg
BCR-327	Unalloyed zinc - Trace elements Disc of 80 mm diameter and 20 mm thickness Certified values Cd 301.4 mg/kg Fe 144.0 mg/kg Pb 409.4 mg/kg Indicative value for Cu	0.7 kg

Zinc alloys

BCR-351	Zinc alloy ZnAl4 - Trace elements Disc of 80 mm diameter and 20 mm thickness Certified values Al 43.55 x 10 ³ mg/kg Mg 131.0 mg/kg Tl 0.74 mg/kg Cu 12.13 mg/kg Pb 4.50 mg/kg In <0.2 mg/kg Sn <1 mg/kg Indicative values for Cd, Ni	0.7 kg
---------	---	--------

Non ferrous metals and alloys

Code	Product	Unit
BCR-352	Zinc alloy ZnAl4 - Trace elements Disc of 80 mm diameter and 20 mm thickness Certified values Al..... 41.50×10^3 mg/kg In.....3.02 mg/kg Sn.....3.0 mg/kg Cd.....2.88 mg/kg Mg.....283.0 mg/kg Tl.....3.2 mg/kg Cu.....31.26 mg/kg Ni.....6.74 mg/kg Indicative value for Pb	0.7 kg
BCR-353	Zinc alloy ZnAl4 - Trace elements Disc of 80 mm diameter and 20 mm thickness Certified values Al..... 39.5×10^3 mg/kg In.....2.55 mg/kg Sn.....5.6 mg/kg Cd.....10.44 mg/kg Mg.....452.5 mg/kg Tl.....3.95 mg/kg Cu.....100.0 mg/kg Pb.....24.4 mg/kg	0.7 kg
BCR-354	Zinc alloy ZnAl4 - Trace elements Disc of 80 mm diameter and 20 mm thickness Certified values Al..... 37.27×10^3 mg/kg In.....9.8 mg/kg Pb.....30.8 mg/kg Cd.....29.7 mg/kg Mg.....602 mg/kg Sn.....14.1 mg/kg Cu.....312.3 mg/kg Ni.....83.1 mg/kg Tl.....11.01 mg/kg	0.7 kg
BCR-355	Zinc alloy ZnAl4 - Trace elements Disc of 80 mm diameter and 20 mm thickness Certified values Al..... 34.43×10^3 mg/kg In.....24.6 mg/kg Pb.....56.9 mg/kg Cd.....58.1 mg/kg Mg.....786 mg/kg Sn.....29.1 mg/kg Cu.....1035 mg/kg Ni.....268 mg/kg Tl.....23.25 mg/kg	0.7 kg
BCR-356	Zinc alloy ZnAl4Cu1 - Trace elements Disc of 80 mm diameter and 20 mm thickness Certified values Al..... 44.34×10^3 mg/kg Fe.....31.5 mg/kg Ni.....3.43 mg/kg Cd.....0.73 mg/kg In.....<0.2 mg/kg Pb.....9.87 mg/kg Cu..... 3.944×10^3 mg/kg Mg.....132.3 mg/kg Tl.....0.79 mg/kg Indicative value for Sn	0.7 kg
BCR-357	Zinc alloy ZnAl4Cu1 - Trace elements Disc of 80 mm diameter and 20 mm thickness Certified values Al..... 42.27×10^3 mg/kg In.....3.30 mg/kg Sn.....3.51 mg/kg Cd.....2.83 mg/kg Mg.....273 mg/kg Tl.....2.76 mg/kg Cu..... 5.849×10^3 mg/kg Ni.....9.82 mg/kg Fe.....25.7 mg/kg Pb.....13.8 mg/kg	0.7 kg
BCR-359	Zinc alloy ZnAl4Cu1 - Trace elements Disc of 80 mm diameter and 20 mm thickness Certified values Al..... 37.11×10^3 mg/kg In.....15.5 mg/kg Sn.....16.93 mg/kg Cd.....29.8 mg/kg Mg.....557 mg/kg Tl.....13.34 mg/kg Cu..... 9.89×10^3 mg/kg Ni.....92.6 mg/kg Fe.....119.7 mg/kg Pb.....36.2 mg/kg	0.7 kg
BCR-360	Zinc alloy ZnAl4Cu1 - Trace elements Disc of 80 mm diameter and 20 mm thickness Certified values Al..... 34.27×10^3 mg/kg In.....29.8 mg/kg Pb.....73.9 mg/kg Cd.....59.5 mg/kg Mg.....705 mg/kg Sn.....33.0 mg/kg Cu..... 12.34×10^3 mg/kg Ni.....267 mg/kg Tl.....25.9 mg/kg	0.7 kg
BCR-361	Zinc alloy ZnAl4Cu1 - Trace elements Disc of 80 mm diameter and 20 mm thickness Certified values Al..... 40.68×10^3 mg/kg Fe.....10.34 mg/kg Sn.....46.3 mg/kg Cu..... 7.98×10^3 mg/kg Pb.....5.31 mg/kg Tl.....37.4 mg/kg Indicative values for Cd, In	0.7 kg

Code	Product	Unit
Lead base standards		
BNM 010	Tin-lead solder - granulated powder Year of issue: 1991 Certified values Sn 63.40 ± 0.07 % Ni 0.0021 ± 0.0002 % Fe (0.0020) % Pb 36.47 ± 0.17 % Sb 0.0488 ± 0.0008 % In (< 0.001) % Bi 0.0245 ± 0.0010 % Ag (0.014) % Zn (< 0.0001) % Cd 0.0016 ± 0.0002 % As (0.012) % Cu 0.0417 ± 0.0014 % Au (< 0.001) % (Values in parenthesis are indicative values)	100 g
ERM-EB101 Lead alloy PbCaSnAl - disc (50 mm x 30 mm) (BAM-101) Year of issue: 1999 Certified values Ca 0.1436 ± 0.0016 % Al 0.0257 ± 0.0006 % Bi 0.0165 ± 0.0007 % Sn 0.293 ± 0.007 % Ag 0.00288 ± 0.00007 % Cu 0.00173 ± 0.00018 %		
ERM-EB102	Lead alloy PbCaSn - disc (50 mm x 30 mm) (BAM-102) Year of issue: 1999 Certified values Ca 0.0705 ± 0.0011 % Al 0.0124 ± 0.0004 % Bi 0.0148 ± 0.0005 % Sn 0.895 ± 0.011 % Ag 0.00248 ± 0.00007 % Cu 0.00109 ± 0.00007 %	disc
ERM-EB103	Lead-alloy, PbSb1.6 - disc (50 mm x 30 mm) (BAM-103) Year of issue: 2006 Certified values Sn 0.183 ± 0.026 %* As 0.097 ± 0.004 %* Ni 3.02 ± 0.27 mg/kg* Bi 0.0158 ± 0.0004 %* Se 0.0180 ± 0.0010 %* Cd 0.20 ± 0.08 mg/kg* Cu 9.7 ± 0.9 mg/kg* Ag 0.0066 ± 0.0006 %* S (5.4 ± 1.2 mg/kg*) Sb 1.64 ± 0.06 %* Tl 15.2 ± 0.7 mg/kg* Te (1.9 ± 0.6 mg/kg*) (Values in parenthesis are indicative values)	disc
BCR-286A	Lead (electrolytically refined) - Trace elements Block of 60 mm x 60 mm x 12 mm Certified values Bi 21.5 mg/kg Tl 2.47 mg/kg Sb 0.099 mg/kg Zn <0.1 mg/kg	600 g
BCR-286B	Lead (electrolytically refined) - Trace elements Chips in bottles containing about 160 g	160 g
BCR-287A	Lead (thermally refined) -Trace elements Block of 60 mm x 60 mm x 12 mm Certified values Ag 15.2 mg/kg Cu 0.98 mg/kg Tl 0.73 mg/kg Bi 67.3 mg/kg Sb 0.040 mg/kg Zn <0.1 mg/kg Cd 0.356 mg/kg Sn <0.05 mg/kg	600 g
BCR-287B	Lead (thermally refined) -Trace elements Chips in bottles containing about 160 g	160 g
BCR-288A	Lead with added impurities - Trace elements Block of 60 mm x 60 mm x 12 mm Certified values As 55.7 mg/kg Cu 19.3 mg/kg Tl 2.26 mg/kg Bi 215.8 mg/kg Se <0.2 mg/kg Zn 8.2 mg/kg Cd 33.3 mg/kg Te 32.8 mg/kg	600 g
BCR-288B	Lead with added impurities - Trace elements Chips in bottles containing about 160 g	160 g
BCR-055	Lead (refined) - Oxygen Disc of 30 mm diameter x 9 mm thickness Certified value O 1.0 mg/kg	71 g

Gold base standards

BAM-501	Potassium Dicyanoaurate (I) Year of issue: 1997 Certified value Au 682.23 ± 0.25 g/kg	6 g
---------	--	-----

Non ferrous metals and alloys

Code	Product	Unit
Nickel base standards		
BCR-099	Nickel - Oxygen and nitrogen 25 cubes of nominally 2 g each Certified values O 8.4 mg/kg N 1.1 mg/kg	50 g
Titanium base standards		
BCR-059A	Titanium alloy Ti6Al4V - Oxygen and nitrogen 26 mm diameter x 9 mm discs Certified values O 1750 mg/kg N 172 mg/kg	21 g
BCR-059B	Titanium alloy Ti6Al4V - Oxygen and nitrogen 25 cubes of 0.2 g	5 g
BCR-089	Titanium alloy Ti 6Al 4V - Aluminium, Vanadium Cylinder of 40 mm diameter and 20 mm height Certified values Al 59.7 g/kg V 39.76 g/kg	disc
BCR-090A	Titanium - Impurities Cylinder of 40 mm diameter and 20 mm height Certified values B 28.2×10^{-3} g/kg Cu 0.513 g/kg Mo 0.488 g/kg Co 0.501 g/kg Fe 0.563 g/kg Ni 0.667 g/kg Cr 0.533 g/kg Mn 0.314 g/kg Indivativ values for Nb, Sn, W, Zr	disc
BCR-090B	Titanium - Impurities Cubes of about 0.2 g	25 g
BCR-024B	Titanium - Oxygen and nitrogen 25 cubes of 0.4 g Certified values N 117 mg/kg O 608 mg/kg	10 g
BCR-024C	Titanium - Oxygen, nitrogen 25 cubes of 0.2 g	5 g
BCR-318	Titanium - Hydrogen 7 mm diameter x 1 mm discs (approx. 100 discs of at least 0.2 g each) Certified value H 12.2 mg/kg	20 g
Zirconium base standards		
BCR-275	Zirconium alloy Zircaloy-4 - Oxygen, nitrogen and carbon 10 discs each 13 mm diameter x 1 mm Certified values O 1670 mg/kg N 39.0 mg/kg C 113 mg/kg	9 g
BCR-276	Zirconium alloy Zircaloy-4 - Oxygen, nitrogen and carbon 4.5 mm diameter x 2 mm (approx. 100 cylinders) Certified values O 1540 mg/kg N 41.0 mg/kg C 108 mg/kg	20 g
Tungsten base standards		
BCR-102	Tungsten carbide powder - Oxygen Powder sealed under argon Certified value O 185 mg/kg	2-3 g

Bauxite and Fluorspar

Code	Product	Unit
BCS CRM395	Bauxite Certified values % Fe ₂ O ₃ 16.3 % Na ₂ O (0.02) % Pb (0.0028) % SiO ₂ 1.24 % K ₂ O (0.02) % Sr (0.0023) % Al ₂ O ₃ 52.4 % Cr (0.0453) % Zn (0.0043) % TiO ₂ 1.93 % Cu (0.0021) % L.O.I. 27.8 % CaO 0.05 % Mn (0.0042) % MgO 0.02 % Ni (0.0034) % (Values in parenthesis are indicative values)	100 g
NIST-69B	Bauxite, Arkansas - Constituents Certified values Al ₂ O ₃ 48.8 % MgO 0.085 % TiO ₂ 1.90 % CaO 0.13 % MnO 0.110 % V ₂ O ₅ 0.028 % Cr ₂ O ₃ 0.011 % P ₂ O ₅ 0.118 % ZnO 0.0035 % Fe ₂ O ₃ 7.14 % SO ₃ 0.551 % ZrO ₂ 0.29 % K ₂ O 0.068 % SiO ₂ 13.43 % LOI* 27.2 % Indicative values for Co, BaO, Na ₂ O LOI* Loss on Ignition	60 g
NIST-600	Bauxite, Australian - Constituents Certified values Al ₂ O ₃ 40.0 % MnO 0.013 % V ₂ O ₅ 0.060 % CaO 0.22 % Na ₂ O 0.022 % ZnO 0.003 % Cr ₂ O ₃ 0.024 % P ₂ O ₅ 0.039 % ZrO ₂ 0.060 % Fe ₂ O ₃ 17.0 % SO ₃ 0.155 % LOI* 20.5 % K ₂ O 0.23 % SiO ₂ 20.3 % MgO 0.055 TiO ₂ 1.31 % LOI* Loss on Ignition	90 g
NIST-696	Bauxite, Surinam - Constituents Certified value Al ₂ O ₃ 54.5 % MgO 0.012 % TiO ₂ 2.64 % CaO 0.018 % MnO 0.004 % V ₂ O ₅ 0.072 % Cr ₂ O ₃ 0.047 % P ₂ O ₅ 0.050 % ZnO 0.0014 % Fe ₂ O ₃ 8.70 % SiO ₂ 3.79 % ZrO ₂ 0.14 % K ₂ O 0.009 % SO ₃ 0.15 % LOI* 29.9 % Indicative values for Co, BaO, Na ₂ O LOI* Loss on Ignition	60 g
NIST-697	Bauxite, Dominican - Constituents Certified value Al ₂ O ₃ 45.8 % MgO 0.18 % TiO ₂ 2.52 % CaO 0.71 % MnO 0.41 % V ₂ O ₅ 0.063 % Cr ₂ O ₃ 0.100 % P ₂ O ₅ 0.97 % ZnO 0.037 % Fe ₂ O ₃ 20.0 % SiO ₂ 6.81 % ZrO ₂ 0.065 % K ₂ O 0.062 % SO ₃ 0.077 % LOI* 22.1 % Indicative values for Co, BaO, Na ₂ O LOI* Loss on Ignition	60 g
NIST-698	Bauxite, Jamaican - Constituents Certified value Al ₂ O ₃ 48.2 % MgO 0.058 % TiO ₂ 2.38 % CaO 0.62 % MnO 0.38 % V ₂ O ₅ 0.064 % Cr ₂ O ₃ 0.080 % P ₂ O ₅ 0.37 % ZnO 0.029 % Fe ₂ O ₃ 0.013 % SiO ₂ 0.69 % ZrO ₂ 0.0.061 % K ₂ O 0.010 % SO ₃ 0.143 % LOI* 27.3 % Indicative values for Co, BaO, Na ₂ O LOI* Loss on Ignition	60 g
BCS CRM392	Fluorspar Certified values SiO ₂ 0.67 % BaO 0.37 % Pb 0.18 % CaF ₂ 97.2 % S 0.12 % CaO 0.52 % CO ₂ 0.48 %	100 g
NIST-79A	Fluorspar, customs grade - Calcium fluoride Certified value CaF ₂ 297.39 %	120 g
NIST-180	Fluorspar, high grade - Calcium fluoride Certified value CaF ₂ 98.80%	120 g

Ores

Ores

Code	Product	Unit
BCS CRM308	Grecian chrome ore Certified values FeO 15.3 % TiO ₂ (0.16) % Na ₂ O (0.04) % SiO ₂ 4.25 % CaO 0.34 % K ₂ O (0.01) % Al ₂ O ₃ 19.4 % MgO 16.4 % L.O.I. (0.9) % Cr ₂ O ₃ 41.5 % MnO (0.14) % (Values in parenthesis are indicative values)	100 g
BCR-010	Tin ore powder - Tin Certified value Sn 765.9 g/kg	225 g
BCS CRM355	Tin ore Certified values Sn 31.42 % Zn 0.059 % Si 7.14 % Fe 17.08 % Pb 0.012 % Ti 0.37 % Cu 0.085 % S 0.5 % Al 4.12 % As 0.14 % W 0.35 % Ca 2.63 % Bi 0.015 % Ni 0.004 % F 2.07 %	100 g
BCR-109	Zinc ore concentrate - Trace elements Certified values Cd 4.61 g/kg Fe 145.1 g/kg Pb 7.38 g/kg Cu 9.46 g/kg Hg 0.00096 g/kg F 0.081 g/kg Mg 0.20 g/kg	200 g
BCR-110	Zinc ore concentrate - Trace elements Certified values Cd 10.51 g/kg Fe 5.46 g/kg Pb 97.8 g/kg Cu 16.28 g/kg Hg 0.1484 g/kg F 0.055 g/kg Mg 1.36 g/kg	200 g
NIST-181	Lithium ore, Spodumene - Lithium oxide Certified value Li ₂ O 6.39 %	45 g
NIST-182	Lithium ore, Petalite - Lithium oxide Certified value Li ₂ O 4.34 %	45 g
NIST-183	Lithium ore, Lepidolite - Lithium oxide Certified value Li ₂ O 4.12 %	45 g
NIST-330	Copper ore mill heads - Copper, molybdenum and rhenium Certified values Cu 0.84 % Mo 0.018 % Re 0.30 mg/kg Indicative values for Ag, Au	100 g
NIST-25D	Manganese ore - Constituents Certified values Mn 51.78 % K ₂ O 0.93 % TiO ₂ 0.13 % Al ₂ O ₃ 5.32 % P ₂ O ₅ 0.25 % Available oxygen 14.28 % Fe ₂ O ₃ 3.92 % SiO ₂ 2.52 % Indicative values for BaO, CaO, Moisture	100 g
BCS CRM176/2	Manganese ore Certified values Fe 6.86 % MgO 0.04 % K ₂ O 1.3 % SiO ₂ 2.53 % S 0.018 % BaO 0.19 % Al ₂ O ₃ 5.2 % P 0.087 % Pb 0.01 % TiO ₂ 0.3 % Mn 47.5 % As ₂ O ₃ 0.22 % CaO 0.09 % Na ₂ O 0.11 %	100 g
NCS DC47004	Manganese ore - Constituents Certified values Cu 0.013 % Zn 0.027 % MnO ₂ 67.25 % Fe 1.22 % Al ₂ O ₃ 2.20 % Na ₂ O 0.044 % Mn 45.39 % BaO 0.68 % SiO ₂ 16.16 % Ni 0.019 % CaO 1.06 % TiO ₂ 0.063 % P 0.054 % K ₂ O 1.00 % S 0.007 % MgO 0.64 % Indicative value for Mn(SiO ₃) ²⁻	100 g

Code	Product	Unit
NCS DC47005	Manganese ore - Constituents Certified values Cu 0.014 % Zn 0.029 % MnO ₂ 54.38 % Fe..... 2.24 % Al ₂ O ₃ 3.00 % Na ₂ O..... 0.048 % Mn 36.99 % BaO..... 0.47 % SiO ₂ 22.24 % Ni 0.019 % CaO..... 3.60 % TiO ₂ 0.10 % P 0.081 % K ₂ O..... 0.46 % S 0.013 % MgO..... 1.44 % Indicative value for Mn(SiO ₃ ²⁻)	100 g
NCS DC47006	Manganese ore - Constituents Certified values Cu 0.036 % Zn 0.064 % MnO ₂ 48.01 % Fe..... 11.24 % Al ₂ O ₃ 8.55 % Na ₂ O..... 0.039 % Mn 32.54 % BaO..... 0.18 % SiO ₂ 14.50 % Ni 0.099 % CaO..... 0.083 % TiO ₂ 0.43 % P 0.207 % K ₂ O..... 0.93 % S 0.019 % MgO..... 0.11 % Indicative value for Mn(SiO ₃ ²⁻)	100 g
NCS DC47007	Manganese ore - Constituents Certified values Cu 0.028 % Zn 0.048 % MnO ₂ 36.93 % Fe..... 20.99 % Al ₂ O ₃ 8.97 % Na ₂ O..... 0.030 % Mn 25.00 % BaO..... 0.23 % SiO ₂ 10.46 % Ni 0.073 % CaO..... 0.051 % TiO ₂ 0.54 % P 0.275 % K ₂ O..... 0.72 % S 0.032 % MgO..... 0.10 % Indicative value for Mn(SiO ₃ ²⁻)	100 g
NCS DC47008	Manganese ore - Constituents Certified values Cu 0.009 % Zn 0.018 % Mn(CO ₃ ²⁻)..... 22.46 % Fe..... 1.40 % Al ₂ O ₃ 1.68 % Na ₂ O 0.024 % Mn 22.54 % BaO..... 0.13 % SiO ₂ 14.07 % Ni 0.041 % CaO..... 14.73 % TiO ₂ 0.10 % P 0.043 % K ₂ O..... 0.46 % S 0.21 % MgO..... 3.50 % Indicative value for Mn(SiO ₃ ²⁻)	100 g
NCS DC47009	Manganese ore - Constituents Certified values Cu 0.01 % Zn 0.020 % Mn(CO ₃ ²⁻)..... 15.69 % Fe..... 2.07 % Al ₂ O ₃ 2.49 % Na ₂ O 0.040 % Mn 15.74 % BaO..... 0.15 % SiO ₂ 15.82 % Ni 0.050 % CaO..... 19.78 % TiO ₂ 0.15 % P 0.061 % K ₂ O..... 0.70 % S 0.27 % MgO..... 3.82 % Indicative value for Mn(SiO ₃ ²⁻)	100 g
NIST-RM 8602	Chinese lead ore - Elements and constituents This Reference Material (RM) is a lead ore mined from a skarn lead-zinc-copper deposit, then prepared and certified by the Hubei Geological Research Laboratory (HGRL) in China which is the sole authority for the information provided with this material including reference values and other technical information. It is intended for use in geological and geochemical investigations, particularly in geochemical exploration programs and in studies of ore genesis, and for environmental monitoring in mining areas. Each unit of this RM consists of 100 g of ore pulverized to pass a 100 µm screen (-160 mesh). Reference values for selected elements and gangue constituents are given.	100 g
NIST-RM 8603	Chinese lead ore - Elements and constituents This Reference Material (RM) is a lead ore mined from a skarn lead-zinc-copper deposit, then prepared and certified by the Hubei Geological Research Laboratory (HGRL) in China which is the sole authority for the information provided with this material including reference values and other technical information. It is intended for use in geological and geochemical investigations, particularly in geochemical exploration programs and in studies of ore genesis, and for environmental monitoring in mining areas. Each unit of this RM consists of 100 g of ore pulverized to pass a 100 µm screen (-160 mesh). Reference values for selected elements and gangue constituents are given.	100 g
NIST-277	Tungsten concentrate - Tungsten oxide Certified value WO ₃ 67.4 % Indicative values for Ca, Fe, Mn, Mo, Nb, O, P, Pb, S, Si, Sn, Ta, Ti	100 g
NIST-RM 8607	Chinese tungsten ore - Elements and constituents This Reference Material (RM) is a tungsten ore mined from a greisen tungsten-tin-bismuth deposit, then prepared and certified by the Hubei Geological Research Laboratory (HGRL) in China which is the sole authority for the information provided with this material including reference values and other technical information. It is intended for use in geological and geochemical investigations, particularly in geochemical exploration programs and in studies of ore genesis, and for environmental monitoring in mining areas. Each unit of this RM consists of 100 g of ore pulverized to pass a 100 µm screen (-160 mesh). Reference values for selected elements and gangue constituents are given.	100 g

Refractories and carbides

Code	Product	Unit
NIST-RM 8608	Chinese tungsten ore - Elements and constituents This Reference Material (RM) is a molybdenum ore mined from a skarn molybdenum deposit, then prepared and certified by the Hubei Geological Research Laboratory (HGRL) in China which is the sole authority for the information provided with this material including reference values and other technical information. It is intended for use in geological and geochemical investigations, particularly in geochemical exploration programs and in studies of ore genesis, and for environmental monitoring in mining areas. Each unit of this RM consists of 100 g of ore pulverized to pass a 100 µm screen (-160 mesh). Reference values for selected elements and gangue constituents are given	100 g
NCS DC70008	Tungsten ore - Constituents Certified values for 38 elements, 9 oxides Indicative values for 5 elements	100 g
NIST-886	Gold ore, refractory - Trace elements Certified values S.....1.466 % Al ₂ O ₃8.25 mg/kg Indicative values for C (total), SO ₃	200 g
NIST-1835	Borate ore - Constituents Certified values F.....0.348 % Fe ₂ O ₃1.141 % Na ₂ O3.484 % Al ₂ O ₃3.474 % MgO.....3.411 % SO ₃1.477 % B ₂ O ₃18.739 % MnO.....0.0333 % SrO.....0.9418 % BaO.....0.0497 % K ₂ O.....1.261 % TiO ₂0.1332 % CaO.....21.622 % SiO ₂18.408 % LOI*.....25.72 %	60 g

Refractories and carbides

Code	Product	Unit
BAM-S002	Tungsten metal powder Certified values Al.....29.4 ± 0.9 mg/kg Fe.....53 ± 5 mg/kg Na.....41 ± 5 mg/kg Ca.....46 ± 4 mg/kg K.....40.0 ± 1.8 mg/kg Ni.....29 ± 4 mg/kg Co.....45 ± 6 mg/kg Mg.....38.8 ± 2.7 mg/kg P.....(7.2 ± 1.3) mg/kg Cr.....47.0 ± 1.4 mg/kg Mn.....16.7 ± 1.9 mg/kg Si.....106 ± 10 mg/kg Cu.....28.4 ± 2.9 mg/kg Mo.....59 ± 4 mg/kg Sn.....42 ± 6 mg/kg (Values in parenthesis are indicative values)	100 g
NIST-76A	Burnt refractory (Al ₂ O ₃ -40%) - Constituents Certified values Al ₂ O ₃38.7 % Li ₂ O.....0.042 % SiO ₂54.9 % CaO.....0.22 % MgO.....0.52 % SrO.....0.037 % Fe ₂ O ₃1.60 % Na ₂ O.....0.07 % TiO ₂2.03 % K ₂ O.....1.33 % P ₂ O ₅0.120 % Indicative value for LOI* * Loss on Ignition	75 g
NIST-77A	Burnt refractory (Al ₂ O ₃ - 60%) - Constituents Certified values Al ₂ O ₃60.2 % Li ₂ O.....0.025 % SiO ₂35.0 % CaO.....0.05 % MgO.....0.38 % SrO.....0.009 % Fe ₂ O ₃1.00 % Na ₂ O.....0.037 % TiO ₂2.66 % K ₂ O.....0.090 % P ₂ O ₅0.092 % Indicative value for LOI* * Loss on Ignition	75 g
NIST-78A	Burnt refractory (Al ₂ O ₃ - 70%) - Constituents Certified values Al ₂ O ₃71.7 % Li ₂ O.....0.12 % SiO ₂19.4 % CaO.....0.11 % MgO.....0.70 % SrO.....0.25 % Fe ₂ O ₃1.2 % Na ₂ O.....0.078 % TiO ₂3.22 % K ₂ O.....1.22 % P ₂ O ₅1.3 % Indicative value for LOI* * Loss on Ignition	75 g
NIST-198	Silica brick - Constituents Certified values Al ₂ O ₃0.16 % Li ₂ O.....0.001 % P ₂ O ₅0.022 % CaO.....2.71 % MgO.....0.07 % TiO ₂0.02 % Fe ₂ O ₃0.66 % MnO.....0.008 % LOI*.....0.21 % K ₂ O.....0.017 % Na ₂ O.....0.012 %	45 g

* Loss on Ignition

Refractories and carbides

Code	Product	Unit
NIST-199	Silica brick - Constituents Certified values Al ₂ O ₃ 0.48 % Li ₂ O 0.002 % P ₂ O ₅ 0.015 % CaO 2.41 % Na ₂ O 0.015 % TiO ₂ 0.06 % Fe ₂ O ₃ 0.74 % MgO 0.13 % LOI* 0.17 % K ₂ O 0.094 % MnO 0.007 %	45 g
	* Loss on Ignition	
NIST-887	Cemented carbide (W83-Co25-Ti4) - Cobalt Certified value Co 10.35 % Indicative value for C	100 g
NIST-888	Cemented carbide (W64-Co25-Ti4) - Cobalt and tantalum Certified values Co 24.7 % Ta 4.77 % Indicative value for C	100 g
NIST-889	CEM carbide (W75-Co9-Ta5-Ti4) - Cobalt, tantalum and titanium Certified values Co 9.50 Ta 4.60 Ti 4.03 Indicative value for C	100 g
NIST-889	CEM carbide (W75-Co9-Ta5-Ti4) - Cobalt, tantalum and titanium Certified values Co 9.50 Ta 4.60 Ti 4.03 Indicative value for C	100 g
NIST-889	CEM carbide (W75-Co9-Ta5-Ti4) - Cobalt, tantalum and titanium Certified values Co 9.50 Ta 4.60 Ti 4.03 Indicative value for C	100 g
BCS CRM352/1	Tungsten carbide Certified values C _{Total} 6.154 % Fe 0.0029 % C _{Free} 0.036 % O 0.11 %	100 g
NIST-276B	Tungsten carbide - Carbon Certified value C 6.10 % Indicative values for C (free), O, N	75 g
BCS CRM359	Nitrogen bearing silicon carbide Certified values C _{Total} 23.46 % Mn (<0.01) % Ca 0.108 % Si _{Total} 67.6 % Ni (0.014) % Mg (<0.01) % Al _{Total} 0.118 % O (0.532) % Na (<0.01) % C _{Free} (0.061) % N (7.84) % K (<0.01) % Si _{Free} (0.325) % Ti 0.022 % Fe _{Total} 0.175 % V (0.027) % (Values in parenthesis are indicative values)	100 g
BCS CRM360	Sialon bonded silicon carbide Certified values C _{Total} 23.53 % Mn (<0.01) % Ca 0.115 % Si _{Total} 60.8 % Cr (<0.01) % Mg (<0.02) % Al _{Total} 6.52 % Ni (0.013) % Na (<0.01) % C _{Free} (0.085) % O (4.03) % K (<0.01) % Si _{Free} (0.538) % N (4.77) % Fe _{Total} (0.19) % Ti 0.025 % (Values in parenthesis are indicative values)	100 g

Special materials

Code	Product	Unit
Cement		
BCS CRM353	Sulphate resist portland cement Certified values SiO ₂20.5 % Mn ₂ O ₃0.23 % P ₂ O ₅0.077 % Al ₂ O ₃3.77 % CaO64.8 % SO ₃2.25 % TiO ₂0.16 % MgO2.42 % SrO0.23 % Fe ₂ O ₃4.82 % Na ₂ O (Acid Sol)0.1 % Cl(0.01) % Cr ₂ O ₃(0.02) % K ₂ O (Acid Sol)0.49 % (Values in parenthesis are indicative values)	100 g
BCS CRM354	White Portland Cement Certified values SiO ₂21.8 % Mn ₂ O ₃0.057 % P ₂ O ₅0.12 % Al ₂ O ₃4.85 % CaO70 % SO ₃2.25 % TiO ₂0.04 % MgO0.42 % SrO0.11 % Fe ₂ O ₃0.3 % Na ₂ O (Acid Sol)0.1 % Cl0.005 % Cr ₂ O ₃(0.003 %) K ₂ O (Acid Sol)0.11 % (Values in parenthesis are indicative values)	100 g
BCS CRM372/1	Ordinary portland cement Certified values SiO ₂20.3 % Mn ₂ O ₃0.074 % P ₂ O ₅(0.07) % Al ₂ O ₃5.37 % CaO65.3 % SO ₃2.95 % TiO ₂0.27 % MgO1.31 % SrO(0.05) % Fe ₂ O ₃3.42 % Na ₂ O (Acid Sol)0.1 % Cl(0.008) % Cr ₂ O ₃(0.01) % K ₂ O (Acid Sol)0.75 % (Values in parenthesis are indicative values)	100 g

Special materials

High tech ceramics

Code	Product	Unit
ERM-ED101	Silicon nitride powder (BAM-S001) Certified values Al469 ± 12 mg/kg Mg4.3 ± 0.4 mg/kg N38.1 ± 0.2 % Ca14.1 ± 0.5 mg/kg Na7.59 ± 0.27 mg/kg O(1.91 ± 0.07) % Co43.5 ± 0.8 mg/kg W41.3 ± 1.3 mg/kg β-phase7.43 ± 0.09 % Fe79.5 ± 1.3 mg/kg C0.162 ± 0.024 % (Values in parenthesis are indicative values)	50 g
BAM-S003	Silicon carbide powder Certified values Al372 ± 20 mg/kg Mn1.44 ± 0.17 mg/kg O910 ± 35 mg/kg B63 ± 7 mg/kg Na17.7 ± 0.8 mg/kg N(93 ± 22) mg/kg Ca29.4 ± 1.8 mg/kg Ni32.9 ± 2.7 mg/kg SiO ₂ free(600 ± 148) mg/kg Cr3.5 ± 0.4 mg/kg Ti79 ± 4 mg/kg Si _{free}(481 ± 223) mg/kg Cu1.5 ± 0.4 mg/kg V41.4 ± 2.8 mg/kg C _{total}29.89 ± 0.07 % Fe149 ± 10 mg/kg Zr25.2 ± 2.0 mg/kg Mg6.3 ± 0.6 mg/kg C _{free}493 ± 79 mg/kg (Values in parenthesis are indicative values)	bottle
NIST-154C	Titanium dioxide - Purity This Standard Reference Material® (SRM®) is intended primarily for use in the evaluation of techniques employed in the assay of titanium dioxide in the paint and ceramic industries. Certified value Titanium dioxide99.951 % Informational values for impurities	90 g
BCS CRM201A	Nepheline Syenite Certified values SiO ₂(57.3) % MgO(0.025) % Mn ₂ O ₃(0.007) % Al ₂ O ₃(23.54) % Na ₂ O(7.53) % SrO(0.43) % TiO ₂(0.05) % K ₂ O(8.9) % L.O.I.(0.76) % Fe ₂ O ₃(0.12) % P ₂ O ₅(0.025) % CaO(1.07) % BaO(0.37) % (Values in parenthesis are indicative values)	100 g
BCS CRM202A	Plaster (Gypsum) Certified values SiO ₂(1.38) % CaO(37.4) % P ₂ O ₅(<0.01) % Al ₂ O ₃(0.33) % MgO(0.39) % SrO(0.33) % TiO ₂(0.03) % Na ₂ O(<0.03) % SO ₃(53) % Fe ₂ O ₃(0.1) % K ₂ O(0.1) % (Values in parenthesis are indicative values)	100 g

Coals

Code	Product		Unit
BCS CRM203A	Talc		100 g
	SiO ₂(59.7) %	CaO(0.25) %	P ₂ O ₅(0.13) %
	Al ₂ O ₃(0.3) %	MgO(32.08) %	L.O.I.(6.78) %
	TiO ₂(<0.01) %	Na ₂ O(0.02) %	
	Fe ₂ O ₃(0.22) %	K ₂ O(0.005) %	
	(Values in parenthesis are indicative values)		
BCS CRM204A	Zircon		100 g
	SiO ₂(37.6) %	CaO(0.15) %	P ₂ O ₅(0.77) %
	Al ₂ O ₃(0.74) %	MgO(0.012) %	SnO ₂(1.69) %
	TiO ₂(2.22) %	Na ₂ O(0.014) %	ZrO ₂ +HfO ₂(53.8) %
	Fe ₂ O ₃(0.18) %	K ₂ O(0.017) %	L.O.I.(0.5) %
	(Values in parenthesis are indicative values)		

Multielement glass for XRF analysis – type A – type B

BAM-S005A	Multielement glass for XRF analysis (soda lime glass) - disc 39 mm x 5 mm (26 g - 30 g)	disc
Certified values		
	Arsenic (III) oxide	132 ± 8 mg/kg
	Barium oxide.....	115 ± 9 mg/kg
	Cadmium oxide.....	..62 ± 4 mg/kg
	Cerium (IV) oxide	105 ± 6 mg/kg
	Chloride	247 ± 33 mg/kg
	Cobalt oxide.....	.49.4 ± 2.4 mg/kg
	Chromium (III) oxide	15.6 ± 2.4 mg/kg
	Copper (II) oxide	112 ± 5 mg/kg
	Iron (III) oxide	422 ± 11 mg/kg
	Manganese (II) oxide	124 ± 5 mg/kg
	Molybdenum (VI) oxide	343 ± 12 mg/kg
	Nickel (II) oxide59.0 ± 2 mg/kg
	Lead (II) oxide	202 ± 8 mg/kg
	Antimony (III) oxide	132 ± 7 mg/kg
	Selenium.....	19.6 ± 1.7 mg/kg
(Values in parenthesis are indicative values)		
BAM-S005B	Multielement glass for XRF analysis (soda lime glass) - disc 39 mm x 5 mm (26 g - 30 g)	disc
Certified values		
	Arsenic (III) oxide	132 ± 8 mg/kg
	Barium oxide.....	115 ± 5 mg/kg
	Cadmium oxide.....	..62 ± 3 mg/kg
	Cerium (IV) oxide	105 ± 5 mg/kg
	Chloride	247 ± 24 mg/kg
	Cobalt oxide.....	.49.4 ± 2.3 mg/kg
	Chromium (III) oxide	15.2 ± 1.2 mg/kg
	Copper (II) oxide	112 ± 4 mg/kg
	Iron (III) oxide	422 ± 10 mg/kg
	Manganese (II) oxide	124 ± 5 mg/kg
	Molybdenum (VI) oxide	343 ± 12 mg/kg
	Nickel (II) oxide59.0 ± 1.9 mg/kg
	Lead (II) oxide	202 ± 7 mg/kg
	Antimony (III) oxide	132 ± 6 mg/kg
(Values in parenthesis are indicative values)		

Coals

Code	Product		Unit
BCR-182	Steam coal - Elements, ash and gross calorific value		15 g
	A medium volatile coal of high ash content.		
	Certified values		
	C732.9 g/kg	Mn195 mg/kg	Zn.....33.3 mg/kg
	Cd0.057 mg/kg	N.....16.36 g/kg	Ash.....122.7 g/kg
	Cl3.70 g/kg	Se.....0.68 mg/kg	GCV*.....29.68 MJ/kg
	Hg0.040 mg/kg	V.....24.3 mg/kg	
	Indicative values for Al, As, B, Br, Ce, Co, Cr, Cu, Fe, H, K, La, Mn, Na, Ni, Pb, Rb, Sc, Th, Ti, Vol. matter		
	*Gross calorific value		
BCR-331	Low volatile steam coal - Sulphur		20 g
	Certified value		
	S4.99 g/kg		
BCR-332	High volatile industrial coal - Sulphur		20 g
	Certified value		
	S9.61 g/kg		
BCR-333	Coking steam coal - Sulphur		20 g
	Certified value		
	S13.44 g/kg		

Coals

Code	Product	Unit
BCR-335	Flame coal - Sulphur Certified value S..... 50.8 g/kg	20 g
BCR-336	High volatile steam coal - Sulphur Certified value S..... 32.90 g/kg	20 g
BCR-460	Coal - Fluorine Certified value F..... 225 mg/kg Indicative value for Cl	40 g
NCS FC82007	Coal - Fluorine Certified value F..... 248±12 µg/g	50 g
NCS FC82008	Coal - Fluorine Certified value F..... 864±16 µg/g	50 g
NCS FC82009	Coal - Fluorine Certified value F..... 1496±20 µg/g	50 g
NIST-2682B	Coal (sub-bituminous) - Sulphur, mercury and heat of combustion Certified values S..... 0.4917 % Hg..... 108.8 µg/kg Indicative values for chlorine, ash content, gross calorific value	50 g
NIST-2683B	Coal (bituminous) - Sulphur Certified value S..... 1.955 % Indicative value for Furnace ash	50 g
NIST-2684B	Coal (bituminous) - Sulphur, mercury and heat of combustion Certified values S..... 3.076 % Heat of combustion 28.56 MJ/kg Hg 97.4 µg/kg Indicative values for Ash content, Gross calorific value	50 g
NIST-2692b	Coal (bituminous) - Sulphur, mercury and gross calorific value Certified values S..... 1.170 % ± 0.020 % Chlorine..... 1651 ± 28 mg/kg Hg 133.3 µg/kg ± 4.1 µg/kg Reference values Ash Content (mass fraction) 7.90 % Gross Calorific Value 32.81 MJ/kg	50 g
NIST-2693	Coal (bituminous) - Sulphur and mercury Certified values S..... 0.4571 % Hg..... 37.3 µg/kg Cl..... 369.6 mg/kg	50 g
NIST-2685B	Coal - Sulphur, mercury and heat of combustion This Standard Reference Material (SRM®) is intended primarily for use in the evaluation of techniques employed in the determination of sulfur, mercury, ash content, and calorific value (MJ·kg⁻¹) in coal and materials of a similar matrix. NIST-2685b consists of 50 g of bituminous coal ground to pass a 250 µm (60 mesh) sieve, homogenized, and packaged in an amber glass bottle. Certified values S..... 4.730 % Cl..... 517 mg/kg Hg 146.2 µg/kg Indicative values for the ash content, gross calorific value and elements	50 g
NIST-1632C	Coal (bituminous) - Elements Certified values H 5.11% Ba 41.1 mg/kg Selenium 1.326 mg/kg S 1.462 % Co 3.48 mg/kg Sodium 298.8 mg/kg Potassium 0.1100 % Mn 13.04 mg/kg Strontium 63.8 mg/kg Chlorine 0.1139 % Hg 0.0938 mg/kg Thorium 1.40 mg/kg Sb 0.461 mg/kg Rb 7.52 mg/kg Zinc 12.1 mg/kg Indicative values for Carbon (total), N , Si, Al, As, B, Br, Ca, Cd, Ce, Cr, Cs, Co, Eu, Fe, Hf, Mg, Ni, Pb, Sm, Sc,Ti, U, V, Ash	50 g

Code	Product	Unit
NIST-1635	Coal - Trace elements Certified values As 0.42 mg/kg Cd 0.03 mg/kg Cr 2.5 mg/kg Cu 3.6 mg/kg F 25.9 mg/kg Fe 0.239 % Mn 21.4 mg/kg Ni 1.74 mg/kg Pb 1.9 mg/kg S 0.3616 % Se 0.9 mg/kg Th 0.62 mg/kg U 0.24 mg/kg V 5.2 mg/kg Zn 4.7 mg/kg Indicative values for Al, Ce, Co, Eu, Ga, Hf, Hg, Na, Sb, Sc, Ti	75 g
NCS FC28001J	Coal (bitumite) - Elements and properties Coal type: Bitumite Certified values C 78.35 % H 4.28 % N 1.36 % Total S 0.50 % Ash 9.63 % Volatile matter 24.43 % Relative density (20°C) 1.41 Gross calorific value 31.47 MJ/kg	50 g
NCS FC28002H	Coal - Elements and properties Coal type: Bitumite Certified values C 57.84 % H 3.77 % N 1.08 % Total S 1.5 % Ash 28.69 % Volatile matter 27.27 % Relative density (20°C) 1.57 Calorific value 23.63 MJ/kg	50 g
NCS FC28003F	Coal - Elements and properties Coal type: Anthracite Certified values C 78.22 % H 0.95 % N 0.23 % Total S 0.28 % Ash 16.29 % Volatile matter 6.32 % Relative density (20°C) 1.95 Calorific value 26.44 MJ/kg	50 g
NCS FC28004E	Coal - Elements and properties Coal type: Anthracite Certified values C 66.8 % H 1.43 % N 0.74 % Total S 0.99 % Ash 28.13 % Volatile matter 4.81 % Relative density (20°C) 1.89 Calorific value 23.8 MJ/kg	50 g
NCS FC28005C	Coal (anthracite) - Elements and properties Coal type: Anthracite Certified values C 80.73 % H 2.74 % N 1.06 % S 1.86 % Ash 11.63 % Volatile matter 6.94 % Relative density (20°C) 1.57 Gross calorific value 30.84 MJ/kg	50 g
NCS FC28008D	Coal - Elements and properties Coal type: Bitumite Certified values C 66.71 % H 3.63 % N 1.19 % Total S 3.54 % Ash 21.65 % Volatile matter 18.34 % Relative density (20°C) 1.54 Calorific value 27.06 MJ/kg	50 g
NCS FC28009E	Coal - Elements and properties Coal type: Bitumite Certified values C 54.1 % H 3.02 % N 0.97 % Total S 4.09 % Ash 32.78 % Volatile matter 17.14 % Relative density (20°C) 1.73 Calorific value 21.5 MJ/kg	50 g
NCS FC28011C	Coal - Elements and properties Certified values C 72.79 % H 2.12 % N 0.52 % Total S 2.20 % Ash 20.41 % Volatile matter 6.45 % Density (20°C) 1.72 Gross calorific value 27.27 MJ/kg	50 g

Coals

Code	Product	Unit
NCS FC28017	Coal - Elements and properties Coal type: Anthracite Certified values C 73.64 % Ash 21.01 % H 0.91 % Volatile matter 5.9 % N 0.23 % Relative density (20 °C) 1.97 Total S..... 0.21 % Calorific value..... 25.01 MJ/kg	50 g
NCS FC28101	Coal - Elements and properties Coal type: Anthracite Certified values C 90.27 % Ash 3.95 % H 3.01 % Volatile matter 6.64 % N 0.60 % Relative density (20 °C) 1.47 Total S..... 0.20 % Calorific value..... 34.34 MJ/kg	50 g
NCS FC28102	Coal - Elements and properties Coal type: Anthracite Certified values C 87.47 % Ash 6.46 % H 2.86 % Volatile matter 7.90 % N 0.60 % Relative density (20 °C) 1.5 Total S..... 0.19 % Calorific value..... 33.1 MJ/kg	50 g
NCS FC28103	Coal - Elements and properties Coal type: Anthracite Certified values C 81.55 % Ash 10.51 % H 3.33 % Volatile matter 9.45 % N 1.30 % Relative density (20 °C) 1.47 Total S..... 0.36 % Calorific value..... 31.8 MJ/kg	50 g
NCS FC28104	Coal - Elements and properties Coal type: Anthracite Certified values C 81.60 % Ash 10.09 % H 3.52 % Volatile matter 11 % N 1.34 % Relative density (20 °C) 1.45 Total S..... 0.41 % Calorific value..... 32.04 MJ/kg	50 g
NCS FC28105	Coal - Elements and properties Coal type: Anthracite Certified values C 81.54 % Ash 9.61 % H 3.70 % Volatile matter 12.21 % N 1.16 % Relative density (20 °C) 1.43 Total S..... 1.06 % Calorific value..... 32.31 MJ/kg	50 g
NCS FC28106	Coal - Elements and properties Coal type: Bitumite Certified values C 79.09 % Ash 8.56 % H 4.95 % Volatile matter 31.92 % N 1.38 % Relative density (20 °C) 1.35 Total S..... 1.72 % Calorific value..... 32.98 MJ/kg	50 g
NCS FC28107	Coal - Elements and properties Coal type: Bitumite Certified values C 79.89 % Ash 10.41 % H 3.80 % Volatile matter 15.30 % N 1.12 % Relative density (20 °C) 1.43 Total S..... 0.67 % Calorific value..... 31.64 MJ/kg	50 g
NCS FC28108	Coal - Elements and properties Coal type: Bitumite Certified values C 72.94 % Ash 13.68 % H 4.46 % Volatile matter 30.84 % N 1.26 % Relative density (20 °C) 1.42 Total S..... 0.57 % Calorific value..... 29.9 MJ/kg	50 g

Code	Product	Unit
NCS FC28109	Coal - Elements and properties Coal type: Anthracite Certified values C 79.42 % Ash 11.98 % H 3.28 % Volatile matter 11.30 % N 1.09 % Relative density (20 °C) 1.49 Total S 0.58 % Calorific value 30.66 MJ/kg	50 g
NCS FC28110	Coal - Elements and properties Coal type: Bitumite Certified values C 75.96 % Ash 8.42 % H 4.56 % Volatile matter 32.94 % N 1.33 % Relative density (20 °C) 1.41 Total S 0.87 % Calorific value 30.92 MJ/kg	50 g
NCS FC28111	Coal - Elements and properties Coal type: Bitumite Certified values C 60.24 % Ash 25.19 % H 3.73 % Volatile matter 28.39 % N 1.04 % Relative density (20 °C) 1.57 Total S 1.28 % Calorific value 24.35 MJ/kg	50 g
NCS FC28112	Coal - Elements and properties Coal type: Bitumite Certified values C 78.75 % Ash 8.08 % H 5.01 % Volatile matter 33.70 % N 1.31 % Relative density (20 °C) 1.33 Total S 2.10 % Calorific value 33.04 MJ/kg	50 g
NCS FC28113	Coal - Elements and properties Coal type: Bitumite Certified values C 74.80 % Ash 7.06 % H 4.47 % Volatile matter 33.40 % N 1.02 % Relative density (20 °C) 1.41 Total S 0.27 % Calorific value 30.03 MJ/kg	50 g
NCS FC28114	Coal - Elements and properties Coal type: Bitumite Certified values C 76.36 % Ash 4.66 % H 4.54 % Volatile matter 33.07 % N 1.08 % Relative density (20 °C) 1.4 Total S 0.20 % Calorific value 30.73 MJ/kg	50 g
NCS FC28115	Coal - Elements and properties Coal type: Bitumite Certified values C 77.44 % Ash 6.38 % H 4.42 % Volatile matter 32.22 % N 1.21 % Relative density (20 °C) 1.41 Total S 0.42 % Calorific value 31.05 MJ/kg	50 g
NCS FC28116	Coal - Elements and properties Coal type: Bitumite Certified values C 78.68 % Ash 6.08 % H 4.59 % Volatile matter 32.34 % N 1.34 % Relative density (20 °C) 1.39 Total S 0.54 % Calorific value 31.82 MJ/kg	50 g
NCS FC59001	Coke - Sulphur, ash and volatile matter Certified values Total Sulphur 0.63±0.03 % Ash 7.22±0.06 % Volatile matter 1.39±0.08 %	60 g
NCS FC59002	Coke - Sulphur, ash and volatile matter Certified values Total Sulphur 0.47±0.02 % Ash 12.62±0.08 % Volatile matter 1.5±0.2 %	60 g
NCS FC82004	Coal - Chlorine Certified value Cl 0.010±0.002 %	50 g

Coals

Code	Product	Unit
NCS FC82005	Coal - Chlorine Certified value Cl.....0.057±0.003 %	50 g
NCS FC82006	Coal - Chlorine Certified value Cl.....0.110±0.006 %	50 g
NCS FC82001	Coal - Arsenic and phosphorous Certified value As.....15±1 µg/g P0.031±0.002%	50 g
NCS FC82002	Coal - Arsenic and phosphorous Certified value As.....34±2 µg/g P0.007±0.001%	50 g
NCS FC82003	Coal - Arsenic and Phosphorus Certified value As.....51±3 µg/g P0.092±0.005%	50 g
NCS FC62001	Bituminous Coal Certified values M _{ad}4.51% V _{ad}22.89% S _{t,d}1.53% A _{ad}21.17% V _d23.97% Q _{gr,ad}23.91 MJ/kg A _d22.17% S _{t,ad}1.46% Q _{gr,d}25.04 MJ/kg M _{ad} , A _{ad} , V _{ad} , S _{t,ad} and Q _{gr,ad} mean moisture, ash, volatile material, total sulfur and calorific value of air dried sample; A _d , V _d , S _{t,d} and Q _{gr,d} mean ash, volatile material, total sulfur and calorific value of dried sample.	20 g
NCS FC62002	Bituminous Coal Certified values M _{ad}2.27% V _{ad}4.21% S _{t,d}2.95% A _{ad}25.82% V _d4.31% Q _{gr,ad}23.48 MJ/kg A _d26.42% S _{t,ad}2.88% Q _{gr,d}24.02 MJ/kg M _{ad} , A _{ad} , V _{ad} , S _{t,ad} and Q _{gr,ad} mean moisture, ash, volatile material, total sulfur and calorific value of air dried sample; A _d , V _d , S _{t,d} and Q _{gr,d} mean ash, volatile material, total sulfur and calorific value of dried sample.	20 g
NIST-2775	Foundry coke - Sulphur Certified value S.....0.5816 % Indicative values for ash, volatile matter, carbon, hydrogen and nitrogen	50 g
NIST-2776	Foundry coke - Sulphur, ash and volatile matter Certified values S.....0.825 ± 0.016 % Indicative values for ash, volatile matter, carbon, hydrogen, nitrogen	50 g
NIST-2718	Green petroleum coke - Trace elements Certified values Al.....16.5 mg/kg Fe290 mg/kg S.....47032 mg/kg Ca174 mg/kg Ni139.1 mg/kg V.....302 mg/kg Indicative values for C, Co, H, N, Na, Si, Ash, Volatile matter, Gross calorific value	50 g
NIST-2719	Calcinated petroleum coke - Trace elements Certified values Al.....58.9 mg/kg Fe201.6 mg/kg S.....8877 mg/kg Ca57.7 mg/kg Ni204 mg/kg V.....58.6 mg/kg Indicative values for C, Co, H, N, Na, Si, Ash, Volatile content, Gross calorific value	50 g



LGC Standards, Queens Road, Teddington, Middlesex TW11 0LY, UK
www.lgcstandards.com

Excellence through measurement